

KREATECH DIAGNOSTICS CATALOGUE 2013 - 2014



FISH4U™



NEW CUSTOM PROBE SERVICE FOR YOU FROM KREATECH

Our flexible custom probe service FISH4U™ gives you access to your probe of choice. We can develop completely new probe designs to your specifications or you can just request an existing DNA probe from our portfolio to be labeled with an alternative color from our range.

Custom made probes

- Leading-edge probe design
- Designed to meet your specifications

REPEAT-FREE™ POSEIDON™ FISH DNA Probes

- A clearer background
- Greater signal intensity

Check out our custom probe service

Contact customer service by mail:

customerservice@kreatech.com

or visit our website:

www.kreatech.com

Disclaimer FISH4U:

The Customer acknowledges that no licenses or other rights are provided to the Customer under any patents, patent applications, trade secrets or other proprietary rights of Kreatech, except the right to conduct research.

Kreatech retains ownership of the Product and Kreatech is free to commercialize the Product to any third party except if customer has rights or desires to receive rights to such probe as designed, developed and produced by Kreatech.

If the Customer desires to use the Product for commercial purposes and/or acquire exclusive rights in the Product, the Customer agrees, in advance of such use, to negotiate in good faith with Kreatech to establish the terms of such an agreement.

WELCOME TO THE KREATECH DIAGNOSTICS 2013-2014 CATALOGUE

As a leader in molecular diagnostics we can offer you the broadest portfolio of Fluorescent In Situ Hybridization (FISH) probes available. Our REPEAT-FREE™ POSEIDON™ FISH probes and the unique ULS™-labeling technology are used for diagnostic- and (clinical) research applications worldwide.

We continuously strive to increase our level of customer- and technical services, as we are committed to delivering high-performance innovative products. Building on collaborations with leading scientists and global partners we continue to invest in advancing diagnostics and science by continuing to expand the portfolio of our REPEAT-FREE™ POSEIDON™ FISH probes.

To meet the ever-growing demand of high-quality DNA FISH probes throughout the scientific and diagnostic communities, we have

made our expertise in developing REPEAT-FREE™ POSEIDON™ FISH probes available through our custom design program – FISH4U™. This custom probe service allows you to get access to DNA FISH probes that are tailored to your application.

Check our website at www.kreatech.com on a regular basis for up to date information on our products. On our website you can also find the KREATECH chromosome finder: an easy-to-use tool that will guide you through our range of FISH-probes.

We look forward to serving you in the coming years. Your trusted partner in molecular diagnostics and microarray labeling applications.

KREATECH Diagnostics




TABLE OF CONTENTS

| | Page |
|-------------------------------|---|
| | I Catalogue 2013-2014 – New Products |
| Molecular Cytogenetics | 7 Oncology – Hematology DNA Probes |
| | 48 Oncology – Solid Tumor DNA Probes |
| | 74 Oncology – Chromogenic <i>In Situ</i> Hybridization |
| | 76 Preimplantation Genetic Screening |
| | 78 Prenatal DNA Probes |
| | 80 Microdeletion DNA Probes |
| | 92 Sub-Telomere DNA Probes |
| | 94 Satellite Enumeration DNA Probes |
| | 96 Whole Chromosome DNA Probes |
| | 98 Arm Specific / Band Specific DNA Probes |
| | 100 Mouse DNA Probes |
| | 101 Cell Culture Media |
| | 103 Pretreatment Kits / Reagents |
| | 104 Equipment |
| ULS Labeling Products | 107 Array CGH |
| | 108 DNA Products |
| | 109 Gene Expression |
| | 110 MicroRNA |
| | 111 General Nucleic Acid Labeling |
| | 113 Protein Labeling / Antibody Labeling Kits |
| | 114 Supporting Products |
| Appendix | 115 Product Index CE marking & ISO certificates Ordering Information |

REPEAT-FREE POSEIDON FISH DNA PROBES

Kreatech Diagnostics offers the next generation REPEAT-FREE™ POSEIDON™ DNA probes developed with the use of the REPEAT-FREE (RF) technology.

This technology is based on a proprietary subtractive hybridization specifically removing all repetitive elements which are dispersed throughout the human genome. Eliminating these repeat sequences leads to a more specific binding kinetics and makes the need to use C₀t1 DNA for pre-annealing obsolete. In addition, our DNA probes are labeled with Kreatech's proprietary ULS™ linkage system leading to a very homogenous and consistent labeling of fluorescent dyes.

This results in FISH probes that are brighter and give a clearer background.

The signal to noise ratio of RF POSEIDON FISH probes in comparison to standard probes today is 1.5 to 3 times greater.

Our recommended Pretreatment Kits with Ready-to-Use reagents makes the pretreatment of different tissue and cytological samples easier and guarantees optimal results in combination with the RF POSEIDON FISH probes.


Traditional FISH is very time consuming with standard hybridizations being on average 16 hours in duration. With the use of RF POSEIDON FISH probes, hybridization times may be cut 75% to 4 hours.

This hybridization flexibility coupled with optimal pretreatment and the probes being available in a Ready-to-Use format eliminates time consuming steps and optimizes workflow.

The POSEIDON™ Probe Maps

All locus specific probes are provided with a map showing the genome region and/or genes the probe cover. The STS markers indicating the most proximal and distal end covered by the probe are added as well. Those markers can be used to

describe in situ hybridization results in accordance with the International System for Human Cytogenetic Nomenclature (ISCN, 2013). For more information see e.g. Genome Browser at genome.ucsc.edu or Ensemble at www.ensembl.org



Fluorophores and filter recommendations

Please check out the website www.kreatech.com for recommended filters covering the excitation and emission wavelengths as specified below.

When ordering a filter from a supplier please make sure to either have a cube or slider available for your specific brand and type of microscope.

Single Filters

| Fluorophore | Color | Ex/Em | Recommended Filter | Excitation Filter | Dichromatic Mirror | Emission Filter |
|--------------------|---|---------|--------------------|-------------------|--------------------|-----------------|
| DAPI | DAPI | | DAPI / UV | 360/20 | 400 LP | 425 LP |
| PlatinumBright™405 |  Dark Blue | 410/455 | Blue | 405/10 | 425 LP | 460/50 |
| PlatinumBright415 |  Blue | 429/470 | Aqua | 436/20 | 455 LP | 480/30 |
| PlatinumBright495 |  Green | 495/517 | FITC / Green | 480/30 | 505 LP | 535/40 |
| PlatinumBright505 |  Green | 500/528 | FITC / Green | 480/30 | 505 LP | 535/40 |
| PlatinumBright530 |  Gold | 531/561 | Gold | 520/18 | 545 LP | 565/20 |
| PlatinumBright547 |  Light Red | 547/565 | Cy3 | 535/50 | 565 LP | 610/75 |
| PlatinumBright550 |  Red | 550/580 | TRITC | 540/25 | 565 LP | 605/55 |
| PlatinumBright570 |  Red | 570/591 | Texas Red | 560/40 | 595 LP | 645/75 |
| PlatinumBright590 |  Dark Red | 587/612 | Texas Red | 560/40 | 595 LP | 645/75 |
| PlatinumBright647 |  Far Red | 647/665 | Cy5 / Far red | 620/60 | 660 LP | 700/75 |

Dual/Triple Filters

| Fluorophore | Color | Recommended Filter | Excitation Filter | Dichromatic Mirror | Emission Filter |
|-------------------------------|----------------|----------------------------|----------------------------|------------------------------|----------------------------|
| PlatinumBright495/550 | Green/Red | Dual FITC/Cy3 | 485/10 563/12 | 500-550 580 LP | 521/20 595 LP |
| PlatinumBright 415/495/550 | Blue/Green/Red | DAPI/FITC/Cy3 | 430/12 490/12 562/12 | 446-478 503-548 579 LP | 463/10 523/15 585 LP |
| DAPI / PlatinumBright 495/550 | DAPI/Green/Red | DAPI/FITC/ TRITC or Cy3 | 396/8 482/8 554/14 | 440-472 498-540 573 LP | 450/10 516/20 604/22 |

CATALOGUE 2013-2014 – NEW PRODUCTS

Table of Contents

| <i>Description</i> | <i>Cat#</i> | <i>Page</i> |
|--|-------------|-------------|
| ONCOLOGY - HEMATOLOGY DNA PROBES | | II |
| CML | | II |
| ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | KBI-10013 | II |
| AML | | II |
| ON MECOM / RUNX1 t(3;21) Fusion | KBI-10310 | II |
| NUP98 (11p15) Break | KBI-10311 | III |
| ALL | | III |
| ON MLL/AFF1 t(4;11) Fusion | KBI-10404 | III |
| MM | | IV |
| ON MAFB/IGH@ t(14;20) Fusion | KBI-10510 | IV |
| ONCOLOGY – SOLID TUMOR DNA PROBES | | V |
| Lung Cancer | | V |
| ON ROS1 (6q22) Break | KBI-10752 | V |
| ON RET (10q11) Break | KBI-10753 | V |
| FGFR1 (8p11) / SE 8 (D8Z1) | KBI-12754 | VI |
| Sarcoma | | VI |
| ON COL1A1/PDGFB t(17;22) DC, S-Fusion | KBI-10742 | VI |
| ON EWSR1/NFATC2 t(20;22) DC, S-Fusion | KBI-10751 | VII |
| Lymphoma (tissue) | | VII |
| ON MYC (8q24) Break, TC (tissue) | KBI-10749 | VII |
| ON BCL2/IGH@ t(14;18) Fusion (tissue) | KBI-10755 | VIII |
| MICRODELETION DNA PROBES | | VIII |
| MD GATA4 (8p23) / SE 8 | KBI-40118 | VIII |



* All Microdeletion DNA Probes are available in a 5 or a 10 test format.

CML

ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES

Chronic Myeloid Leukemia (CML) is characterized by the formation of the BCR/ABL fusion gene as a result of the reciprocal translocation t(9;22)(q34;q11). The BCR/ABL fusion gene is found on the derivative chromosome 22, called the Philadelphia (Ph) chromosome. The same translocation is also observed in Acute Lymphocytic Leukemia (ALL) and in Acute Myeloid Leukemia (AML). This chimeric BCR/ABL gene encodes a constitutively activated protein tyrosine kinase with profound effects on cell cycle, adhesion, and apoptosis. Understanding this process has led to the development of the drug imatinib mesylate (Gleevec™). Breakpoints in the BCR gene region can occur in different regions, predominately in a major breakpoint cluster region (M-BCR) but can also occur in a minor breakpoint cluster region (m-BCR) or micro breakpoint cluster region (μ-BCR). Further research indicates that CML patients with different BCR-ABL transcripts respond differently to treatment with Gleevec.

Cat# KBI-10013 Mm-BCR/ABL t(9;22), DC, S-Fusion, ES



Mm-BCR/ABL probes hybridized to patient material showing t(9;22) with M-BCR (1F1r1R1G). Mm-BCR/ABL probes hybridized to patient material showing t(9;22) with m-BCR (2F1R1G).

Literature:
Dewald et al., 1998, Blood 91; 3357-3365.
Huntly et al., 2003, Blood 102; 1160-1168.
Sharma et al., 2010, Ann Hematol, 89, 241-7.
Tkachuk et al., 1990, Science 250, 559-56.
Kolomietz et al., 2001. Blood 97; 3581-3588.

| Ordering information | Gene Region | Tests | Cat# |
|---|-------------|-------|-----------|
| ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | green/red | 10 | KBI-10013 |



AML

ON MECOM / RUNX1 t(3;21) Fusion

The MECOM (EVI1, 3q26) / RUNX1 (AML1, 21q22) translocation, t(3;21), is consistently found in blastic crisis of chronic myelogenous leukemia (CML) and myelodysplastic syndrome-derived leukemias. The translocation produces MECOM / RUNX1 chimeric transcription factor and is thought to play important roles in acute leukemic transformation of hemopoietic stem cells.

The MECOM/RUNX1 t(3;21)(q26;q22) specific DNA Probe is optimized to detect the reciprocal translocation t(3;21) in a dual-color, dual-fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10310 MECOM / RUNX1 t(3;21) Fusion



MECOM / RUNX1 Fusion probe hybridized to patient material showing t(3;21) (2F1R1G). Images kindly provided by Dr. Mohr, Dresden.



| Ordering information | Gene Region | Tests | Cat# |
|---------------------------------|-------------|-------|-----------|
| ON MECOM / RUNX1 t(3;21) Fusion | red/green | 10 | KBI-10310 |

NUP98 (11p15) Break

Nucleoporin 98kDa gene (NUP98) rearrangements have been identified in a wide range of hematologic malignancies, including acute myeloid leukemia (AML), acute lymphoblastic leukemia (ALL), chronic myeloid leukemia in blast crisis (CML-bc), myelodysplastic syndrome (MDS) and bilineage/ biphenotypic leukemia. The NUP98 gene is highly promiscuous with regard to its recombination spectrum, as at least 28 different partner genes have been identified for NUP98 rearrangements, all forming in-frame fusion genes. Patients with NUP98 gene rearrangements have an aggressive clinical course and the outcome of treatment is disappointing.

The NUP98 (11p15) Break Probe is optimized to detect translocations involving the NUP98 gene region at 11p15 in a dual-color assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10311 NUP98 (11p15) Break



NUP98 (11p15) Break Probe hybridized to AML patient sample showing a rearrangement of 11p15 involving the NUP98 gene (1F1R1G).

Images kindly provided by Prof. Manuel R. Teixeira, Porto.

Literature:

Gough et al, 2011, Blood 118; 62 47-6257.
Nebral et al, 2005, Haematologica 90; 74 6-752.
Romana et al, 2006, Leukemia 20; 696-70 6.

| Ordering information | Gene Region | Tests | Cat# |
|----------------------|-------------|-------|-----------|
| NUP98 (11p15) Break | red/green | 10 | KBI-10311 |



ALL

ON MLL/AFF1 t(4;11) Fusion

The t(4;11) MLL/AFF1 is the most frequently (approximately 66% according to Meyer et al.) observed translocation involving the MLL gene resulting in Acute Lymphoblastic Leukemia (ALL). The MLL/AFF1 translocation results in the generation of fusion proteins MLL/AFF1 and AFF1/MLL; both seem to have leukemogenic properties. Furthermore, EVI1 (3q26) is one of the targets of the MLL oncoproteins, which increased expression correlates with unfavorable prognosis in Acute Myeloid Leukemia. Patients with ALL and the MLL/AFF1 translocation are associated with a high risk of treatment failure.

The MLL/AFF1 t(4;11) Fusion Probe is optimized to detect translocations involving the MLL and AFF1 gene regions at 4q21-22 and 11q23 in a dual-color, fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10404 MLL/AFF1 t(4;11) Fusion



MLL/AFF1 t(4;11) Fusion Probe. MLL/AFF1 t(4;11) Fusion Probe.
Standard t(4;11) 2 Fusion, 1 Red, ins(4;11) 1F2R1G.
1 Green (2F1R1G).

Literature

Harrison CJ et al, 2010, Br J Haem, 151; 132-142.
Arai S et al, 2011, Blood, 117; 6304 - 6314.
Meyer C et al, 2009, Leukemia, 23; 1490-1499.

| Ordering information | Gene Region | Tests | Cat# |
|----------------------------|-------------|-------|-----------|
| ON MLL/AFF1 t(4;11) Fusion | red/green | 10 | KBI-10404 |



MM

ON MAFB/IGH@ t(14;20) Fusion

The immunoglobulin heavy chain (IGH@) gene at 14q32 is an important cause of genetic deregulation in multiple myeloma (MM). Among the known fusion partners for the IGH gene, reciprocal translocation with the MAFB gene at 20q12 is relatively rare in MM (~2% occurrence). However, the MAFB/IGH@ t(14;20) translocation is associated with poor prognosis in multiple myeloma patients.

The MAFB/IGH@ t(14;20) Fusion Probe is optimized to detect the reciprocal translocation t(14;20) in a dual-color, dual-fusion assay on metaphase/interphase spreads and bone marrow cells.

Cat# KBI-10510 MAFB/IGH@ t(14;20) Fusion



The MAFB/IGH@ t(14;20) Fusion Probe hybridized to patient material showing a complex pattern with a t(14;20) translocation. Images kindly provided by Erasmus Medical Center, Rotterdam..

Literature:

Boersma-Vreugdenhil GR et al, 2004, Br J Haematol, 126, 355-363.
Bergsagel PL et al, 2005, JCO, 23, 6333-6338.

Ordering information

| | Gene Region | Tests | Cat# |
|------------------------------|-------------|-------|-----------|
| ON MAFB/IGH@ t(14;20) Fusion | red/green | 10 | KBI-10510 |



Lung Cancer

ON ROS1 (6q22) Break


Translocations involving the ROS1 (repressor of silencing 1) gene at chromosome 6q22 can increase expression of the gene by fusion with SLC34A2 (4p15), but also with other fusion partners. Elevated expression is observed in non-small cell lung cancer (NSCLC), where the success of tyrosine kinase-based therapeutics is based on inhibiting the activity of these fusion genes. The fusion of ROS1 to the GOPC (FIG; 6q22) gene, by deletion of a 240 kb DNA fragment, also results in activation of a fusion gene.

The ROS1 (6q22) Break Probe is optimized to detect translocations involving the ROS1 gene region at the 6q22 locus, as well as the 240 kb deletion forming the ROS1-GOPC fusion gene, in a dual-color assay on formalin- fixed paraffin-embedded tissue samples.

Cat# KBI-10752 ROS1 (6q22) Break



Hybridization of ROS1 (6q22) Break Probe (KBI-10752) to a tissue section harboring a ROS1 rearrangement.



Hybridization of ROS1 (6q22) Break Probe (KBI-10752) to a cell line harboring a GOPC-ROS1 rearrangement (deletion of red signal).

Literature:

Charest et al., Genes Chromosomes Cancer, 2003, 37: 58-71.
Rikova et al., Cell, 2007, 131: 1190-120.
Rimkunas et al., Clin. Can. Res., 2012, 18: 4449-4457.
Takeuchi et al., Nat. Med., 2012, 18: 378-381.
Gu et al., PLoS, 2011, 6: e15640.



| Ordering information | Gene Region | Tests | Cat# |
|----------------------|-------------|-------|-----------|
| ON ROS1 (6q22) Break | green/red | 10 | KBI-10752 |

ON RET (10q11) Break

Pericentric inversion of chromosome 10 involving the RET (ret proto-oncogene) gene at chromosome 10q11 is known to increase expression of the RET gene by fusion with KIF5B (10p11). Translocations with other fusion partners have also been described. Elevated expression of RET is observed in non-small cell lung cancer (NSCLC), in which the function of tyrosine kinase-based therapeutics is based upon the inhibition of such fusion proteins. Translocations involving RET have also been described in thyroid carcinomas.

The RET (10q11) Break Probe is optimized to detect translocations involving the RET gene region at the 10q11 locus..

Cat# KBI-10753 RET (10q11) Break



Hybridization of RET (10q11) Break Probe (KBI-10753) to a tissue section.

Literature:

Chen et al., Cancer Genet Cytogenet, 2007, 178: 128-134.
Kohno et al., Nat Med, 2012, 18: 375-377.
Takeuchi et al., Nat Med, 2012, 18: 378-381.



| Ordering information | Gene Region | Tests | Cat# |
|----------------------|-------------|-------|-----------|
| ON RET (10q11) Break | red/green | 10 | KBI-10753 |

FGFR1 (8p11) / SE 8 (D8Z1)

Amplification of the fibroblast growth factor receptor type 1 gene (FGFR1) has been observed in numerous cancer types including lung cancer (especially squamous cell carcinoma) and breast cancer. With the development of new therapeutic strategies, FGFR1 amplification could act as a valuable biomarker for R&D and provide an attractive tool for clinical stratification².

The FGFR1 (8p11) / SE8 Amplification probe is optimized to detect amplification involving the FGFR1 gene region at 8p11 in a dual-color assay on paraffin embedded tissue sections.

Cat# KBI-12754 / KBI-14754 FGFR1 (8p11) / SE 8 (D8Z1)



FGFR1 gene locus amplification in FFPE NSCLC tissue.

Literature:

Weiss et al., 2010, Sci. Transl. Med. 2(62): 62ra93.
Brooks et al., 2012, Clin. Cancer Res. 18(7): 1855-62.

Ordering information

| | Gene Region | Tests | Cat# |
|----------------------------|-------------|-------|-----------|
| FGFR1 (8p11) / SE 8 (D8Z1) | red/green | 20 | KBI-12754 |
| FGFR1 (8p11) / SE 8 (D8Z1) | red/green | 50 | KBI-14754 |



Sarcoma

ON COL1A1/PDGFB t(17;22) DC, S-Fusion

The diagnosis of primary soft tissue and bone tumors is often challenging as they are relatively rare. The misdiagnosis between dermatofibroma (DF) and dermatofibrosarcoma protuberans (DFSP) or giant cell fibroblastoma (GCF) might result in improper primary management. DFSP and GCF have in most cases diagnosed today a translocation involving the COL1A1 (collagen, type I, alpha 1) gene at 17q21 and the PDGFB (platelet-derived growth factor beta polypeptide) gene at 22q13. Also, a supernumerary ring chromosome derived from the translocation r(17;22) can be present.

The COL1A1/PDGFB t(17;22) Dual-Color Single-Fusion Probe is optimized to detect the t(17;22)(q21;q13) involving the COL1A1 (17q21) and PDGFB (22q13) gene regions in dual-color, single-fusion assay on paraffin embedded tissue sections.

Cat# KBI-10742 COL1A1/PDGFB t(17;22) DC, S-Fusion



Interphase FISH result of ON COL1A1/PDGFB Fusion probe hybridized to dermatofibrosarcoma protuberans tumor tissue. Co-localisation and amplification of the fusion gene is well visible.

Literature:

Maire et al, 2007, Arch Dermatol, 143; 203-210.
Labropoulos et al, 2007, Biologics, 1; 347-353.
Patel et al, 2008, Hum Path, 39; 184-193.
Sandberg, 2003, Cancer Genet Cytogenet, 140; 1-12.

Ordering information



| | Gene Region | Tests | Cat# |
|---------------------------------------|-------------|-------|-----------|
| ON COL1A1/PDGFB t(17;22) DC, S-Fusion | red/green | 10 | KBI-10742 |

ON EWSR1/NFATC2 t(20;22) DC, S-Fusion

Ewing's sarcoma is the second most frequent primary bone cancer. In most cases a translocation involving the EWSR1 gene at 22q12 and the FLI1 gene at 11q24 is observed. Several other translocation partners of the ETS gene family can also be involved. The first non-ETS family translocation partner described is the NFATC2 gene (nuclear factor of activated T-cells, cyto-plasmic, calcineurin-dependent 2) at 20q13.

The EWSR1/NFATC2 single fusion probe is best used to analyze the specific trans-locations of the EWSR1 and NFATC2 gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat# KBI-10751 EWSR1/NFATC2 t(20;22) DC, S-Fusion



Interphase FISH result of the EWSR1/NFATC2 t(20;22) DC, S-Fusion probe.

Literature:

Szuhai et al, 2009, Clin Cancer Res, 15; 2259-2268.
Zucman-Rossi et al, 1998, PNAS, 95; 11786-11791.
Bernstein et al, 2006, Oncologist, 11; 503-519.

| Ordering information | Gene Region | Tests | Cat# |
|---------------------------------------|-------------|-------|-----------|
| ON EWSR1/NFATC2 t(20;22) DC, S-Fusion | red/green | 10 | KBI-10751 |

Lymphoma (tissue)



ON MYC (8q24) Break, TC (tissue)

Rearrangements of the proto oncogene C-myc (or MYC) have been consistently found in tumor cells of patients suffering from Burkitt's lymphoma. In cases with the common t(8;14) chromosomal translocation, the c-myc gene is translocated to chromosome 14 and rearranged with the immunoglobulin heavy chain genes; the breakpoint occurs 5' to the c-myc gene and may disrupt the gene itself separating part of the first untranslated exon from the remaining two coding exons. In Burkitt's lymphoma showing the variant t(2;8) or t(8;22) translocations, the genes coding for the k and l immunoglobulin light chain are translocated to chromosome 8.

The MYC (8q24) break-apart probe is optimized to detect rearrangements involving the 8q24 locus in a triple-color, split assay on formalin fixed paraffin embedded tissue.

In addition Kretech has developed a probe for the specific use on cell material (KBI-10611).

Cat# KBI-10749 MYC (8q24) Break, TC (tissue)



MYC (8q24) Break, TC (tissue) probe hybridized to patient material showing a 8q24 distal break (1GB1R1GBR).
Image kindly provided by N. Van Acker, HistoGeneX, Antwerp.

Literature:

Fabris et al, 2003, Genes Chromosomes Cancer 37; 261-269.
Hummel et al, 2006, N Engl J Med 354; 2419-30.



| Ordering information | Gene Region | Tests | Cat# |
|----------------------------------|----------------|-------|-----------|
| ON MYC (8q24) Break, TC (tissue) | red/green/blue | 10 | KBI-10749 |

ON BCL2/IGH@ t(14;18) Fusion (tissue)

Follicular lymphoma is a mature B-Cell lymphoma, characterized by the presence of the t(14;18) translocation that juxtaposes the BCL2 locus on chromosome 18q21 to the immunoglobulin H (IGH) locus on chromosome 14q32, resulting in the overexpression of the antiapoptotic protein BCL2.

The BCL2/IGH@ t(14;18) Fusion Probe is optimized to detect the reciprocal translocation t(14;18) in a dual-color, dual-fusion assay on formalin fixed paraffin embedded tissue samples. In addition KREATECH has developed a probe for the specific use on cell material (KBI-10606).

Cat# KBI-10755 BCL2/IGH@ t(14;18) Fusion (tissue)



BCL2/IGH@ t(14;18) Fusion Probe hybridized to paraffin embedded lymph node material.

Image kindly provided by Philippa May, Imperial College, Hammersmith Hospital, London.

Literature:

Taniwaki M et al, 1995, Blood, 86; 1481-1486.
Poetsch M et al, 1996, J Clin Oncol, 14; 963-969.

Ordering information

| Gene Region | Tests | Cat# |
|---------------------------------------|-----------|--------------|
| ON BCL2/IGH@ t(14;18) Fusion (tissue) | red/green | 10 KBI-10755 |



MICRODELETION

MD GATA4 (8p23) / SE 8

The deletion of GATA4 (8p23) is found in patients with congenital heart disease. Besides the deletion of the region, duplications are found of the region flanked by low copy repeats 8p-OR-REPD (distal) and -REPP (proximal). These recurrent deletions are associated with a spectrum of anomalies, including congenital diaphragmatic hernia, developmental delay and neuropsychiatric findings. GATA4 is expressed in adult heart, epithelium and gonads. During fetal development, GATA4 is expressed in yolk sac endoderm and cells involved in heart formation.

The MD GATA4 (8p23) / SE 8 probe is optimized to detect deletions of the GATA4 gene region at 8p23 in a dual color assay on metaphase/interphase spreads, blood smears and bone marrow cells. The Chromosome 8 Satellite Enumeration (SE) probe is included to facilitate chromosome identification.

Cat# KBI-40118 GATA4 (8p23) / SE 8



GATA4 (8p23) / SE 8 Probe hybridized to patient material showing a deletion of the GATA4 (8p23) region (1R2G).

Image kindly provided by Dr. Marie-France Portnoi, Service de génétique et embryologie médicales Hôpital Armand-Trousseau, Paris.

Literature:

Bhatia et al, 1999, Prenat Diagn, 19, 863-867.
Giorda et al, 2007, Hum Mut, 28, 459-468.
Wat et al, 2009, Am J Med Genet Part A, 149A, 1661-1677.

Ordering information

| Gene Region | Tests | Cat# |
|------------------------|-----------|--------------|
| MD GATA4 (8p23) / SE 8 | red/green | 5 KBI-45118 |
| MD GATA4 (8p23) / SE 8 | red/green | 10 KBI-40118 |

ONCOLOGY – HEMATOLOGY DNA PROBES

From the 25,000 genes in the human genome, approximately 350 genes have been causally linked to the development of cancer. Variant or aberrant function of these so-called cancer genes may result from changes in genome copy number (through amplification, deletion, chromosome loss, or duplication), changes in gene and chromosome structure (through chromosomal translocation, inversion, or other rearrangements that lead to chimeric transcripts or deregulated gene expression) and point mutations (including base substitutions, deletions, or insertions in coding regions and splice sites).

The vast majority (90%) of cancer genes are mutated or altered through chromosomal aberrations in somatic tissue, about 10% are altered in the germ line, thereby transmitting heritable cancer susceptibility through successive generations. In addition to high resolution chromosome banding and advanced chromosomal imaging technologies, chromosome aberrations in cancer cells can be analyzed with an increasing number of large-scale, comprehensive genomic and molecular genetic technologies – including Fluorescence *In Situ* Hybridization (FISH).

The REPEAT-FREE™ POSEIDON™ Hematology DNA Probes are direct labeled, Ready-to-Use in hybridization buffer and available in a 10 test kit. The hematology probes are designed for use on interphase, metaphase chromosomes from cultured peripheral blood cells or cultured bone marrow samples. Most of the lymphoma probes are also optimized for use on FFPE tissue material and can be found at the product section solid tumor probes.

Table of Contents


| Description | Cat# | Page |
|---|-----------|-------|
| CML | | |
| ON FIP1L1-CHIC2-PDGFR (4q12) Del, Break | KBI-10003 | 14 |
| ON PDGFRB (5q33) Break | KBI-10004 | 15 |
| ON BCR/ABL t(9;22) Fusion | KBI-10005 | 11,35 |
| ON BCR/ABL t(9;22) Fusion | KBI-12005 | 11,35 |
| ON BCR/ABL t(9;22) TC, D-Fusion | KBI-10006 | 11,35 |
| ON FIP1L1-CHIC2-PDGFR (4q12) Del, Break, TC | KBI-10007 | 15 |
| ON BCR/ABL t(9;22,) DC, S-Fusion, ES | KBI-10008 | 12,35 |
| ON BCR/ABL t(9;22) DC, S-Fusion | KBI-10009 | 12,35 |
| ON p53 (17p13) / MPO (17q22) "ISO 17q" | KBI-10011 | 13 |
| ON JAK2 (9p24) Break | KBI-10012 | 13 |
| ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | KBI-10013 | II |
| CLL | | |
| ON DLEU (13q14) / 13qter | KBI-10102 | 17 |
| ON ATM (11q22) / SE 11 | KBI-10103 | 18 |
| ON GLI (12q13) / SE 12 | KBI-10104 | 20 |
| ON 6q21 / SE 6 | KBI-10105 | 18 |
| ON C-MYC (8q24) / SE 8 | KBI-10106 | 19 |
| ON ATM (11q22) / GLI (12q13) | KBI-10108 | 21 |
| ON 6q21 / MYC (8q24) | KBI-10117 | 21 |
| ON hTERC (3q26) / 3q11 | KBI-10110 | 19 |
| ON p53 (17p13) / SE 17 | KBI-10112 | 17 |
| ON p53 (17p13) / SE 17 | KBI-12112 | 17 |
| ON DLEU (13q14) / p53 (17p13) | KBI-10113 | 20 |
| ON p53 (17p13) / ATM (11q22) | KBI-10114 | 22 |
| MDS | | |
| ON MDS 7q- (7q22; 7q36) | KBI-10202 | 25 |
| ON MDS 20q- (PTPRT 20q12) / 20q11 | KBI-10203 | 26 |
| ON EVI t(3;3); inv(3) (3q26) Break | KBI-10204 | 26 |
| ON EVI t(3;3); inv(3) (3q26) Break, TC | KBI-10205 | 27 |
| ON MDS 7q- (7q22; 7q36) / SE 7 TC | KBI-10207 | 25 |
| ON hTERT (5p15) / 5q31 | KBI-10208 | 23 |
| ON MDS 5q- (5q31; 5q33) | KBI-10209 | 24 |
| ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC | KBI-10210 | 24 |
| AML | | |
| ON AML/ETO t(8;21) Fusion | KBI-10301 | 28 |
| ON PML/RARA t(15;17) Fusion | KBI-10302 | 29 |
| ON PML/RARA t(15;17) Fusion | KBI-12302 | 29 |

| Description | Cat# | Page |
|---------------------------------|-----------|----------|
| ON MLL (11q23) Break | KBI-10303 | 29,35 |
| ON CBFb t(16;16), inv(16) Break | KBI-10304 | 32 |
| ON RARA (17q21) Break | KBI-10305 | 31 |
| ON DEK / NUP214 t(6;9) Fusion | KBI-10306 | 32 |
| ON MLL/MLLT1 t(11;19) Fusion | KBI-10307 | 30 |
| ON MLL/MLLT3 t(9;11) Fusion | KBI-10308 | 30 |
| ON MLL/MLLT4 t(6;11) Fusion | KBI-10309 | 31 |
| ON MECOM / RUNX1 t(3;21) Fusion | KBI-10310 | II |
| ON NUP98 (11p15) Break | KBI-10311 | III |
| ALL | | |
| ON TEL/AML t(12;21) Fusion | KBI-10401 | 33 |
| ON p16 (9p21) / 9q21 | KBI-10402 | 34 |
| ON ETV6 (TEL) (12p13) Break | KBI-10403 | 34 |
| ON MLL/AFF1 t(4;11) Fusion | KBI-10404 | III |
| Multiple Myeloma | | |
| ON MM 11q23 / DLEU (13q14) | KBI-10502 | 38 |
| ON MM 1q21 / 8p21 | KBI-10503 | 39 |
| ON MM 15q22 / 6q21 | KBI-10504 | 40 |
| ON MM 1q21 / SRD (1p36) | KBI-10507 | 41 |
| ON MM 15q22 / 9q34 | KBI-10508 | 40 |
| ON MM 19q13 / p53 (17p13) | KBI-10509 | 39 |
| ON MAFB/IGH@ t(14;20) Fusion | KBI-10510 | IV |
| Lymphoma related probes | | |
| ON IGH (14q32) Break | KBI-10601 | 22,41,44 |
| ON FGFR3/IGH t(4;14) Fusion | KBI-10602 | 38,47 |
| ON MYC/IGH t(8;14) Fusion | KBI-10603 | 42 |
| ON BCL1/IGH t(11;14) Fusion | KBI-10604 | 43 |
| ON MYEOV/IGH t(11;14) Fusion | KBI-10605 | 36,47 |
| ON BCL2/IGH t(14;18) Fusion | KBI-10606 | 43 |
| ON BCL6 (3q27) Break | KBI-10607 | 45 |
| ON MALT (18q21) Break | KBI-10608 | 45 |
| ON CCND1 (BCL1;11q13) Break | KBI-10609 | 46 |
| ON MAF/IGH t(14;16) Fusion | KBI-10610 | 37 |
| ON MYC (8q24) Break, TC | KBI-10611 | 47 |
| ON BCL2 (18q21) Break | KBI-10612 | 44 |
| ON MAFB/IGH@ t(14;20) Fusion | KBI-10510 | IV |

Translocation, Dual-Fusion Assay

Dual-fusion, dual-color FISH assays for translocation utilize large probes that span 2 breakpoints or flanking regions on the different chromosomes. Dual-fusion, dual-color FISH is optimal for detection of low levels of nuclei possessing a simple balanced translocation, as it greatly reduces the number of normal background nuclei with an abnormal signal pattern.

Translocation, Dual-Fusion Assay




Expected signal pattern:

In normal intact cells, two separate red and two separate green individual signals will be visible, whereas a reciprocal translocation will generate two fused red/green signals (often appearing as single yellow signals), accompanied by one red and one green signal (representing the normal chromosomes).

Translocation, Break-Apart or Split Assay

FISH using dual-color, break-apart probes is very useful in the evaluation of genes known to have multiple translocation partners; the differently colored probes hybridize to targets on opposite sides of the breakpoint of the affected gene.

Translocation, Break or Split Assay



Expected signal pattern:

In normal cells two sets of red/green-fused signals (representing the two alleles) will be visible. In an abnormal diploid cell, in which one allele has been split by a translocation, a separated red and green signal will be visible in addition to the normal fused signal.

Chronic Myeloproliferative Disorders (CMPD)


Chromosomal translocations in chronic myeloproliferative diseases (CMPD) almost invariably correlated with expression of constitutively activated fusion tyrosine kinases. The hallmark of these diseases is CML, where the BCR/ABL activated tyrosine kinase results from the balanced reciprocal Philadelphia chromosome translocation t(9;22).

Chronic Myelogenous Leukemia (CML) - BCR/ABL t(9;22)


CML is a malignant chronic myeloproliferative disorder (MPD) of the hematopoietic stem cell. All CML have a t(9;22) resulting from fusion of the 3' ABL region at 9q34 with the 5' BCR region at 22q11. This chimeric BCR/ABL gene encodes a constitutively activated protein tyrosine kinase with profound effects on cell cycle, adhesion, and apoptosis. Understanding this process has led to the development of the drug imatinib mesylate (Gleevec™), the first in a new class of genetically targeted agents, this is a major advance in cancer treatment. Several different approaches are used to analyze the BCR/ABL t(9;22)(q34;q11) by FISH each providing different details about this translocation.

Interpretation guidelines for POSEIDON™ BCR/ABL Probes


BCR/ABL t(9;22) Dual-color, Dual-Fusion Probe




BCR/ABL t(9;22) Triple-color, Dual-Fusion Probe




BCR/ABL t(9;22) Dual-color, Single-Fusion Extra Signal Probe



BCR/ABL t(9;22) Dual-color, Single-Fusion Probe



Mm-BCR/ABL t(9;22) Dual-Color, Single-Fusion, Extra Signal




BCR/ABL Product Family

The Philadelphia chromosome is an abnormally short chromosome 22 that is one of the two chromosomes involved in a translocation with chromosome 9. This translocation t(9;22)(q34;q11) takes place in a single bone marrow cell and, through the process of clonal expansion, gives rise to the leukemia.

ABL and BCR are genes on chromosomes 9 and 22, respectively. The ABL gene encodes a tyrosine kinase enzyme whose activity is tightly controlled. In the formation of the Ph translocation, two fusion genes are generated: BCR-ABL on the Ph chromosome and ABL-BCR on the chromosome 9 participating in the translocation. The BCR-ABL gene encodes a protein with deregulated tyrosine kinase activity.

The presence of this protein in the CML cells is strong evidence of its pathogenetic role. The efficacy in CML of a drug that inhibits the BCR-ABL tyrosine kinase has provided the final proof that the BCR-ABL oncoprotein is the unique cause of CML. The POSEIDON™ portfolio contains now 4 different probes for BCR/ABL to suit all needs for the detection of t(9;22) by FISH:



ON BCR/ABL t(9;22), Fusion


The BCR/ABL t(9;22) Fusion is optimized to detect the t(9;22) (q34;q11) reciprocal translocation in a dual-color, dual-fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

This probe will also detect cryptic insertions of ABL into the BCR region not detectable by karyotyping and therefore described as Ph-negative.

ON BCR/ABL t(9;22), TC, D-Fusion



The BCR/ABL t(9;22), TC, D-Fusion probe is a triple-color, dual-fusion probe build from the same regions as the dual-color, dual-fusion probe, but the proximal BCR region is labeled in blue. Using the triple-color probe allows to distinguish between the derivative chromosome 22, the Philadelphia chromosome, which will be observed as purple (red/blue) color, while the derivative chromosome 9 will show a yellow (red/green) signal.

Cat.# KBI-10005 BCR/ABL t(9;22), Fusion



BCR/ABL t(9;22) Fusion probe hybridized on patient material showing t(9;22) (q34;q11) reciprocal translocation (2RG1R1G).
Image kindly provided by Monika Conchon, São Paulo.

Cat.# KBI-10006 BCR/ABL t(9;22), TC, D-Fusion



BCR/ABL t(9;22), TC, D-Fusion probe hybridized on patient material showing translocation of distal BCR (1BG1RB1R1G).
Image kindly provided by Prof. Siebert, Kiel.

Literature:

Morris et al, 1990, Blood, 76: 1812-1818.
Dewald et al, 1998, Blood, 91: 3357-3365.
Kolomietz et al, 2001, Blood, 97: 3581-3588.
Huntly et al, 2003, Blood, 102: 1160-1168.
Tkachuk et al., 1990, Science 250, 559-562.

Literature:

Morris et al, 1990, Blood, 76: 1812-1818.
Dewald et al, 1998, Blood, 91: 3357-3365.
Kolomietz et al, 2001, Blood, 97: 3581-3588.
Huntly et al, 2003, Blood, 102: 1160-1168.
Tkachuk et al., 1990, Science 250, 559-562.

| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| ON BCR/ABL t(9;22) Fusion | red/green | 10 | KBI-10005 |
| ON BCR/ABL t(9;22) Fusion | red/green | 20 | KBI-12005 |

| Ordering information | Color | Tests | Cat# |
|---------------------------------|----------------|-------|-----------|
| ON BCR/ABL t(9;22) TC, D-Fusion | red/green/blue | 10 | KBI-10006 |



ON BCR/ABL t(9;22), DC, S-Fusion, ES

A single-fusion assay is preferably used for the initial screening of CML and ALL patients. Proximal to the breakpoints on chromosome 9q34, this probe will provide an extra signal on the der(9q34) in case of a t(9;22). The Philadelphia chromosome, der(22q), is visualized by the fusion signal.

ON BCR/ABL t(9;22), DC, S-Fusion



A simple dual-color, single-fusion assay is preferably used for the initial screening of CML and ALL patients. The Philadelphia chromosome, der(22q), is visualized by a fusion signal while the der(9q) shows no signal.

Cat.# KBI-10008 BCR/ABL t(9;22), Dual-Color, Single-Fusion, Extra Signal



BCR/ABL t(9;22), DC, S-Fusion, ES probe hybridized to patient material showing t(9;22)translocation (1RG1r1R1G).
Material kindly provided by Dr. Balogh, Budapest.

Cat.# KBI-10009 BCR/ABL t(9;22) Dual-Color, Single-Fusion



BCR/ABL t(9;22), DC, S-Fusion probe hybridized to patient material showing t(9;22)translocation (1RG1R1G).
Material kindly provided by Dr. Balogh, Budapest.

Literature:

Morris et al, 1990, Blood, 76: 1812-1818.
Dewald et al, 1998, Blood, 91: 3357-3365.
Kolomietz et al, 2001, Blood, 97: 3581-3588.
Huntly et al, 2003, Blood, 102: 1160-1168.
Tkachuk et al., 1990, Science 250, 559-562.

Literature:

Morris et al, 1990, Blood, 76: 1812-1818.
Dewald et al, 1998, Blood, 91: 3357-3365.
Kolomietz et al, 2001, Blood, 97: 3581-3588.
Huntly et al, 2003, Blood, 102: 1160-1168.
Tkachuk et al., 1990, Science 250, 559-562.

| Ordering information | Color | Tests | Cat# |
|-------------------------------------|-----------|-------|-----------|
| ON BCR/ABL t(9;22) DC, S-Fusion, ES | red/green | 10 | KBI-10008 |

| Ordering information | Color | Tests | Cat# |
|---------------------------------|-----------|-------|-----------|
| ON BCR/ABL t(9;22) DC, S-Fusion | red/green | 10 | KBI-10009 |



CML secondary chromosomal changes

ON p53 (17p13) / MPO (17q22) "ISO 17q"

Isochromosome 17q is the most common isochromosome in cancer. It plays an important role in tumor development and progression. Hematologic malignancies such as chronic myeloid leukemia (CML) with isochromosome 17q carry a poor prognosis. Isochromosome 17q is the most common chromosome abnormality in primitive neuroectodermal tumors and medulloblastoma. Isochromosome 17q is, by convention, symbolized as i(17q).

The p53 (17p13) / MPO (17q22) "ISO 17q" probe is optimized to detect copy numbers of the p53 gene region at 17p13 and MPO gene region at 17q22. In case of i(17q) a signal pattern of three red signals for MPO (17q22) and one signal for p53 at 17p13 is expected.

Cat.# KBI-10011 p53 (17p13) / MPO (17q22) "ISO 17q"



p53 (17p13) / MPO (17q22) "ISO 17q" probe hybridized to peripheral blood of a CLL patient with an isochromosome 17 (3R1G).
Image kindly provided by Dr. Lana Harder, Kiel.

Literature:

Becher et al, 1990, Blood, 75: 1679-1683.
Fioretos et al, 1999, Blood, 94: 225-232.



| Ordering information | Color | Tests | Cat# |
|--|-----------|-------|-----------|
| ON p53 (17p13) / MPO (17q22) "ISO 17q" | green/red | 10 | KBI-10011 |

ON JAK2 (9p24) Break

Janus Kinase 2 (JAK2) is a tyrosine kinase involved in cytokine signaling. Mutations and translocations involving the JAK2 gene region are observed in myeloproliferative neoplasms. The common JAK2617V>F point mutation and translocations results in constitutive activation of JAK2. Translocations are described with the following fusion partners: PCM1, BCR, ETV6 (TEL), SSBP2 and 3q21. Patients with the JAK2617V>F point mutation can also exhibit a numerical gain of the gene.

The JAK2 (9p24) Break probe is optimized to detect translocations involving the JAK2 gene region at region 9p24 in a dual-color, split assay on metaphase/interphase spreads. The JAK2 (9p24) Break probe can not be used to detect point mutations, and it has not been optimized to detect gene amplifications.

Cat# KBI-10012 JAK2 (9p24) Break



JAK2 (9p24) Break probe hybridized to bone marrow sample (2RG).

Literature:

Najfeld V et al, 2007, Exp Hematol, 35, 1668-1676.
Smith C et al, 2008, Hum Pathol, 39, 795-810.
Poitras J et al, 2008, Genes Chromosomes Cancer, 47, 884-889.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON JAK2 (9p24) Break | red/green | 10 | KBI-10012 |

SE 8 (D8Z1)


SE 7 (D7Z1) / SE 8 (D8Z1)

Gain of chromosome 8 is the most common secondary chromosomal aberration in CML (approx. 34%).

Cat.# KBI-20008 SE 8 (D8Z1)

Cat.# KBI-20031 SE 7 (D7Z1) / SE 8 (D8Z1)

See description under Satellite Enumeration probes on page 94.





SE 7 / SE 8 showing trisomy 8.
Material kindly provided by Dr. Balogh, Budapest.

Other Myeloproliferative Diseases:

ON FIP1L1-CHIC2-PDGFRα (4q12) Del, Break

Idiopathic hypereosinophilic syndrome (HES) and chronic eosinophilia leukemia (CEL) represent the most recent additions to the list of molecularly defined chronic myeloproliferative disorders. A novel tyrosine kinase that is generated from fusion of the Fip1-like 1 (FIP1L1) and PDGFRα (PDGFRA) genes has been identified as a therapeutic target for imatinib mesylate in hypereosinophilic syndrome (HES). This fusion results from an approximately 800-kb interstitial chromosomal deletion that includes the cysteine-rich hydrophobic domain 2 (CHIC2) locus. The FIP1L1-CHIC2-PDGFRα probe is optimized to detect the CHIC2 deletion at 4q12 associated with the FIP1L1/PDGFRα fusion in a Dual-Color, split assay. It also allows the detection of translocation involving the FIP1L1 and PDGFRA region. However, chromosome 4 polyploidy may provide additional signals not associated with a translocation involving 4q12.

Cat.# KBI-10003 FIP1L1-CHIC2-PDGFRα (4q12) Del, Break



FIP1L1-CHIC2-PDGFRα (4q12) Del, Break probe hybridized to a normal interphase/metaphase (2RG).

Literature:

Cools et al, N Engl J Med, 2003, 348, 1201-1214.
Godlib et al, Blood, 2004, 103, 2879-2891.

| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| SE 8 (D8Z1) | red/green | 10 | KBI-20008 |
| SE 7 (D7Z1) / SE 8 (D8Z1) | red/green | 10 | KBI-20031 |

| Ordering information | Color | Tests | Cat# |
|--|-----------|-------|-----------|
| ON FIP1L1-CHIC2-PDGFRα (4q12) Del, Break | red/green | 10 | KBI-10003 |

ON FIP1L1-CHIC2-PDGFR A (4q12) Del, Break, TC



The FIP1L1-CHIC2-PDGFR A triple-color probe is optimized to detect the CHIC2 deletion at 4q12 associated with the FIP1L1/PDGFR A fusion in a triple-color, split assay. It also allows the detection of translocation involving the FIP1L1 and PDGFR A region.

ON PDGFR B (5q33) Break

PDGFR B activation has been observed in patients with chronic myelomonocytic leukemia/atypical chronic myeloid leukemia and has been associated with 11 translocation partners, the best known is the ETV6 gene on 12p13, causing a t(5;12) translocation. Cytogenetic responses are achieved with imatinib in patients with PDGFR B fusion positive, BCR/ABL negative CMPDs.



The PDGFR B probe is optimized to detect translocations involving the PDGFR B region at 5q33 in a dual-color, split assay.

Cat.# KBI-10007 FIP1L1-CHIC2-PDGFR A (4q12) Del, Break, Triple-Color



FIP1L1-CHIC2-PDGFR A (4q12) Del, Break, TC probe hybridized to a normal metaphase (2BRG).

Cat.# KBI-10004 PDGFR B (5q33) Break



PDGFR B (5q33) Break probe hybridized to a normal metaphase (2RG).

Literature:

Cools et al, 2003, N Engl J Med, 348: 1201-1214.
Griffin et al, 2003, PNAS, 100: 7830-7835.
Gotlib et al, 2004, Blood, 103; 2879-2891.

Literature:

Wlodarska et al, 1997, Blood, 89: 1716-1722.
Wilkinson et al, 2003, Blood, 102: 4287-419.



| Ordering information | Color | Tests | Cat# |
|---|----------------|-------|-----------|
| ON FIP1L1-CHIC2-PDGFR A (4q12) Del, Break, TC | red/green/blue | 10 | KBI-10007 |

| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON PDGFR B (5q33) Break | red/green | 10 | KBI-10004 |

ON FGFR1 (8p12) Break

FGFR1 has been implicated in the tumorigenesis of haematological malignancies, where it is frequently involved in balanced chromosomal translocations, including cases of chronic myeloid leukaemia (BCR-FGFR1 fusion) and the 8p11 myeloproliferative syndrome/stem cell leukaemia-lymphoma syndrome, which is characterized by myeloid hyperplasia and non-Hodgkin's lymphoma with chromosomal translocations fusing several genes, the most common being a fusion between ZNF198 and FGFR1.

Cat.# KBI-10737 FGFR1 (8p12) Break



FGFR1 (8p12) Break probe hybridized to patient material showing a break at 8p12 (1RG1R1G).

Literature:

Smedley et al, 1998, Hum Mol Genet. 7; 627-642.
Sohal et al, 2001, Genes Chrom. Cancer 32; 155-163.
Kwak et al, J Clin Oncol., 27(26):4247-53.

Chronic Lymphocytic Leukemia (CLL)

CLL accounts for about 30% of all leukemias in Europe and the USA. Distinct clonal chromosomal abnormalities can be identified in up to 90% of CLL cases of the B-cell lineage. By FISH the most common chromosomal changes in CLL and their frequencies have been identified as shown in the table below.

| | |
|-------------|-----|
| Del(13q14) | 55% |
| Del(11q) | 18% |
| Trisomy 12q | 16% |
| Del(17p) | 7% |
| Del(6q) | 6% |
| Trisomy 8q | 5% |
| t(14q32) | 4% |
| Trisomy 3q | 3% |



| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON FGFR1 (8p12) Break | red/green | 10 | KBI-10737 |

ON DLEU (13q14) / 13qter

Deletions of chromosome 13q14 have been reported not only in CLL but in a variety of human tumors, including other types of lymphoid and myeloid tumors, as well as prostate, head and neck, and non-small cell lung cancers. The deletion of 13q may be limited to a single locus (13q14), or accompanied with the loss of a larger interstitial region of the long arm of chromosome 13. A minimal critical region of 400 kb has been described containing the DLEU1, DLEU2 and RFP2 genes.

The DLEU (13q14) specific DNA probe is optimized to detect copy numbers of the DLEU gene region at 13q14. The 13qter (13q34) region is included to facilitate chromosome identification.

Cat.# KBI-10102 DLEU (13q14) / 13qter



DLEU (13q14) / 13qter probe hybridized to patient material showing a 13q14 deletion (1R2G).

Image kindly provided by Dr. Dastugue, Toulouse.

Literature:

Wolf et al, 2001, Hum Mol Genet, 10: 1275-1285.
Corcoran et al, 1998, Blood, 91: 1382-1390.



| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON DLEU (13q14) / 13qter | red/green | 10 | KBI-10102 |

ON p53 (17p13) / SE 17

The p53 tumor suppressor gene at 17p13, has been shown to be implicated in the control of normal cellular proliferation, differentiation, and apoptosis. Allelic loss, usually by deletion, and inactivation of p53 have been reported in numerous tumor types and are associated with poor prognosis in CLL.

The p53 (17p13) specific DNA probe is optimized to detect copy numbers of the p53 gene region at 17p13. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

Cat.# KBI-10112 p53 (17p13) / SE 17



p53 (17p13) / SE 17 probe hybridized to patient material showing a 17p13 deletion at the p53 gene region (1R2G).

Literature:

Amiel A et al, 1997, Cancer Gener.Cytogenet., 97; 97-100.
Drach J et al, 1998, Blood, 92; 802-809.



| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON p53 (17p13) / SE 17 | red/green | 10 | KBI-10112 |
| ON p53 (17p13) / SE 17 | red/green | 20 | KBI-12112 |

ON ATM (11q22) / SE 11

Chromosome 11q22.3-23.1 deletions involving the ataxia-telangiectasia mutated (ATM) locus are detected at diagnosis in 15-20% of cases of B-cell chronic lymphocytic leukemia (CLL) and are associated with a relatively aggressive disease. Loss of the 11q22-23 region, where the ataxia-telangiectasia mutated (ATM) gene is located, is also one of the most frequent secondary chromosomal aberrations in mantle cell lymphoma.

The ATM (11q22.3) specific DNA probe is optimized to detect copy numbers of the ATM gene region at region 11q22.3. The chromosome 11 satellite enumeration (SE 11) at D11Z1 probe is included to facilitate chromosome identification.

Cat.# KBI-10103 ATM (11q22) / SE 11



ATM (11q22) / SE 11 hybridized to patient material showing a 11q22 deletion at the ATM gene region (1R2G).
Image kindly provided by Dr. Wenzel, Basel.

Literature:

Doehner et al, 1997, Blood, 89: 2516-2522.
Bigoni et al, 1997, Leukemia, 11: 1933-1940.



| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON ATM (11q22) / SE 11 | red/green | 10 | KBI-10103 |

ON 6q21 / SE 6

Deletions affecting the long arm of chromosome 6 (6q) are among the most commonly observed chromosomal aberrations in lymphoid malignancies and have been identified as an adverse prognostic factor in subsets of tumors including CLL. A minimal deletion region has been identified within a 2-megabase segment of 6q21, between D6S447 and D6S246. The SEC63 gene is located within this critical region.

The 6q21 specific DNA probe is optimized to detect copy numbers of 6q at region 6q21. The chromosome 6 satellite enumeration probe (SE 6) at D6Z1 is included to facilitate chromosome identification.

Cat.# KBI-10105 6q21 / SE 6



6q21 / SE 6 probe hybridized to a normal metaphase (2R2G).

Literature:

Sherratt et al, 1997, Chromosome Res, 5: 118-124.
Zhang et al, 2000, Genes Chrom Cancer, 27: 52-58.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON 6q21 / SE 6 | red/green | 10 | KBI-10105 |

ON C-MYC (8q24) / SE 8

The C-MYC gene produces an oncogenic transcription factor that affects diverse cellular processes involved in cell growth, cell proliferation, apoptosis and cellular metabolism. The C-MYC oncogene has been shown to be amplified in many types of human cancer such as bladder, breast and cervical. Amplification at 8q24 including C-MYC is also observed in 5% of CLL patients. C-MYC is also the prototype for oncogene activation by chromosomal translocation.

The C-MYC (8q24) specific DNA probe is optimized to detect copy numbers of the C-MYC gene region at 8q24. The chromosome 8 satellite enumeration probe (SE 8) at D8Z1 is included to facilitate chromosome identification.

Cat.# KBI-10106 C-MYC (8q24) / SE 8



C-MYC (8q24) / SE 8 hybridized to a normal metaphase (2R2G).

Literature:

Greil et al, 1991, Blood, 78: 180-191.

Note:

This probe should not be used to detect translocations involving C-MYC.



| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON C-MYC (8q24) / SE 8 | red/green | 10 | KBI-10106 |

ON hTERC (3q26) / 3q11

Amplification of the 3q26-q27 has a high prevalence in cervical, prostate, lung, and squamous cell carcinoma. This amplification can also be found to a lesser extent in CLL patients. The minimal region of amplification was refined to a 1- to 2-Mb genomic segment containing several potential cancer genes including hTERC, the human telomerase RNA gene.

The hTERC (3q26) specific DNA probe is optimized to detect copy numbers of the hTERC gene region at region 3q26. The 3q11 region probe is included to facilitate chromosome identification.

Cat.# KBI-10110 hTERC (3q26) / 3q11



hTERC (3q26) / 3q11 probe hybridized to a normal interphase/metaphase (2R2G).

Literature:

Arnold et al, 1996, Genes Chrom Cancer, 16: 46-54.

Soder et al, 1997, Oncogene, 14: 1013-1021.



| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON hTERC (3q26) / 3q11 | red/green | 10 | KBI-10110 |

ON GLI (12q13) / SE 12

Trisomy 12 is the most common numerical chromosomal aberration in patients with B-cell chronic lymphocytic leukemia (B-CLL). Partial trisomy 12 of the long arm of chromosome 12 consistently includes a smaller region at 12q13-15 and has been observed in CLL and several other tumors. A number of loci located close to either MDM2 or CDK4/SAS, including the genes GADD153, **GLI**, RAP1B, A2MR, and IFNG, were found to be coamplified.

The GLI (12q13) specific DNA probe is optimized to detect copy numbers of the GLI gene region at region 12q13. The chromosome 12 satellite enumeration probe (SE 12) D12Z3 is included to facilitate chromosome identification.

Cat.# KBI-10104 GLI (12q13) / SE 12



GLI (12q13) / SE 12 hybridized to patient material showing GLI (12q13) amplification (3R2G).
Image kindly provided by Dr. Wenzel, Basel.

Literature:

Merup et al, 1997, Eur J Haematol, 58: 174-180.
Dierlamm et al., 1997, Genes Chrom Cancer, 20: 155-166.

| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON GLI (12q13) / SE 12 | red/green | 10 | KBI-10104 |



CLL probe combinations

Most of the CLL probes are also available as combinations covering two critical loci in one hybridization. This is of particular interest if patient material is limited. The disadvantage is that no internal control is added. The following CLL probe combinations are available:

ON DLEU (13q14) / p53 (17p13)

Deletion of DLEU at 13q14 indicates a rather good prognosis, deletion of p53 at 17p13 is associated with poor prognosis.

Cat.# KBI-10113 DLEU (13q14) / p53 (17p13)



DLEU (13q14) / p53 (17p13) hybridized to a normal metaphase (2R2G).

Literature:

Amiel A et al, 1997, Cancer Genet.Cytogenet, 97: 97-100.
Drach J et al, 1998, Blood, 92: 802-809.
Stilgenbauer S et al, 1998, Oncogene, 16: 1891 – 1897.
Wolf S et al, 2001, Hum. Molec. Genet., 10: 1275-1285.

| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON DLEU (13q14) / p53 (17p13) | red/green | 10 | KBI-10113 |



ON ATM (11q22) / GLI (12q13)

Deletion of ATM at 11q22-q23 indicates a rather poor prognosis, amplification of GLI at 12q13 is associated with an intermediate prognosis.

ON 6q21 / MYC (8q24)



Deletion of 6q21 indicates an intermediate risk category, amplification of MYC at 8q24 is associated with poor prognosis.

Cat.# KBI-10108 ATM (11q22) / GLI (12q13)



ATM (11q22) / GLI (12q13) hybridized to a normal metaphase (2R2G).

Cat.# KBI-10117 6q21 / MYC (8q24)



6q21 / MYC (8q24) hybridized to a normal interphase/metaphase (2R2G).

Literature:

Döhner H et al, 1997, Blood, 7; 2516-2522.
Boultonwood J, 2001, J. Clin. Pathol., 54; 512-516.
Dierlamm J et al, 1998, Genes Chromosomes Cancer, 20; 155-166.
Döhner H et al, 1999, J. Molec. Med., 77; 266-281.

| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON ATM (11q22) / GLI (12q13) | red/green | 10 | KBI-10108 |

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON 6q21 / MYC (8q24) | red/green | 10 | KBI-10117 |

ON p53 (17p13) / ATM (11q22)

Deletion of p53 and ATM are both indicating poor prognosis in CLL.

Other relevant CLL probes:

ON IGH (14q32) Break

Cat.# KBI-10601 IGH (14q32) Break

See description under Lymphoma on page 44.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON IGH (14q32) Break | red/green | 10 | KBI-10601 |



SE 12 (D12Z3)

Cat.# KBI-20012 SE 12 (D12Z3)

See description under Satellite Enumeration probes on page 94.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| SE 12 (D12Z3) | red/green | 10 | KBI-20012 |

Cat.# KBI-10114 p53 (17p13) / ATM (11q22)



p53 (17p13) / ATM (11q22) hybridized to a normal metaphase (2R2G).

| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON p53 (17p13) / ATM (11q22) | red/green | 10 | KBI-10114 |



Myelodysplastic Syndromes (MDS)

The myelodysplastic syndromes (MDS) are a heterogeneous group of hematopoietic disorders characterized in most patients by peripheral blood cytopenia with hypercellular bone marrow and dysplasia of the cellular elements. Cytogenetic studies play a major role in confirmation of diagnosis and prediction of clinical outcome in MDS, and have contributed to the understanding of its pathogenesis. Clonal chromosomal abnormalities are detected by routine karyotyping techniques in 40%–70% of cases of de novo MDS, and 95% of cases of therapy-related MDS.

ON hTERT (5p15) / 5q31

The hTERT / 5q31 dual-color probe can be used to detect deletions involving band 5q31 in MDS and AML. The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31. The hTERT gene region at 5p15 is included to facilitate chromosome identification.

Cat.# KBI-10208 hTERT (5p15) / 5q31



hTERT (5p15) / 5q31 probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:
Zhao et al, 1997, PNAS, 94; 6948-6053.
Horrigan et al, 2000, Blood, 95; 2372-2377.


| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON hTERT (5p15) / 5q31 | red/green | 10 | KBI-10208 |

ON MDS 5q- (5q31; 5q33)

The presence of del(5q), either as the sole karyotypic abnormality or as part of a more complex karyotype, has distinct clinical implications for myelodysplastic syndromes (MDS) and acute myeloid leukemia. Interstitial 5q deletions are the most frequent chromosomal abnormalities in MDS and are present in 10% to 15% of MDS patients. Two different critical regions are described, one at 5q31-q33 containing the CSF1R and RPS14 gene regions, characteristic for the '5q-' syndrome, and a more proximal located region at 5q13-q31 containing the CDC25C and EGR1 gene regions.

The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31 and the CSF1R/RPS14 gene region at 5q33 simultaneously in a dual-color assay.

Cat.# KBI-10209 MDS 5q- (5q31; 5q33)



MDS 5q- (5q31; 5q33) probe hybridized to patient material showing a 5q33 deletion (1R2G).
Image kindly provided by Dr. Mohr, Dresden.

Literature:


Boultonwood J e.a., Blood 2002; 99: 4638-4641.
Zhao N e.a., PNAS 1997; 94: 6948-6953.
Wang e.a., Haematologica 2008; 93: 994-1000 .
Ebert BL e.a., Nature 2008; 451: 335-339.
Mohamedali A and Mufti GJ, Brit J Haematol 2008; 144: 157-168.

| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON MDS 5q- (5q31; 5q33) | red/green | 10 | KBI-10209 |

ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC

The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31 and the CSF1R/RPS14 gene region at 5q33 simultaneously in a dual-color assay. The triple-color probe provides in addition to the two critical regions a control in blue targeting the hTERT gene region at 5p15.

Cat.# KBI-10210 MDS 5q- (5q31; 5q33) / hTERT (5p15), Triple-Color



MDS 5q- (5q31; 5q33) / hTERT (5p15) probe hybridized to a normal metaphase (2R2G2B).

Literature:

Boultonwood J e.a., Blood 2002; 99: 4638-4641.
Zhao N e.a., PNAS 1997; 94: 6948-6953.
Wang e.a., Haematologica 2008; 93: 994-1000.
Ebert BL e.a., Nature 2008; 451: 335-339.
Mohamedali A and Mufti GJ, Brit J Haematol 2008; 144: 157-168.



| Ordering information | Color | Tests | Cat# |
|---|----------------|-------|-----------|
| ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC | red/green/blue | 10 | KBI-10210 |

ON MDS 7q- (7q22; 7q36)

Loss of a whole chromosome 7 or a deletion of the long arm, del(7q), are recurring abnormalities in malignant myeloid diseases. Most deletions are interstitial and there are two distinct deleted segments of 7q. The majority of patients have proximal breakpoints in bands q11-22 and distal breakpoints in q31-36. The CCAAT displacement protein CUX1 gene region is located in the 7q22 critical region.

The 7q- specific DNA probe is optimized to detect copy number of 7q at 7q22 and at 7q36 simultaneously in a dual-color assay.

Cat.# KBI-10202 MDS 7q- (7q22; 7q36)



MDS 7q- (7q22; 7q36) hybridized to patient material showing a 7q36 deletion (1R2G).

Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

LeBeau et al., 1996, Blood, 88: 1930-1935.
Doehner et al, 1998, Blood, 92: 4031-4035.



| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON MDS 7q- (7q22; 7q36) | red/green | 10 | KBI-10202 |

ON MDS 7q- (7q22; 7q36) / SE 7 TC

The 7q- specific DNA probe is optimized to detect copy number of 7q at 7q22 and at 7q36 simultaneously in a dual-color assay.

The chromosome 7 satellite enumeration probe (SE 7) at D7Z1 in blue is included to facilitate chromosome identification.

Cat.# KBI-10207 MDS 7q (7q22; 7q36) / SE 7, Triple-Color



MDS 7q (7q22; 7q36) / SE 7 TC probe hybridized to a normal metaphase (2R2G2B).

Literature:

LeBeau et al., 1996, Blood, 88: 1930-1935.
Doehner et al, 1998, Blood, 92: 4031-4035.



| Ordering information | Color | Tests | Cat# |
|----------------------------------|----------------|-------|-----------|
| ON MDS 7q- (7q22; 7q36) / SE7 TC | red/green/blue | 10 | KBI-10207 |

ON MDS 20q- (PTPRT 20q12) / 20q11

Acquired deletions of the long arm of chromosome 20 are found in several hematologic conditions and particularly in the myeloproliferative disorders (MPD) and myelodysplastic syndromes and acute myeloid leukemia (MDS/AML). A minimal critical region deleted in MPD and MDS has been identified at 20q12 which includes a protein tyrosine phosphatase receptor gene (PTPRT).

The 20q- (PTPRT, 20q12) specific DNA probe is optimized to detect copy numbers of 20q at region 20q12. A 20q11 region specific probe is included to facilitate chromosome identification.

Cat.# KBI-10203 MDS 20q- (PTPRT 20q12) / 20q11



MDS 20q- (PTPRT 20q12) / 20q11 probe hybridized to patient material showing 20q- deletion (1R2G).
Material kindly provided by Labdia Labordiagnostik, Vienna.

Literature:

Bench et al, 2000, Oncogene, 19: 3902-3913.
Asimakopoulos et al, 1994, Blood, 84: 3086-3094.



| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON 20q- (PTPRT 20q12) / 20q11 | red/green | 10 | KBI-10203 |

ON EVI t(3;3); inv(3) (3q26) Break

The inv(3)(q21;q26) is a recurrent cytogenetic aberration of myeloid malignancy associated with fusion of EVI1 and RPN1 and a poor disease prognosis. Genomic breakpoints in 3q26 are usually located proximal to the EVI1 locus, spanning a region of several hundred kilobases. Other recurrent and sporadic rearrangements of 3q26 also cause transcriptional activation of EVI1 including the translocations t(3;3)(q21;q26) and t(3;21)(q26;q22). Breakpoints in the latter rearrangements span a wider genomic region of over 1 megabase encompassing sequences distal to EVI1 and neighboring gene MDS1.

The EVI t(3;3) inv(3) Break, dual-color probe is optimized to detect the inversion of chromosome 3 involving the EVI1 gene region at 3q26 in a dual-color, split assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat.# KBI-10204 EVI t(3;3); inv(3) (3q26) Break



EVI t(3;3); inv(3) (3q26) Break probe hybridized to patient material showing a rearrangement involving the EVI gene region at 3q26 (1RG1R1G). Image kindly provided by Dr. Reed, London.

Literature:


Levy et al, 1994, Blood, 83: 1348-1354.
Wieser et al, 2003, Haematologica, 88: 25-30.
Melo et al, 2007, Leukemia, 22, 434-437.

Note: In t(3;3) the breakpoint cluster can span 1.3 Mb region. The described probe set will therefore provide false negative results in cases with very distal breakpoints.

| Ordering information | Color | Tests | Cat# |
|------------------------------------|-----------|-------|-----------|
| ON EVI t(3;3); inv(3) (3q26) Break | red/green | 10 | KBI-10204 |

The EVI t(3;3) inv(3) Break, triple-color probe is optimized to detect the inversion of chromosome 3 involving the EVI gene region at 3q26 in a dual-color, split assay on metaphase/interphase spreads, blood smears and bone marrow cells. By using a third color breakpoint variations can also be easily observed.

The diagram illustrates the structural rearrangement in the human chromosome 3p26.2 region. On the left, the normal configuration is shown with a blue chromosome and a red band at 3p26.2. On the right, the inverted configuration is shown with a blue chromosome and a green band at 3p26.2. The diagram includes labels for the RH103072, RH67219, EVI1, RH1348, RH10606, RH111834, and RH123089 genes, and the GAP and TERC regions. The inverted region is labeled 'inv(3) breakpoint region' and the normal region is labeled 't(3;3) breakpoint region'.



Literature:

Levy et al, 1994, Blood, 83: 1348-1354.
Wieser et al, 2003, Haematologica, 88: 25-30.
Melo et al, 2007, Leukemia, 22, 434-437.

27

Acute Myeloid Leukemia (AML)



At least 80% of patients with acute myeloid leukemia (AML) have an abnormal karyotype. Cytogenetic analysis provides some of the strongest prognostic information available, predicting outcome of both remission induction and postremission therapy. Abnormalities which indicate a good prognosis include t(8;21), inv(16), and t(15;17). Patients with AML that is characterized by deletions of the long arms or monosomies of chromosomes 5 or 7; by translocations or inversions of chromosome 3, t(6;9), t(9;22); or by abnormalities of chromosome 11q23 have particularly poor prognoses with chemotherapy.

ON AML/ETO t(8;21) Fusion

t(8;21)(q21;q22) is the most frequently observed karyotypic abnormality associated with acute myeloid leukemia (AML), especially in FAB M2. As a consequence of the translocation the AML1 (CBFA2, RUNX1) gene in the 21q22 region is fused to the ETO (MTG8, RUNX1T) gene in the 8q21 region, resulting in one transcriptionally active gene on the 8q-derivative chromosome.

The AML/ETO t(8;21)(q21;q22) specific DNA probe is optimized to detect the reciprocal translocation t(8;21) in a dual-color, dual-fusion assay.

Cat.# KBI-10301 AML/ETO t(8;21) Fusion



AML/ETO t(8;21) Fusion probe hybridized to a normal metaphase (2R2G).

Literature:

Sacchi et al, 1995, Genes Chrom Cancer, 79: 97-103.
Hagemeijer et al, 1998, Leukemia, 12: 96-101.



| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| ON AML/ETO t(8;21) Fusion | red/green | 10 | KBI-10301 |

ON PML/RARA t(15;17) Fusion

A structural rearrangement involving chromosomes 15 and 17 in acute promyelocytic leukemia (APL) was first recognized in 1977. The critical junction is located on the der(15) chromosome and consists of the 5' portion of PML fused to virtually all of the RARA gene. The PML/RARA fusion protein interacts with a complex of molecules known as nuclear co-repressors and histone deacetylase. This complex binds to the fusion protein and blocks the transcription of target genes. Other less common variant translocations fuse the RARA gene on 17q21 to the PLZF, NPM, NUMA, and STAT5b genes, respectively.

The PML/RARA t(15;17) specific DNA probe is optimized to detect the reciprocal translocation t(15;17) (q24;q21) in a dual-color, dual-fusion assay.

Cat.# KBI-10302 PML/RARA t(15,17) Fusion



PML/RARA t(15,17) Fusion probe hybridized to a normal metaphase (2R2G).

Literature:

Schad et al, 1994, Mayo Clin Proc, 69: 1047-1053.
Brockman et al, 2003, Cancer Genet Cytogenet, 145: 144-151.



| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON PML/RARA t(15,17) Fusion | red/green | 10 | KBI-10302 |
| ON PML/RARA t(15,17) Fusion | red/green | 20 | KBI-12302 |

ON MLL (11q23) Break

The human chromosome band 11q23 is associated with a high number of recurrent chromosomal abnormalities including translocations, insertions, and deletions. It is involved in over 20% of acute leukemias. The MLL (Myeloid-Lymphoid Leukemia or Mixed-Lineage Leukemia) gene, named for its involvement in myeloid (usually monoblastic) and lymphoblastic leukemia, and less commonly in lymphoma is located in the 11q23 breakpoint region. Leukemias involving the MLL gene usually have a poor prognosis.

The MLL (11q23) break probe is optimized to detect translocations involving the MLL gene region at 11q23 in a dual-color split assay.

Cat.# KBI-10303 MLL (11q23) Break



MLL (11q23) Break probe hybridized to patient material showing a translocation at 11q23 (1R1G1G).

Literature:

Kobayashi et al, 1993, Blood, 81: 3027-3022
Martinez-Climent et al, 1995, Leukemia, 9: 1299-1304.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON MLL (11q23) Break | red/green | 10 | KBI-10303 |

ON MLL/MLLT1 t(11;19) Fusion

One of the relatively frequently observed translocations (around 10 %) in human Acute Myeloid Leukemia (AML) and Acute Lymphoblastic Leukemia (ALL) involves the genes MLL and MLLT1 (aka ENL) at 11q23 and 19p13. The MLL/MLLT1 translocation results in the generation of fusion protein that retains the MLL N-terminus, including both an A-T hook domain and a region similar to mammalian DNA methyltransferase. There are several breakpoints within the MLLT1 gene described, without clear differences in clinicohematologic features. Patients with AML and the MLL/MLLT1 translocation carry a poor prognosis, but noninfant children with ALL and MLL/MLLT1 fusion may have a favorable prognosis.

The MLL/MLLT1 Fusion probe is optimized to detect translocations involving the MLL and MLLT1 gene regions at 11q23 and 19p13 in a dual-color, fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10307 MLL/MLLT1 t(11;19) Fusion



MLL/MLLT1 t(11;19) Fusion probe hybridized to patient material showing t(11;19) translocation (2RG1R1G).
Image kindly provided by Dr. Mohr, Dresden.

Literature:

Mitterbauer-Hohdanner G et al, 2004, Eur J Clin Invest, 34; 12-24.
Meyer C et al, 2009, Leukemia, 23; 1490-1499.
Fu JF et al, 2007, Am J Clin Pathol, 127; 24-30.



| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON MLL/MLLT1 t(11;19) Fusion | red/green | 10 | KBI-10307 |

ON MLL/MLLT3 t(9;11) Fusion

Chromosomal rearrangements involving the mixed lineage leukemia (MLL) gene at 11q23 are frequently observed in adult and childhood acute leukemia and are, in general, associated with poor prognosis. However, children with Acute Myeloid Leukemia (AML) carrying the t(9;11) MLL/MLLT3 (aka AF9) translocation have been described to be more sensitive to chemotherapy than patients with other 11q23 rearrangements.

The MLL/MLLT3 Fusion probe is optimized to detect translocations involving the MLL and MLLT3 gene regions at 11q23 and 9p21 in a dual-color fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10308 MLL/MLLT3 t(9;11) Fusion



MLL/MLLT3 t(9;11) Fusion probe hybridized to patient material showing t(9;11) translocation (2RG1R1G).
Image kindly provided by Dr. Mohr, Dresden.

Literature:

Palle J et al, 2005, Br J Haematol, 129; 189-198.
Meyer C et al, 2009, Leukemia, 23; 1490-1499.
Cavazzini F et al, 2006, Haematologica, 91; 381-5.
Balgobind BV et al, 2009, Blood, 114; 2489-96.
Keefe JG et al, 2010, J Mol Diagn, 12; 441-452.



| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON MLL/MLLT3 t(9;11) Fusion | red/green | 10 | KBI-10308 |

ON MLL/MLLT4 t(6;11) Fusion

One of the relatively frequently observed translocations in human Acute Myeloid Leukemia (AML) involves the genes MLL and MLLT4 (aka AF6) at 11q23 and 6q27. The MLL/MLLT4 translocation results in the generation of fusion protein that retains the MLL N-terminus, including both an A-T hook domain and a region similar to mammalian DNA methyltransferase. The breakpoint region of the MLLT4 gene is located within intron 1 and downstream of the initiation codon. In all age groups and all phenotypes of leukemia, the MLL/MLLT4 translocation carries a poor prognosis.

The MLL/MLLT4 t(6;11) Fusion probe is optimized to detect translocations involving the MLL and MLLT4 gene regions at 11q23 and 6q27 in a dual-color, fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat# KBI-10309 MLL/MLLT4 t(6;11) Fusion





MLL/MLLT4 t(6;11) Fusion probe hybridized to patient material showing 47,XX,t(6;11)(q27;q23),+der(6)t(6;11)(q27;q23).
Image kindly provided by Dr. Mohr, Dresden.

ON RARA (17q21) Break

This break apart probe can detect the numerous types of recurrent rearrangement of the RAR α (Retinoid acid receptor, alpha) gene with various gene partners (e.g., PML, NPM, MLL, FIP1L1, NuMA1, PLZF, amongst the others), leading to the formation of different reciprocal fusion proteins. The importance of retinoid metabolism in acute promyelocytic leukemia (APL) is highlighted by the numerous recent studies, but the different leukemogenic functions of the RAR α fusion proteins in the neoplastic myeloid development still has to be defined, as well as the distinct clinical outcome of the patients with the variant forms of APL.

Cat.# KBI-10305 RARA (17q21) Break



RARA (17q21) Break probe hybridized to patient material showing a translocation at 17q21 (1RG1R1G).

Literature:

Mitterbauer-Hohdanner G et al, 2004, Eur J Clin Invest, 34; 12-24.
Meyer C et al, 2009, Leukemia, 23; 1490-1499.

| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON MLL/MLLT4 t(6;11) Fusion | red/green | 10 | KBI-10309 |



| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON RARA (17q21) Break | red/green | 10 | KBI-10305 |

ON CBFB t(16;16); inv(16) Break

Inv(16)(p13;q22) and t(16;16)(p13;q22) are recurring chromosomal rearrangements in AML. In both the inversion and translocation, the critical genetic event is the fusion of the CBFB gene at 16q22 to the smooth muscle myosin heavy chain (MYH11) at 16p13. A deletion of between 150 and 350 kb centromeric to the p-arm inversion breakpoint cluster region can be observed in some patients containing the 5' portion of the myosin heavy chain (MYH11) gene.

The CBFB t(16;16) inv(16) break probe is optimized to detect the inversion of chromosome 16 involving the CBFB gene region at 16q22 in a dual-color, split assay.

Cat.# KBI-10304 CBFB t(16;16); inv(16) Break



CBFB t(16;16); inv(16) Break probe hybridized to a normal metaphase (2RG).

Literature:

Dauwerse et al, 1993, Hum.Mol.Genet., 2: 1527-1534.
Marlton et al, 1995, Blood, 85: 772-779.



| Ordering information | Color | Tests | Cat# |
|---------------------------------|-----------|-------|-----------|
| ON CBFB t(16;16); inv(16) Break | red/green | 10 | KBI-10304 |

ON DEK / NUP214 t(6;9) Fusion

The chromosomal translocation t(6;9) (p22;q34) is associated with a specific subtype of acute myeloid leukemia (AML) and constitutes 0.5% to 4% of all AML cases. The translocation results in a fusion between the DEK oncogene (6p22) and the nucleoporin 214 kDa (NUP214 at 9q34; previously known as CAN). The exact mechanism by which the fusion protein DEK-NUP214 contributes to leukemia development has not been identified. Patients with t(6;9) AML have a very poor prognosis. The currently available chemotherapy does not seem to improve overall survival. However, accurate diagnosis is crucial because these patients may benefit from early allogeneic stem cell transplant.

The DEK / NUP214 t(6;9) specific DNA Probe has been optimized to detect the reciprocal translocation t(6;9) in a dual-color, dual-fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

Cat.# KBI-10306 DEK / NUP214 t(6;9) Fusion probe



DEK / NUP214 t(6;9) Fusion probe hybridized to patient material showing a t(6;9)(p22;q34) translocation (2RG1R1G).
Image kindly provided by Dr. Stevens-Kroef, UMC St. Radboud, Nijmegen.

Literature:

Von Lindern et al, 1992, Mol. Cell. Biol., 12: 1687-1697.
Ageberg et al, 2008, Gen. Chrom. Canc., 47: 276-287.
Chi et al, 2008, Arch. Pathol. Lab. Med., 132: 1835-1837.

| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON DEK / NUP214 t(6;9) Fusion | red/green | 10 | KBI-10306 |

Acute Lymphoblastic Leukemia (ALL)



Acute lymphocytic leukemia, also called acute lymphoblastic leukemia, is a type of cancer that starts from white blood cells in the bone marrow. A number of recurring cytogenetic abnormalities are associated with distinct immunologic phenotypes of ALL and characteristic outcomes. The ETV6 (TEL) / RUNX1 (AML1) fusion arising from the translocation t(12;21)(p13;q22) has been associated with a good outcome; the BCR/ABL fusion of t(9;22)(q34;q11), rearrangements of the MLL gene (11q23), and abnormalities of the short arm of chromosomes 9 involving the tumor suppressor genes p16 (INK4A) have a poor prognosis.

ON TEL/AML t(12;21) Fusion

The t(12;21), a cryptic translocation rarely observed by conventional cytogenetics, was first identified by fluorescence in situ hybridization (FISH). In ALL blasts, this translocation fuses the 5' part of the TEL (ETV6) gene with almost the entire AML1 (CBFA2) gene, producing the chimeric transcript ETV6-CBFA2. The t(12;21) (p13;q22) has also been identified as the most frequent chromosomal abnormality in childhood ALL, affecting 20% to 25% of B-lineage cases.

The TEL/AML t(12;21) specific DNA probe is optimized to detect the reciprocal translocation t(12;21) (p13;q22) in a dual-color, dual-fusion assay.

Cat.# KBI-10401 TEL/AML t(12;21) Fusion



TEL/AML t(12;21) Fusion probe hybridized to patient material showing t(12;21)translocation (2RG1R1G).
Material kindly provided by Dr. Balogh, Budapest.



Literature:
Romana et al, 1995, Blood, 85: 3662-3670.

| Ordering information | Color | Tests | Cat# |
|----------------------------|-----------|-------|-----------|
| ON TEL/AML t(12;21) Fusion | red/green | 10 | KBI-10401 |

ON ETV6 (TEL) (12p13) Break

ETV6 (TEL) gene is the abbreviation for -ETS variant 6- gene. It encodes an ETS family factor which functions as a transcriptional repressor in hematopoiesis and in vascular development. The gene is located on chromosome 12p13, and is frequently rearranged in human leukemias of myeloid or lymphoid origins. Also systematic deletion of the normal ETV6 allele in patients with ETV6-AML1 fusions can be found.

Cat.# KBI-10403 ETV6 (TEL) (12p13) Break



ETV6 (TEL) (12p13) Break probe hybridized to patient material showing a translocation involving the ETV6 region at 12p13 (1RG1R1G). Image kindly provided by Magret Ratjen, Kiel.

Literature:

Golub et al, 1995, PNAS 92; 4917-4921.
Ford et al, 2001, Blood 98; 558-564.



| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON ETV6 (TEL) (12p13) Break | red/green | 10 | KBI-10403 |

ON p16 (9p21) / 9q21

Hemizygous deletions and rearrangements of chromosome 9, band p21 are among the most frequent cytogenetic abnormalities detected in pediatric acute lymphoblastic leukemia (ALL). This deletion includes loss of the p16 (CDKN2A)/p15 (CDKN2B) genes, which are cell cycle kinase inhibitors and important in leukemogenesis.

The p16 (9p21) specific DNA probe is optimized to detect copy numbers of the p16 (INK4A) gene region at region 9p21. The 9q21 region probe is included to facilitate chromosome identification.

Cat.# KBI-10402 p16 (9p21) / 9q21



p16 (9p21) / 9q21 hybridized on patient material showing an isochromosome 9. Image kindly provided by Dr. Wenzel, Basel.

Literature:

Dreyling et al, 1995, Blood, 86: 1931-1938.
Southgate et al, 1995, Br J Cancer, 72: 1214-1218.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON p16 (9p21) / 9q21 | red/green | 10 | KBI-10402 |

ON MLL (11q23) Break

The human chromosome band 11q23 is associated with a high number of recurrent chromosomal abnormalities including translocations, insertions, and deletions. It is involved in over 20% of acute leukemias. The MLL (Myeloid-Lymphoid Leukemia or Mixed-Lineage Leukemia) gene, named for its involvement in myeloid (usually monoblastic) and lymphoblastic leukemia, and less commonly in lymphoma is located in the 11q23 breakpoint region. Leukemias involving the MLL gene usually have a poor prognosis.

The MLL (11q23) break probe is optimized to detect translocations involving the MLL gene region at 11q23 in a dual-color split assay.

Cat.# KBI-10303 MLL (11q23) Break



MLL (11q23) Break probe hybridized to patient material showing a translocation at 11q23 (1RG1R1G).

Literature:

Kobayashi et al, 1993, Blood, 81: 3027-3022.

Martinez-Climent et al, 1995, Leukemia, 9: 1299-1304.

ON BCR/ABL t(9;22)

The t(9;22) BCR/ABL translocation is present in about 5% of pediatric and up to 50% of adult ALL cases, and is associated with poor prognosis.

Cat.# KBI-10005 ON BCR/ABL t(9;22) Fusion

Cat.# KBI-12005 ON BCR/ABL t(9;22) Fusion

Cat.# KBI-10006 ON BCR/ABL t(9;22) TC, D-Fusion

Cat.# KBI-10008 ON BCR/ABL t(9;22) DC, S-Fusion, ES

Cat.# KBI-10009 ON BCR/ABL t(9;22) DC, S-Fusion

Cat.# KBI-10013 ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES

See description under CML on page 11 and 12.

| Ordering information | Color | Tests | Cat# |
|--------------------------------------|----------------|-------|-----------|
| ON BCR/ABL t(9;22) Fusion | red/green | 10 | KBI-10005 |
| ON BCR/ABL t(9;22) Fusion | red/green | 20 | KBI-12005 |
| ON BCR/ABL t(9;22) TC, D-Fusion | red/green/blue | 10 | KBI-10006 |
| ON BCR/ABL t(9;22) DC, S-Fusion, ES | red/green | 10 | KBI-10008 |
| ON BCR/ABL t(9;22) DC, S-Fusion | red/green | 10 | KBI-10009 |
| Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | red/green | 10 | KBI-10013 |

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON MLL (11q23) Break | red/green | 10 | KBI-10303 |

Multiple Myeloma (MM)

The cytogenetic picture in multiple myeloma (MM) is highly complex, from which non-random numerical and structural chromosomal changes have been identified. Specifically, translocations involving the immunoglobulin heavy chain gene (IGH) at 14q32 and either monosomy or deletions of chromosome 13 have been observed in a significant number of patients. More recently several additional deletions or amplifications have been identified in MM which are currently investigated in large patient studies.



Note: Multiple Myeloma is a cancer of plasma cells. Analysis of such cells is hampered by their low frequency. Enrichment of plasma cells using CD138 is highly recommended.

ON MYEOV/IGH t(11;14) Fusion

The most common chromosomal translocation in multiple myeloma (MM) is t(11;14), resulting in up-regulation of cyclin D1. In MM the breakpoints are scattered within a 360-kb region between CCND1 and MYEOV. This breakpoint is more proximal than the t(11;14) breakpoints observed in mantle cell lymphoma or other leukemias. Patients with MM who have t(11;14)(q13;q32) seem to have an aggressive clinical course.

The MYEOV/IGH t(11;14)(q13;q32) specific DNA probe is optimized to detect the reciprocal translocation t(11;14) in a dual-color, dual-fusion assay.

Cat.# KBI-10605 MYEOV/IGH t(11;14) Fusion



MYEOV/IGH t(11;14) Fusion probe hybridized to MM patient material showing t(11;14) translocation (2RG1R1G).

Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

Janssen et al., 2000, Blood, 95: 2691-2698.
Fonseca et al, 2002, Blood, 99: 3735-3741.


| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON MYEOV/IGH t(11;14) Fusion | red/green | 10 | KBI-10605 |

ON MAF/IGH t(14;16) Fusion

Chromosome translocations involving the immunoglobulin heavy chain gene (IgH) on 14q32 are a fundamental event in the pathogenesis of many B-cell malignancies. It often is preceded by a stable pre-malignant tumor called Monoclonal Gammopathy of Undetermined Significance (MGUS), which can sporadically progress to Multiple Myeloma (MM). One of the recurrent primary rearrangements involving the immunoglobulin heavy chain (IgH) locus on chromosome 14q32 identified in MGUS and MM tumors is the MAF/IgH t(14;16) translocation. Following MGUS appearance, the pathogenesis of multiple myeloma (MM) is thought to involve at least two pathways, which generate hyperdiploid (HRD) or nonhyperdiploid (NHRD) tumors, respectively.

The MAF/IgH is mainly present in NHRD tumors, providing important information on MM patient sub-types. Since these translocations are caused by aberrant IgH switch recombination, and possibly by aberrant somatic hypermutation in germinal center B cells, they provide information of an early and perhaps initiating event of transformation.

Cat.# KBI-10610 MAF/IGH t(14;16) Fusion



MAF/IGH t(14;16) Fusion probe hybridized to patient material showing a deletion of the MAF gene region at 16q23 (2R1G).

Literature:
Chesi et al, 1998, Blood 91; 4457-4463.
Sawyer et al, 1998, Blood 92; 4269-4278.



| Ordering information | Color | Tests | Cat# |
|----------------------------|-----------|-------|-----------|
| ON MAF/IGH t(14;16) Fusion | red/green | 10 | KBI-10610 |

ON FGFR3/IGH t(4;14) Fusion

The t(4;14) translocation is undetectable by conventional cytogenetics. The breakpoints on chromosome 4 occur within an approximately 113-kb region located in small part of a conserved gene cluster including the transforming acidic coiled-coil protein 3 (TACC3), fibroblast growth factor receptor 3 (FGFR3), and multiple myeloma SET domain-containing protein (MMSET). The translocation is indicative for poor survival and poor response to chemotherapy.

The FGFR3/IGH t(4;14)(p16;q32) specific DNA probe is optimized to detect the reciprocal translocation t(4;14) in a dual-color, dual-fusion assay.

Cat.# KBI-10602 FGFR3/IGH t(4;14) Fusion



FGFR3/IGH t(4;14) Fusion probe hybridized to MM patient material showing t(4;14) translocation (2R61R1G).
Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

Chesi et al, 1997, Nat Genet, 16: 260-264.
Finelli et al, 1999, Blood, 94: 724-732.



| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON FGFR3/IGH t(4;14) Fusion | red/green | 10 | KBI-10602 |

ON MM 11q23 / DLEU (13q14)

Hybridization results delineated 11q23 and 11q25 as the most frequently gained regions in MM, which supports a relevant pathogenetic role of genes in this region. Deletions of 13q14 are frequently detected by interphase FISH in patients with newly diagnosed MM, correlate with increased proliferative activity, and represent an independent adverse prognostic feature in MM.

The 11q23 specific DNA probe is optimized to detect copy numbers at 11q23. The DLEU (13q14) specific DNA region is optimized to detect copy numbers of the DLEU gene region at 13q14.

Cat.# KBI-10502 MM 11q23 / 13q14



MM 11q23 / DLEU 13q14 probe hybridized to MM patient material showing a 13q14 deletion (1R2G).
Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

Gonzalez et al, 2004, Haematologica, 89: 1213-1218.
Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.



| Ordering information | Color | Tests | Cat# |
|----------------------------|-----------|-------|-----------|
| ON MM 11q23 / DLEU (13q14) | red/green | 10 | KBI-10502 |

ON MM 1q21 / 8p21

Amplifications of 1q21 are concurrent with dysregulated expression of c-MAF, MMSET/FGFR3, or Del13 and represent an independent unfavorable prognostic factor. Allelic losses of the chromosome 8p21-22 have been reported as a frequent event in several cancers.

The 1q21 specific DNA probe is optimized to detect copy numbers at 1q21. The 8p21 specific DNA region is optimized to detect copy numbers at 8p21.

Cat.# KBI-10503 MM 1q21 / 8p21



MM 1q21 / 8p21 hybridized to a normal metaphase (2R2G).

Literature:

Shaughnessy J., 2005, Hematology, 10 suppl 1: 117-126.
Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON MM 1q21 / 8p21 | red/green | 10 | KBI-10503 |

ON MM 19q13 / p53 (17p13)

P53 gene deletion, which can be identified by interphase FISH in almost a third of patients with newly diagnosed MM, is a novel prognostic factor predicting for short survival of MM patients treated with conventional-dose chemotherapy. Amplification of 19q13 has been reported in a variety of cancer.

The 19q13 specific DNA probe is optimized to detect copy numbers at 19q13. The p53 (17p13) specific DNA region is optimized to detect copy numbers of the p53 gene region at 17p13.

Cat.# KBI-10509 MM 19q13 / p53 (17p13)



MM 19q13 / p53 (17p13) hybridized to patient material showing a p53 (17p13) deletion (1R2G).

Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

Drach et al, 1998, Blood, 92: 802-809.
Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.



| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| ON MM 19q13 / p53 (17p13) | red/green | 10 | KBI-10509 |

ON MM 15q22 / 6q21

Chromosome 6q amplifications encompassing 6q21-22 have been observed in MM including the same region as in CLL. Amplification including band 15q22 has been reported in MM.

The 15q22 specific DNA probe is optimized to detect copy numbers at 15q22. The 6q21 specific DNA region is optimized to detect copy numbers at 6q21.

Cat.# KBI-10504 MM 15q22 / 6q21



MM 15q22 / 6q21 hybridized to MM patient material with gain of both critical regions 6q21 and 15q22.
Image kindly provided by Prof. Jauch, Heidelberg.

Literature:

Cremer et al, 2005,
Genes Chrom Cancer, 44: 194-203.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON MM 15q22 / 6q21 | red/green | 10 | KBI-10504 |

ON MM 15q22 / 9q34

The hyperdiploid subtype in MM is defined by presence of multiple trisomic chromosomes. Combination of the chromosome 9q34 and 15q22 specific regions are important regions to detect the hyperdiploid subtype in MM which is usually associated with a low frequency of IGH translocations.

The 15q22 and 9q34 probe is designed as a dual-color assay to detect amplifications at 15q22 and 9q34.

Cat.# KBI-10508 MM 15q22 / 9q34



MM 15q22 / 9q34 hybridized to a normal interphase/
metaphase (2R2G).

Literature:

Cremer et al, 2005,
Genes Chrom Cancer, 44: 194-203.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON MM 15q22 / 9q34 | red/green | 10 | KBI-10508 |

ON MM 1q21 / SRD (1p36)

Frequent loss of heterozygosity (LOH) on the short arm of chromosome 1 (1p) has been reported in a series of human malignancies.

The 1q21 specific DNA probe is optimized to detect copy numbers at 1q21. The SRD 1p36 specific DNA Probe is optimized to detect copy numbers of 1p at region 1p36 containing the markers D1S2795 and D1S253.

Cat.# KBI-10507 MM 1q21 / SRD (1p36)



MM 1q21 / SRD (1p36) hybridized to a normal metaphase (2R2G).

Literature:

Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.
Shaughnessy J., 2005, Hematology, 10 suppl 1: 117-126.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON 1q21 / SRD (1p36) | red/green | 10 | KBI-10507 |

ON IGH (14q32) Break

Multiple myeloma is characterized by complex rearrangements involving the IgH gene, particularly at the constant locus. The IgH rearrangement provides a useful marker of clonality in B-cell malignancies and amplification of this rearrangement is the method of choice to monitor the residual tumor cells in multiple myeloma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14q32 in a dual-color, split assay.

Cat.# KBI-10601 IGH (14q32) Break



IGH (14q32) Break probe hybridized to patient material showing a partial deletion of 14q32 (1RG1R).

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969.
Gozetti et al, 2002, Cancer Research, 62: 5523-5527.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON IGH (14q32) Break | red/green | 10 | KBI-10601 |

Lymphoma



Lymphoma is a type of cancer that originates in lymphocytes. Following WHO classification there are three large groups: B cell lymphomas, T cell and natural killer cell tumors, and Hodgkin lymphoma. The IGH gene locus at chromosome band 14q32 is very frequently involved in B-cell lymphoma.

ON MYC/IGH t(8;14) Fusion

The translocation t(8;14)(q24;q32) is the characteristic chromosomal aberration of Burkitt's-type of lymphomas. This translocation fuses the C-MYC gene at 8q24 next to the IGH locus at 14q32, resulting in overexpression of the transcription factor C-MYC. Detection of the t(8;14) is aimed to help in the diagnostic process of patients with high-grade B-cell lymphomas because treatment strategies differ between Burkitt and other high-grade lymphomas.

The MYC/IGH t(8;14)(q24;q32) specific DNA probe is optimized to detect the reciprocal translocation t(8;14) in a dual-color, dual-fusion assay.

Cat.# KBI-10603 MYC/IGH t(8;14) Fusion



MYC/IGH t(8;14) Fusion probe hybridized to a normal interphase/metaphase (2R2G).

Literature:

Veronese et al, 1995, Blood, 85: 2132-2138.
Siebert et al, 1998, Blood, 91: 984-990.

| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| ON MYC/IGH t(8;14) Fusion | red/green | 10 | KBI-10603 |

ON BCL1/IGH t(11;14) Fusion

Mantle cell lymphoma is a subtype of non-Hodgkin lymphoma characterized by poor prognosis. Cytogenetically t(11;14) is associated with 75% of mantle cells lymphomas. The translocation breakpoints are scattered within the 120 kb BCL1 region adjacent to CCND1. The translocation leads to overexpression of cyclin D1 due to juxtaposition of the Ig heavy chain gene enhancer on 14q32 to the cyclin D1 gene on 11q13.



The BCL1/IGH t(11;14)(q13;q32) specific DNA probe is optimized to detect the reciprocal translocation t(11;14) in a dual-color, dual-fusion assay.

ON BCL2/IGH t(14;18) Fusion

The t(14;18) chromosomal translocation that results in the juxtaposition of the BCL2 proto-oncogene with the heavy chain JH locus. It a common cytogenetic abnormality in human lymphoma and is observed in about 85% of follicular lymphoma (FL) and up to one-third of diffuse lymphomas (DL). Two breakpoint region clusters (brc) have been identified: a major breakpoint region (mbr) within the 3' untranslated region of the BCL2 proto-oncogene accounting for approximately 60% of the cases and a minor cluster region (mcr) 30 kb 3' of BCL2 accounting for approximately 25% of the breakpoints.



The BCL2/IGH t(14;18)(q21;q32) specific DNA probe is optimized to detect the reciprocal translocation t(18;14), involving either of the two brc in the BCL2 gene in a dual-color, dual-fusion assay. KREATECH has optimized this probe for the specific use on cell material (KBI-10606), or for the use on tissue (KBI-10755).

Cat.# KBI-10604 BCL1/IGH t(11;14) Fusion



BCL1/IGH t(11;14) Fusion probe hybridized to a normal interphase/metaphase (2R2G).

Cat.# KBI-10606 BCL2/IGH t(14;18) Fusion



BCL2/IGH t(14;18) probe hybridized to a normal interphase/metaphase (2R2G).

Literature:

Vaandrager et al, 1996, Blood, 88: 1177-1182.

Literature:

Poetsch et al, 1996, J Clin Oncol, 14: 963-969.

Vaandrager et al, 2000, Genes Chrom Cancer, 27: 85-94.

| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON BCL1/IGH t(11;14) Fusion | red/green | 10 | KBI-10604 |



| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON BCL2/IGH t(14;18) Fusion | red/green | 10 | KBI-10606 |

ON IGH (14q32) Break

Chromosomal rearrangements involving the immunoglobulin heavy chain gene (IGH) at 14q32 are observed in 50% of patients with B-cell non-Hodgkin's lymphoma (NHL) and many other types of Lymphomas. More than 50 translocation partners with IGH have been described. In particular t(8;14), is associated with Burkitt's lymphoma, t(11;14) is associated with Mantle cell lymphoma, t(14;18) is observed in a high proportion of follicular lymphomas and t(3;14) is associated with Diffuse Large B-Cell Lymphoma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14q32 in a dual-color, split assay. Kreatech has developed this probe for the specific use on cell material (KBI-10601), or for the use on tissue (KBI-10729).

Cat.# KBI-10601 IGH (14q32) Break



IGH (14q32) Break probe hybridized to patient material showing a partial deletion of 14q32 (1R1R).

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969.
Gozetti et al, 2002, Cancer Research, 62: 5523-5527.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON IGH (14q32) Break | red/green | 10 | KBI-10601 |

ON BCL2 (18q21) Break

Follicular lymphoma is a mature B-cell lymphoma characterized by the presence of the t(14;18) translocation that juxtaposes the BCL2 locus on chromosome 18q21 to the immunoglobulin H (IGH) locus on chromosome 14q32, resulting in the overexpression of the anti-apoptotic protein BCL2. Next to IGH, other translocation partners to BCL2 are also known (e.g. IGK at 2p11.2 and IGL at 22q11). A break or split assay is therefore best suited to detect rearrangements of the BCL2 gene region at 18q21.

The BCL2 (18q21) Break probe is optimized to detect translocations involving the BCL2 gene region at 18q21 in a dual-color, split assay on metaphase/interphase spreads, bloodsmears and bone marrow cells.

Cat# KBI-10612 BCL2 (18q21) Break



BCL2 (18q21) Break probe hybridized to a normal metaphase.

Literature:

Taniwaki M et al, 1995, Blood, 86; 1481-1486.
Poetsch M et al, 1996, J Clin Oncol, 14; 963- 969.
Einerson R et al, 2005, Am J Clin Pathol, 124; 421-429.



| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON BCL2 (18q21) Break | red/green | 10 | KBI-10612 |

ON BCL6 (3q27) Break

Chromosomal translocations involving band 3q27 with various different partner chromosomes represent a recurrent cytogenetic abnormality in B-cell non-Hodgkin's lymphoma. A FISH strategy using two differently labeled flanking BCL6 probes provides a robust, sensitive, and reproducible method for the detection of common and uncommon abnormalities of BCL6 gene in interphase nuclei. Kreatech has developed this probe for the specific use on cell material (KBI-10607), or for the use on tissue (KBI-10730).

Two different breakpoint regions have been identified; the major breakpoint region (**MBR**) is located within the 5' noncoding region of the BCL6 proto-oncogene, while the atypical breakpoint region (**ABR**) is located approximately 200 kb distal to the BCL6 gene. The BCL6 (3q27) Break probe is designed in a way to flank both possible breakpoints, thereby providing clear split signals in either case.

Cat.# KBI-10607 BCL6 (3q27)



BCL6 (3q27) Break probe hybridized to patient material (1RG1R1G).
Image kindly provided by Prof. Siebert, Kiel.

Literature:

Butler et al, 2002, Cancer Res, 62; 4089-4094.
Sanchez-Izquierdo, 2001, Leukemia, 15; 1475-1484.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON BCL6 (3q27) Break | red/green | 10 | KBI-10607 |

ON MALT (18q21) Break

Low grade malignant lymphomas arising from mucosa associated lymphoid tissue (MALT) represent a distinct clinicopathological entity. The three major translocations seen in MALT lymphomas are t(11;18)(q21;q21)/API2-MALT1, t(14;18)(q32;q21)/IGH-MALT1 and t(1;14)(p22;q32)/IGH-BCL10. A break or split probe for MALT (18q21) is best used to analyze translocation of the MALT gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Kreatech has optimized this probe for the specific use on cell material (KBI-10608), or for the use on tissue (KBI-10731).

Cat.# KBI-10608 MALT (18q21) Break



MALT (18q21) Break probe hybridized to patient material showing a translocation at 18q21 (1RG1RG).

Literature:

Morgan et al, 1999, Cancer Res, 59; 6205-6213.
Dierlamm et al, 2000, Blood, 96; 2215-2218.


| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON MALT (18q21) Break | red/green | 10 | KBI-10608 |

ON CCND1 (BCL1;11q13) Break

Besides the important functions in cellular growth, metabolism, and cellular differentiation, CCND1 (also known as Cyclin D1 or BCL1) can also function as a proto-oncogene, often dysregulated after re-arrangement by translocation. In fact, it can juxtapose into many different gene locus to drive tumorigenic effects. To date, the gene has been found to be rearranged in leukemias, in multiple myelomas (MM), and in some cases of benign parathyroid tumors. More specifically, the chromosomal translocation t(11;14)(q13;q32), involving rearrangement of the CCND1 locus, has been reported to be associated with human lymphoid neoplasia affecting mantle cell lymphomas (MCL).

The rearrangement has been documented in 40-70% of cases of mantle cell lymphoma, whereas it only rarely occurs in other B cell lymphomas. In multiple myeloma, the same translocation t(11;14)(q13;q32) is the most common, with a reported frequency of 15% to 20% of the cases. For this reason, the CCND1 break apart probe KBI-10609 can be considered a very useful tool for routine diagnosis in MCL and Multiple myeloma, to be used in association to the related probes KBI-10604 and KBI-10605 probes that can detect more specifically the translocation t(11;14) in Mantle Cell Lymphoma (KBI-10604) and Multiple Myeloma (KBI-10605).

Cat.# KBI-10609 CCND1 (BCL1; 11q13) Break



CCND1 (BCL1; 11q13) Break probe hybridized to a normal metaphase (2R2G).

Literature:



Vaandrager et al, 1996, Blood, 88 (4); 1177-1182.
Vaandrager et al, 1997, Blood, 89; 349-350.

| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON CCND1 (BCL1;11q13) Break | red/green | 10 | KBI-10609 |

ON MYC (8q24) Break, TC

Rearrangements of the protooncogene C-myc (or MYC) have been consistently found in tumor cells of patients suffering from Burkitt's lymphoma. In cases with the common t(8;14) chromosomal translocation, the c-myc gene is translocated to chromosome 14 and rearranged with the immunoglobulin heavychain genes; the breakpoint occurs 5' to the c-myc gene and may disrupt the gene itself separating part of the first untranslated exon from the remaining two coding exons. In Burkitt's lymphoma showing the variant t(2;8) or t(8;22) translocations, the genes coding for the k and l immunoglobulin light chain are translocated to chromosome 8. The rearrangement takes place 3' to the c-myc gene.

Cat.# KBI-10611 MYC (8q24) Break, TC



MYC (8q24) Break probe hybridized to patient material showing a 8q24 proximal break (1GBR1G1BR).
Image kindly provided by Prof. Siebert, Kiel.

Literature:

Fabris et al, 2003, Genes Chromosomes Cancer 37 ; 261-269.
Hummel et al., 2006, N Engl J Med 354 ;2419-30.

| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON MYC (8q24) Break, TC | red/green | 10 | KBI-10611 |

The MYC (8q24) break-apart probe is optimized to detect rearrangements involving the 8q24 locus in a triple-color, split assay on metaphase/interphase spreads, blood smears and bone marrow cells.

ON FGFR3/IGH t(4;14) Fusion

Cat.# KBI-10602 FGFR3/IGH t(4;14) Fusion

See description under Multiple Myeloma on page 38.

| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON FGFR3/IGH t(4;14) Fusion | red/green | 10 | KBI-10602 |

ON MYEOV/IGH t(11;14) Fusion

Cat.# KBI-10605 MYEOV/IGH t(11;14) Fusion

See description under Multiple Myeloma on page 36.

| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON MYEOV/IGH t(11;14) Fusion | red/green | 10 | KBI-10605 |

ONCOLOGY - SOLID TUMOR DNA PROBES

In solid tumors significantly high levels of chromosome abnormalities have been detected, but distinction between significant (drivers) and irrelevant events (passengers) has been a major challenge. Consequently, the application of cytogenetic technology as diagnostic, prognostic, or therapeutic tools for these malignancies has remained largely underappreciated. The emergence of FISH is particularly useful for solid malignancies, and the spectrum of their application is rapidly expanding to improve efficiency and sensitivity in cancer diagnosis, prognosis, and therapy selection, alone or with other diagnostic methods. The REPEAT-FREE™ POSEIDON™ Solid Tumor DNA Probes are direct labeled, Ready-to-Use in hybridization buffer. The solid tumor probes are designed for the use on Formalin Fixed Paraffin Embedded (FFPE) tissue samples.

Table of Contents

| Description | Cat# | Page |
|---|-----------|--------|
| ON ERBB2, Her-2/Neu (17q12) / SE 17 | KBI-10701 | 49 |
| ON ERBB2, Her-2/Neu (17q12) / SE 17 | KBI-14701 | 49 |
| ON TOP2A (17q21) / SE 17 | KBI-10724 | 49 |
| ON ZNF217 (20q13) / 20q11 | KBI-10733 | 50 |
| ON TOP2A (17q21) / Her-2/neu (17q12) / SE 17 TC | KBI-10735 | 50 |
| ON FGFR1 (8p12) Break | KBI-10737 | 16, 51 |
| ON p16 (9p21) / 9q21 (tissue) | KBI-10710 | 52 |
| ON CC hTERT (3q26) C-MYC (8q24) / SE 7 TC | KBI-10704 | 52 |
| ON EGFR, Her-1 (7p11) / SE 7 | KBI-10702 | 54 |
| ON hTERT (5p15) / 5q31 (tissue) | KBI-10709 | 55 |
| ON C-MET (7q31) / SE 7 | KBI-10719 | 55 |
| ON ALK/EML4 t(2;2); inv(2) Fusion | KBI-10746 | 54 |
| ON ALK (2p23) Break | KBI-10747 | 53, 73 |
| ON ROS1 (6q22) Break | KBI-10752 | V |
| ON RET (10q11) Break | KBI-10753 | V |
| ON FGFR1 (8p11) / SE 8 (D8Z1) | KBI-12754 | VI |
| ON PTEN (10q23) / SE 10 | KBI-10718 | 57 |
| ON AR (Xq12) / SE X | KBI-10720 | 56 |
| ON TMPRSS2-ERG (21q22) Del, Break, TC | KBI-10726 | 57 |
| ON PPARγ (3p25) Break | KBI-10707 | 58 |
| ON MYCN (2p24) / LAF (2q11) | KBI-10706 | 59 |
| ON MLL (11q23) / SE 11 | KBI-10711 | 60 |
| ON SRD (1p36) / SE 1(1qh) | KBI-10712 | 60 |



| Description | Cat# | Page |
|---------------------------------------|-----------|--------|
| ON MDM4 (1q32) / SE 1 | KBI-10736 | 61 |
| ON SYT (18q11) Break | KBI-10713 | 63 |
| ON CHOP (12q13) Break | KBI-10714 | 63 |
| ON FUS (16p11) Break | KBI-10715 | 64 |
| ON FKHR (13q14) Break | KBI-10716 | 64 |
| ON MDM2 (12q15) / SE 12 | KBI-10717 | 65 |
| ON CDK4 (12q13) / SE 12 | KBI-10725 | 66 |
| ON COL1A1/PDGFB t(17;22) DC, S-Fusion | KBI-10742 | VI |
| ON EWSR1 (22q12) Break | KBI-10750 | 62 |
| ON EWSR1/NFATC2 t(20;22) DC, S-Fusion | KBI-10751 | VII |
| ON MYC (8q24) Break, TC (tissue) | KBI-10749 | VII |
| ON BCL2/IGH@ t(14;18) Fusion (tissue) | KBI-10755 | VIII |
| ON AURKA (20q13) / 20q11 | KBI-10721 | 67 |
| ON AURKB (17p13) / SE 17 | KBI-10722 | 68 |
| ON IGH (14q32) Break (tissue) | KBI-10729 | 70 |
| ON BCL6 (3q27) Break (tissue) | KBI-10730 | 71 |
| ON MALT (18q21) Break (tissue) | KBI-10731 | 71 |
| ON CCND1 (11q13) / SE 11 | KBI-10734 | 69 |
| ON p53 (17p13) / SE 17 (tissue) | KBI-10738 | 72 |
| ON ERCC1 (19q13) / ZNF443 (19p13) | KBI-10739 | 66 |
| ON TFE3 (Xp11) Break | KBI-10741 | 72 |
| ON BCL2 (18q21) Break (tissue) | KBI-10745 | 70 |
| ON ALK (2p23) Break | KBI-10747 | 53, 73 |

Breast Cancer

ON ERBB2, Her2/Neu (17q12) / SE 17

The Her2/Neu gene encodes a receptor tyrosine kinase involved in growth factor signaling. Overexpression of this gene is seen in about 20% of invasive breast cancers and is without proper treatment associated with poor survival. Her2 gene amplification is a permanent genetic change that results in this continuous overexpression of Her2. Trastuzumab (commonly known as Herceptin) has been developed to be effective against Her2-positive breast cancer. Her2/Neu amplification is also observed in a variety of other tumors, such as prostate, lung, colon and ovary carcinoma. The Her2/Neu (17q12) specific DNA probe is optimized to detect copy numbers of the Her2/Neu (ERBB2) gene region at region 17q12. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification/enumeration.

Cat.# KBI-10701 ERBB2, Her2/Neu (17q12) / SE 17



ERBB2, Her2/Neu (17q12) / SE 17 probe hybridized to breast tumor tissue showing amplification of Her2/Neu (ERBB2) / SE 17.

Literature:

Pauletti et al, 1996, Oncogene, 13: 63-72.
Xing, et al, 1996, Breast Cancer Res Treat, 39: 203-212.



| Ordering information | Color | Tests | Cat# |
|------------------------------------|-----------|-------|-----------|
| ON ERBB2, Her2/Neu (17q12) / SE 17 | red/green | 10 | KBI-10701 |
| ON ERBB2, Her2/Neu (17q12) / SE 17 | red/green | 50 | KBI-14701 |

ON TOP2A (17q21) / SE 17

The Topoisomerase2A enzyme, which is vital for the cell because of its role in cell replication and repair, catalyzes the relaxation of supercoiled DNA molecules to create a reversible double-strand DNA break. This enzyme is also the target of a number of cytotoxic agents, namely TOP2A inhibitors (anthracyclines, etoposide, teniposide).

The dual-color probe is optimized to detect amplifications (copy numbers) or deletions of the TOP2A gene region at the 17q21. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

Cat.# KBI-10724 TOP2A (17q21) / SE 17



TOP2A (17q21) / SE 17 probe hybridized to breast tissue (2R2G).

Literature:

Järvinen et al, 1999, Genes, Chromosomes and Cancer 26; 142-150.
Järvinen et al, 2000, Am. J. Pathology 156; 639-647.



| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON TOP2A (17q21) / SE 17 | red/green | 10 | KBI-10724 |

ON TOP2A (17q21) / Her2 (17q12) / SE 17 Triple-Color Probe

The presence of both TOP2A amplification and deletion in advanced cancer are associated with decreased survival, and occur frequently and concurrently with Her2 gene amplification.

The TOP2A (17q21)/ Her2 (17q12)/ SE 17 probe is designed as a triple-color assay to detect amplification at 17q12 as well as amplifications or deletions at 17q21. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 in blue is included to facilitate chromosome identification/enumeration.

Cat.# KBI-10735 TOP2A (17q21) / Her2 (17q12) / SE 17



TOP2A (17q21)/ Her2(17q12) / SE 17 TC probe hybridized to breast tumor tissue showing amplification of TOP2A/Her2.

Literature:

Järvinen et al, 1999, Genes, Chromosomes and Cancer 26; 142-150.
Järvinen et al, 2000, Am. J. Pathology 156; 639-647.



| Ordering information | Color | Tests | Cat# |
|----------------------------------|-----------|-------|-----------|
| ON TOP2A (17q21) / Her 2 / SE 17 | red/green | 10 | KBI-10735 |

ON ZNF217 (20q13) / 20q11

Zinc-finger protein 217 (ZNF217) is a Kruppel-like zinc-finger protein located at 20q13.2, within a region of recurrent maximal amplification in a variety of tumor types, and especially breast cancer cell lines and primary breast tumors. Copy number gains at 20q13 are also found in more than 25% of cancers of the ovary, colon, head and neck, brain, and pancreas, often in association with aggressive tumor behavior. ZNF217 is considered a strong candidate oncogene that may have profound effects on cancer progression, which is transcribed in multiple normal tissues, and overexpressed in almost all cell lines and tumors in which it is amplified.

The ZNF217 (20q13) specific DNA probe is optimized to detect copy numbers of 20q at 20q13. The 20q11 probe is included to facilitate chromosome identification.

Cat.# KBI-10733 ZNF217 (20q13) / 20q11



ZNF217 (20q13) / 20q11 probe hybridized to tissue (2R2G).

Literature:

Tanner M et al, 2000, Clin Cancer Res, 6; 1833-1839.
Ginestier C et al, 2006, Clin Cancer Res, 12; 4533-4544.

| Ordering information | Color | Tests | Cat# |
|---------------------------|-----------|-------|-----------|
| ON ZNF217 (20q13) / 20q11 | red/green | 10 | KBI-10733 |

ON FGFR1 (8p12) Break



Translocations affecting the chromosomal locus FGFR1 (8p12) are hallmarks of an atypical stem cell myeloproliferative disorder. These events disrupt the fibroblast growth factor receptor 1 (FGFR1) gene and fuse the FGFR1 C-terminal catalytic domain with unrelated proteins.

FGFR1 expression has been shown to play pivotal roles in mammary development and breast cancer tumorigenesis. It has been shown that FGFR1 amplification is found in up to 10% of breast cancers and is significantly more prevalent in patients older than 50 years of age and in tumors that lack HER2 expression.

Even though the prognostic impact of FGFR1 amplification in breast cancer still remains unclear, the functional data demonstrating that FGFR1 signaling is required for the survival of breast cancer cells harboring FGFR1 amplification and the relatively high prevalence of FGFR1 amplification in breast cancer support the idea that this gene may be a useful therapeutic target for a subgroup of breast cancer patients with FGFR1 gene amplification.

The FGFR1 (8p12) break-apart probe is optimized to detect translocations involving the FGFR1 gene region at 8p12 in a dual-color, split assay on metaphase/interphase spreads and paraffin embedded tissue sections.

Cat.# KBI-10737 FGFR1 (8p12) Break



FGFR1 (8p12) Break probe hybridized to patient material showing a break at 8p12 (1RG1R1G).

Literature:

Smedley et al, 1998, Hum Mol Genet. 7; 627-642.
Sohal et al, 2001, Genes Chrom. Cancer 32; 155-163.
Kwak et al, J Clin Oncol., 27(26):4247-53.

| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON FGFR1 (8p12) Break | red/green | 10 | KBI-10737 |



Bladder Cancer

ON p16 (9p21) / 9q21 (tissue)

Homozygous and hemizygous deletions of 9p21 are the earliest and most common genetic alteration in bladder cancer. The p16 (INK4A) gene has been identified as tumor suppressor gene in this region which is commonly deleted in bladder cancer. The loss of DNA sequences on chromosomal bands 9p21-22 has been documented also in a variety of malignancies including leukemias, gliomas, lung cancers, and melanomas.

The p16 (9p21) specific DNA probe is optimized to detect copy numbers of the p16 gene region at region 9p21. The 9q21 region probe is included to facilitate chromosome identification.

Cat.# KBI-10710 p16 (9p21) / 9q21 (tissue)



p16 (9p21) / 9q21 (tissue) probe hybridized to tissue (2R2G).

Literature:

Stadler et al, 1994, Cancer Res, 54: 2060-2063.
Williamson et al, 1995, Hum Mol Genet, 4: 1569-1577.



| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON p16 (9p21) / 9q21 (tissue) | red/green | 10 | KBI-10710 |

Cervical Cancer

ON CC hTERT (3q26) / C-MYC (8q24) / SE 7 TC

Cervical cancer, a potentially preventable disease, remains the second most common malignancy in women worldwide. The most consistent chromosomal gain in aneuploid tumors of cervical squamous cell carcinoma mapped to chromosome arm 3q, including the human telomerase gene locus (hTERT) at 3q26. High-level copy number increases were also mapped to chromosome 8. Integration of HPV (Human Papilloma Virus) DNA sequences into C-MYC chromosomal regions have been repeatedly observed in cases of invasive genital carcinomas and in cervical cancers. The hTERT (3q26) specific DNA Probe is optimized to detect copy numbers of the hTERT gene region at region 3q26. The C-MYC (8q24) specific DNA probe is optimized to detect copy numbers of the C-MYC gene region at 8q24. The chromosome 7 satellite enumeration probe (SE 7) at D7Z1 is included as ploidy control.

Cat.# KBI-10704 Cervical Cancer hTERT (3q26) / C-MYC (8q24) / SE 7 Triple-Color



CC hTERT (3q26) / C-MYC (8q24) / SE 7 TC probe hybridized to a PAP smear (destained) showing 3q26 and 8q24 amplification. The SE 7 control probe indicates a non-triploid karyotype (2B).

Image kindly provided by Dr. Weimer, Kiel.

Literature:

Xie et al, 2008, Geburtshilfe Frauenheilkunde, 68: 573.
Heselmeyer et al, 1996, PNAS, 93: 479-484.
Herrick et al, 2005, Cancer Res, 65: 1174-1179.

| Ordering information | Color | Tests | Cat# |
|---|-----------|-------|-----------|
| ON CC hTERT (3q26) / C-MYC (8q24) / SE 7 TC | red/green | 10 | KBI-10704 |

Lung Cancer

Lung cancer remains the leading cause of cancer death, annually resulting in more than one million cases worldwide. About 1.3 million new cases are diagnosed each year and prognoses are poor. Non-small cell lung cancer (NSCLC), the most common form (~80%) of lung cancer, has a 5-year survival rate of approximately 15%, mainly of late-stage detection¹.

A personalized medicine approach for treatment of NSCLC is emerging. Promising results have been obtained with specific anaplastic lymphoma kinase or ALK inhibitors like Crizotinib (Xalkori®)² in patients carrying the fusion gene ALK-EML4³.

Literature:



1. World Health Organization (Cancer Fact sheet N°297 2009: <http://www.who.int/mediacentre/factsheets/fs297/en/>)
2. Kwak et al, J Clin Oncol., 27(26):4247-53
3. Koivunen et al, Clin Cancer Res, 2008, 14, 4275-4283

ON ALK (2p23) Break

Translocations of the ALK (anaplastic lymphoma kinase) gene at 2p23 have originally been associated with anaplastic lymphomas, B-cell lymphomas, neuroblastomas and myofibroblastic tumors. To date at least 21 translocation partners have been described, however 80% of the translocations involves the NPM1 gene (5q35). More recently ALK rearrangements have been described in non-small cell lung cancer (NSCLC) cases. Promising results have been obtained with specific anaplastic lymphoma kinase or ALK inhibitors like Crizotinib (Xalkori®) in patients carrying the fusion gene ALK-EML4.

The ALK (2p23) Break probe is optimized to detect translocations involving the ALK gene region at 2p23.

Cat# KBI-10747 ALK (2p23) Break



ALK (2p23) Break probe hybridized to lung adenocarcinoma tissue showing translocation involving the ALK region at 2p23 (1R1G1G). Image kindly provided by Prof. B. Terris, Dr. P.A. Just, Hôpital Cochin, Paris.

Literature:

- Soda et al., Nature, 2007, 448, 561-566.
 Kwak et al, J Clin Oncol., 27(26):4247-53.
 Koivunen et al. Clin Cancer Res, 2008, 14, 4275-4283.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON ALK (2p23) Break | red/green | 10 | KBI-10747 |



ON ALK/EML4 t(2;2); inv(2) Fusion

The inversion in 2p21 and 2p23 leading to a fusion of the kinase domain of ALK (anaplastic lymphoma kinase) and EML4 (echinoderm microtubule associated protein like 4) has been described in 5-7% of non-small cell lung cancer (NSCLC) cases. ALK and EML4 are ~12 MB apart in opposite directions; a simple inversion generates the fusion gene.

Promising results have been obtained with specific anaplastic lymphoma kinase or ALK inhibitors like Crizotinib (Xalkori®) in patients carrying the fusion gene ALK-EML4.

The ALK/EML4 t(2;2); inv(2) Fusion probe is designed as a dual-color assay to detect the fusion of the ALK gene with the EML4 gene by paracentric inversion with breakage and reunion occurring at bands 2p21 and 2p23.

Cat# KBI-10746 ALK/EML4 t(2;2); inv(2) Fusion



ALK/EML4 t(2;2); inv(2) Fusion probe hybridized to lung adenocarcinoma tissue showing ALK-EML4 fusion (2RG1R1G).
Image kindly provided by Prof. B. Terris, Dr. P.A. Just, Hôpital Cochin, Paris.

Literature:

Soda et al., Nature, 2007, 448, 561-566.
Koivunen et al. Clin Cancer Res, 2008, 14, 4275-4283.



| Ordering information | Color | Tests | Cat# |
|-----------------------------------|-----------|-------|-----------|
| ON ALK/EML4 t(2;2); inv(2) Fusion | red/green | 10 | KBI-10746 |

ON EGFR, Her-1 (7p11) / SE 7

Epidermal growth factor receptor (EGFR) is a cell membrane protein, providing signal transduction and cell growth. It is a member of the Her or Erb-B family of type I receptor tyrosine kinases and implicated in the development and progression of non-small cell lung carcinomas (NSCLC), breast, intestine, and other organs. EGFR has been found to act as a strong prognostic indicator in head and neck, ovarian, cervical, bladder and oesophageal cancers. In these cancers, increased EGFR expression was associated with reduced recurrence-free or overall survival.

The EGFR (7p11) specific DNA Probe is optimized to detect copy numbers of the EGFR (Her-1) gene region at region 7p11. The chromosome 7 satellite enumeration (SE) probe is included to facilitate chromosome identification.

Cat.# KBI-10702 EGFR, Her1 (7p11) / SE 7



EGFR, Her1 (7p11) / SE 7 hybridized to colon tissue (2R2G).

Literature:

Wang et al, 1993, Jpn J Hum Genet, 38: 399-406.
Nicholson et al, 2001, Eur J Cancer, 37: 9-15.



| Ordering information | Color | Tests | Cat# |
|------------------------------|-----------|-------|-----------|
| ON EGFR, Her-1 (7p11) / SE 7 | red/green | 10 | KBI-10702 |

ON hTERT (5p15) / 5q31 (tissue)

Amplification of the hTERT gene at band 5p15 has been observed in a variety of cancer, particularly lung cancer, cervical tumors, and breast carcinomas. Several studies have revealed a high frequency of hTERT gene amplification in human tumors, which indicates that the hTERT gene may be a target for amplification during the transformation of human malignancies and that this genetic event probably contributes to a dysregulation of hTERT/ telomerase occurring in a subset of human tumors.

The hTERT (5p15) probe is designed as a dual-color assay to detect amplification at 5p15. The CDC25C/EGR1 (5q31) gene region probe is included as internal control.

Cat.# KBI-10709 hTERT (5p15) / 5q31 (tissue)



hTERT (5p15) / 5q31 (tissue) probe hybridized to paraffine embedded tissue (2R2G).

Literature:

Bryce et al, 2000, Neoplasia, 2; 197-201.
Zhang et al, 2000, Cancer Res, 60; 6230-6235.



| Ordering information | Color | Tests | Cat# |
|---------------------------------|-----------|-------|-----------|
| ON hTERT (5p15) / 5q31 (tissue) | red/green | 10 | KBI-10709 |

ON C-MET (7q31) / SE 7

The C-MET proto-oncogene is a receptor-like tyrosine kinase that drives a physiological cellular program important for development, cell movement, cell repair, cellular growth. Aberrant execution of this program has been associated to neoplastic transformation, invasion and metastasis. Activation of C-MET has been reported in a significant percentage of human cancers including non-small cell lung cancer (NSCLC). And is amplified during the transition between primary tumors and metastasis.

The C-MET (7q31) specific DNA probe is optimized to detect copy numbers of the C-MET gene region at region 7q31. The Chromosome 7 Satellite enumeration probe (SE 7) at D7Z1 is included to facilitate chromosome identification.

Cat.# KBI-10719 C-MET (7q31) / SE 7



Hybridization of MET Amplification Probe (KBI-10719) to a tissue section showing MET amplification.

Literature:

Go et al, 2010, J Thorac Oncol 5: 305-313.
Hara et al, 1998, Lab Invest 78; 1143-1153.
Tsugawa et al., 1998, Oncology 55; 475-481.

| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON C-MET (7q31) / SE 7 | red/green | 10 | KBI-10719 |

Prostate Cancer

Prostate cancer is the most commonly diagnosed, nondermatological malignancy in men; causing death in about 1 out of 35 Western men. Prostate cancer is the second leading cause of cancer death behind lung cancer.¹

Lately, great advances have been made using genomic technologies to develop predictive models that anticipate the risk of developing prostate cancer, prostate cancer progression, and the response of prostate cancer to therapy.²

Literature:



1. American Cancer Society: <http://www.cancer.org/cancer/prostatecancer/detailedguide/prostate-cancer-key-statistics>, revision of 11/22/2010, accessed on Dec 13, 2010.
2. Febbo, P.G., 2009, Cancer, 115: 3046-3057.

ON AR (Xq12) / SE X

The androgen receptor (AR) gene has been identified as a target gene for the Xq12 amplification found in one-third of hormone-refractory prostate cancers. The findings suggest that the AR gene amplification and overexpression is involved in the emergence of prostate cancer.

The AR (Xq12) specific DNA Probe is optimized to detect copy numbers of the AR gene region at region Xq12. The chromosome X satellite enumeration probe (SE X) at DXZ1 is included to facilitate chromosome identification.

Cat# KBI-10720 AR (Xq12) / SE X



AR (Xq12) / SE X probe hybridized to VCaP prostate cancer cell showing highlevel AR gene amplification.

Image kindly provided by Prof. Trapman, Erasmus Medical Centre, Rotterdam.

Literature:

- Visakorpi T et al, 1995, Nat. Genet. 9; 401-406.
Koivisto P et al, 1997, Cancer Res. 57 ; 314-319.



| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON AR (Xq12) / SE X | red/green | 10 | KBI-10720 |

ON PTEN (10q23) / SE 10

The gene 'phosphatase and tensin homolog deleted on chromosome 10' (PTEN), is a tumor suppressor located at chromosome 10q23, that plays an essential role in the maintenance of chromosomal stability, cell survival and proliferation. Loss of PTEN has been found in a wide number of tumors, and his important role is demonstrated by the fact that it is the second most frequently mutated gene after p53. Loss of PTEN significantly correlates with the advanced forms of gliomas, but also of prostate cancer and breast tumors.

The PTEN (10q23) specific DNA probe is optimized to detect copy numbers of the PTEN gene region at region 10q23. The Chromosome 10 Satellite enumeration probe (SE 10) at D10Z1 is included to facilitate chromosome identification.

Cat.# KBI-10718 PTEN (10q23) / SE 10



PTEN (10q23) / SE 10 probe hybridized to prostate cancer material showing deletion of PTEN gene region at 10q23 (1R2G).

Image kindly provided by Portuguese Cancer Inst., Porto.

Literature:

Cairns et al, 1997, Cancer Res, 57 ; 4997-5000.
Hermans et al, 2004, Genes Chrom Cancer, 39; 171-184.



| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON PTEN (10q23) / SE 10 | red/green | 10 | KBI-10718 |

ON TMPRSS2-ERG (21q22) Del, Break, TC

The transmembrane protease serine 2 gene (TMPRSS2) is involved in gene fusions with ERG, ETV1 or ETV4 in prostate cancer. In recent studies it has been reported that the expression of the TMPRSS2-ERG fusion gene is a strong prognostic factor for the risk of prostate cancer recurrence in prostate cancer patients treated by surgery.

The TMPRSS2-ERG rearrangement probe is optimized to detect the deletion between TMPRSS2 and ERG at 21q22 associated with the TMPRSS2-ERG fusion in a triple-color deletion assay. It also detects translocations involving the TMPRSS2 region such as ETV1 t(7;21), or ETV4 t(17;21).

Cat.# KBI-10726 TMPRSS2-ERG (21q22) Del, Break, TC



TMPRSS2-ERG (21q22) rearrangement probe hybridized to prostate carcinoma tissue showing a deletion of the TMPRSS2 (21q22) gene region associated with TMPRSS2-ERG fusion (1RGB 1RB).

Image kindly provided by Dr. Teixeira, Porto.

Literature:

Perner et al, 2006 Cancer Res 66; 8337-8341.
Hermans et al, 2006, Cancer Res 66; 10658-10663.
Attard et al, 2008, Oncogene 27; 253-263.

| Ordering information | Color | Tests | Cat# |
|------------------------------------|-----------|-------|-----------|
| TMPRSS2-ERG (21q22) Del, Break, TC | red/green | 10 | KBI-10726 |

Thyroid Carcinoma


Papillary thyroid carcinoma (PTC) is the most frequent primary carcinoma of the thyroid gland. PTC, conversely, is multifocal and are associated with prior radiation and high iodine intake. The follicular carcinomas is associated with endemic goiter and a diet with low iodine intake.

ON PPAR γ (3p25), Break

Follicular thyroid carcinoma is associated with the chromosomal translocation t(2;3)(q13;p25), fusing PAX8 (2q13) with the nuclear receptor, peroxisome proliferator-activated receptor γ (PPAR γ). PPAR is located in a breakpoint hot spot region, leading to recurrent alterations of this gene in thyroid tumors of follicular origin including carcinomas as well as adenomas with or without involvement of PAX8.

A break or split probe for PPAR γ is best used to analyze translocation of the PPAR γ (3p25) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat.# KBI-10707 PPAR γ (3p25) Break



PPAR γ (3p25) Break probe hybridized to patient material showing a translocation at 3p25 (1RG1R1G).
Image kindly provided by Dr. Valent, Paris.

Literature:
French et al, 2003, Am J Pathol, 162; 1053-1060.
Drieschner et l, 2006, Thyroid, 16; 1091-1096.

| Ordering information | Color | Tests | Cat# |
|--------------------------------|-----------|-------|-----------|
| ON PPAR γ (3p25), Break | red/green | 10 | KBI-10707 |

Neuroblastoma

According to the International Neuroblastoma Risk Grouping (INRG) Biology Committee MYCN remains the only genomic factor to be used currently for treatment stratification. Common data elements to be obtained by all groups include tumor cell ploidy and copy number/ LOH status at chromosome bands 1p36, 11q23, and 17q23-25.


Literature:
Ambros et al, 2006, Advances in Neuroblastoma Research.

ON MYCN (2p24) / LAF (2q11)

Amplification of the human N-myc protooncogene, MYCN, is frequently seen either in extrachromosomal double minutes or in homogeneously staining regions of aggressively growing neuroblastomas. MYCN amplification has been defined by the INRG as > 4-fold MYCN signals compared to 2q reference probe signals.

The MYCN (2p24) specific DNA probe is optimized to detect copy numbers of the MYCN gene region at 2p24. The LAF gene region probe at 2q11 is included to facilitate chromosome identification.

Cat.# KBI-10706 MYCN (2p24) / LAF (2q11)



MYCN (2p24) / LAF (2q11) hybridized to a cell line showing amplification of MYCN on chromosome 13 and 15.
Image kindly provided by Pasteur Workshop 2008, Paris.

Literature:
Shapiro et al, 1993, Am J Pathol, 142: 1339-1346.
Corvi et al, 1994, PNAS, 91: 5523-5527.

| Ordering information | Color | Tests | Cat# |
|-----------------------------|-----------|-------|-----------|
| ON MYCN (2p24) / LAF (2q11) | red/green | 10 | KBI-10706 |



ONCOLOGY - SOLID TUMOR DNA PROBES

ON SRD 1p36 / SE 1(1qh)

Neuroblastomas frequently have deletions of chromosome 1p and amplification of the N-myc oncogene. These deletions tend to be large and extend to the telomere, but a common region within sub-band 1p36.3 is consistently lost in these deletions. Inactivation of a tumor suppressor gene within 1p36.3 is believed to be associated with an increased risk for disease relapse. The 1p36 specific DNA probe has recently been changed to cover the recently described smallest region of consistent deletion (**SRD**) between D1S2795 and D1S253.

The SRD (1p36) specific DNA probe is optimized to detect copy numbers of the 1p36 region on chromosome 1. The chromosome 1 satellite enumeration probe (SE 1) at 1qh is included to facilitate chromosome identification.

Cat.# KBI-10712 SRD (1p36) / SE 1(1qh)



SRD (1p36) / SE 1 probe hybridized to a normal metaphase (2R2G).

Literature:

Caron et al, 1993, Nat Genet, 4: 187-190.
Cheng et al, 1995, Oncogene, 10: 291-297.
White et al, 2005, Oncogene, 24: 2684-2694.



| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON SRD (1p36) / SE1(1qh) | red/green | 10 | KBI-10712 |

ON MLL (11q23) / SE 11

Deletions of the long arm of chromosome 11 (11q) have been noted in primary neuroblastomas. It is assumed that a tumor suppressor gene mapping within 11q23.3 is commonly inactivated during the malignant evolution of a large subset of neuroblastomas, especially those with unamplified MYCN.

The MLL (11q23) specific DNA probe is optimized to detect amplification or deletion involving the MLL gene region at 11q23 in a dual-color assay on metaphase/interphase spreads, blood smears and bone marrow cells. The Chromosome 11 Satellite Enumeration probe (SE 11) at D11Z1 is included to facilitate chromosome identification.

Cat.# KBI-10711 MLL (11q23) / SE 11



MLL (11q23) / SE 11 hybridized to normal interphases (2R2G).

Literature:

Guo et al, 1999, Oncogene, 18: 4948-4957.
Maris et al, 2001, Med Pediatr Oncol, 36: 24-27.


| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON MLL (11q23) / SE11 | red/green | 10 | KBI-10711 |

ON MDM4 (1q32) / SE1

MDM4 (also known as MDMX, murine double minute gene) is a relative of MDM2 that was identified on the basis of its ability to physically interact with p53. MDM4, like MDM2, acts as a key negative suppressor of p53 by interfering with its transcriptional activity. MDM4 amplification and/or overexpression occurs in several diverse tumors. Studies showed an increased MDM4 copy number in 65% of human retinoblastomas compared to other tumors, qualifying MDM4 as a specific chemotherapeutic target for treatment of this tumor.

The MDM4 (1q32) specific DNA probe is designed as a dual-color assay to detect amplification at 1q32. The chromosome 1 Satellite Enumeration (SE 1) probe at 1qh is included to facilitate chromosome identification.

Cat.# KBI-10736 MDM4 (1q32) / SE 1



MDM4 (1q32) / SE 1 probe hybridized to paraffin embedded tissue (2R2G).

Literature:

Riemenschneider et al, 1999, Cancer Res. 59 ; 6091-6096.
Danovi et al, 2004, Mol.Cell.Bio. 24; 5835-5843.

| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON MDM4 (1q32) / SE 1 | red/green | 10 | KBI-10736 |

Sarcoma



Sarcoma is a general class of less common cancers in which the cancer cells arise from or resemble normal cells in the body known as “connective tissues” (fat, muscle, blood vessels, deep skin tissues, nerves, bones, and cartilage). The benign and malignant forms have related karyotypic changes which provide an important resource for identifying the additional genetic changes that occur in the malignant compared with the benign form. In fact, the molecular biology of soft-tissue sarcomas has provided the perfect example of how cytogenetic and molecular approaches can contribute toward a clearer understanding of the development of soft-tissue sarcomas.

ON EWSR1 (22q12) Break

Ewing’s sarcoma is the second most frequent primary bone cancer. In most cases a translocation involving the EWSR1 gene at 22q12 and the FLI1 gene at 11q24 are observed, but several other translocation partners (ERG, ETV1, FEV, and E1A3) can also be involved.

The EWSR1 (22q12) Break Probe is optimized to detect translocations involving the EWSR1 gene region at 22q12 in a dual-color, split assay on metaphase/interphase spreads and paraffin embedded tissue sections.

Cat.# KBI-10750 EWSR1 (22q12) Break



Interphase FISH result of the EWSR1 split probe. Co-localisation indicates intact EWSR1 locus, separation of the green and red signal indicates the presence of a translocation breakpoint.

Literature:
Zucman-Rossi, et al, 1998, PNAS, 95; 11786-11791.
Bernstein et al, 2006, Oncologist, 11; 503-519.

| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON EWSR1 (22q12) Break | red/green | 10 | KBI-10750 |

ON SYT (18q11) Break

The characteristic chromosomal abnormality in synovial sarcoma is t(X;18)(p11.2;q11.2) present in 90% of patients. This translocation results in the fusion of the chromosome 18 SYT gene to either of two distinct genes, SSX1 or SSX2, located on the X chromosome.



A break or split probe for SYT is best used to analyze translocation of the SYT (SS18) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

ON CHOP (12q13) Break

Liposarcoma is one of the most frequent sarcomas in adults, representing 10 to 16 percent of soft tissue sarcomas. Most patients with round cell/myxoid liposarcoma have an acquired t(12;16)(CHOP-FUS) or t(12;22)(CHOP-EWS) translocation, both of which involve the CHOP gene at 12q13. A break or split probe for CHOP is best used to analyze translocation of the CHOP (12q13) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.



The CHOP (12q13) Break probe is optimized to detect translocations involving the CHOP gene region at 12q13 in a dual-color, split assay.

Cat.# KBI-10713 SYT (18q11) Break



SYT (18q11) Break probe hybridized to patient material showing translocation of the SYT (SS18) gene region at 18q11 (1RG1R1G).

Cat.# KBI-10714 CHOP (12q13) Break



CHOP (12q13) Break probe hybridized to a normal metaphase (2RG).

Literature:

Kawai et al, 1998, NEJM, 338; 153-160.
Surace et al, 2004, LabInvest., 84; 1185-1192.

Literature:

Panagopoulos et al, 1994, Cancer Res, 54; 6500-6503.
Schoenmakers et al, 1994, Genomics, 20; 210-222.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON SYT (18q11) Break | red/green | 10 | KBI-10713 |

| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON CHOP (12q13) Break | red/green | 10 | KBI-10714 |



ON FUS (16p11) Break

The FUS gene was originally shown to be rearranged in myxoid liposarcomas harboring a t(12;16)(q13;p11) translocation.

FUS has also been shown to be involved in other recombinations: with ERG in acute myeloid leukemia carrying a t(16;21), with ATF1 in band 12q13 in angiomatoid fibrous histiocytoma, and with CREB3L2 in fibromyxoid sarcoma.

A break or split probe for FUS is best used to analyze translocation of the FUS (16p11) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat.# KBI-10715 FUS (16p11)



FUS (16p11) Break probe hybridized to liposarcoma material.

Literature:

Shing et al, 2003, Cancer Res, 63: 4568-4576.

Storlazzi et al, 2003, Hum. Mol. Genet., 12: 2349-2358.

| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON FUS (16p11), Break | red/green | 10 | KBI-10715 |



ON FKHR (13q14) Break

The t(2;13) is associated with alveolar rhabdomyosarcomas.

This translocation results in the formation of a chimeric transcript consisting of the 5' portion of PAX3, including an intact DNA-binding domain fused to the FKHR gene on chromosome 13. The t(1;13)(p36;q14) also seen in alveolar rhabdomyosarcomas results in the fusion of another member of the PAX family, PAX7 to the FKHR gene on chromosome 13.

A break or split probe for FKHR is best used to analyze translocation of the FKHR (13q14) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat.# KBI-10716 FKHR (13q14) Break



FKHR (13q14) Break probe hybridized to patient material.

Literature:

Barr et al, 1996, Hum. Mol. Genet., 5: 15-21.

Coignet et al, 1999, Genes Chrom. Cancer, 25: 222-229.

| Ordering information | Color | Tests | Cat# |
|-----------------------|-----------|-------|-----------|
| ON FKHR (13q14) Break | red/green | 10 | KBI-10716 |

ON MDM2 (12q15) / SE 12



Well-differentiated liposarcoma/atypical lipomatous tumor and dedifferentiated liposarcoma are among the most common malignant soft tissue tumors presented in older adults. These tumors can be difficult to distinguish from benign lipomatous neoplasms and other high-grade sarcomas.

Amplification of the MDM2 gene has been identified in lipomatous neoplasms. The use of fluorescence in situ hybridization in identifying MDM2 amplification has made the MDM2 amplification probe a valuable diagnostic tool in well-differentiated liposarcomas/atypical lipomatous tumors.

Fibrosarcoma is a rare soft-tissue tumor composed of fascicles of spindled fibroblast-like cells. Gains and high-level amplifications of 12q14–22 were the most common genomic imbalances, and reflected MDM2 amplification, thereby indicating the importance of this gene in the evolution of fibrosarcomas.

The MDM2 (12q15) specific DNA probe is optimized to detect copy numbers of the MDM2 gene region at region 12q15. The Chromosome 12 Satellite Enumeration probe (SE 12) at D12Z3 is included to facilitate chromosome identification.

Cat.# KBI-10717 MDM2 (12q15) / SE 12



MDM2 (12q15) / SE 12 Amplification probe hybridized to patient material showing amplification of the MDM2 gene region at 12q15.

Literature:

Uchida et al., 2010, Cancer Genet Cytogenet 203; 324-327.
Lucas et al, 2010, Am J Surg Pathol 34: 844-851.
Weaver et al, 2008, Mod Pathol 21: 943-949.
Mitchell et al, 1995, Chrom. Res., 3; 261-262.
Reifenberger et al, 1996, Cancer Res., 15; 5141-5145.



| Ordering information | Color | Tests | Cat# |
|------------------------|-----------|-------|-----------|
| ON MDM2 (12q15) / SE12 | red/green | 10 | KBI-10717 |

ON CDK4 (12q13) / SE 12

Amplification of the CDK4 gene region at 12q13-q15 has been observed in several types of cancer, especially in gliomas and sarcomas. CDK4 codes for a cyclin dependent kinase which is involved in controlling progression through the G1 phase of the cell cycle. The oncogenic potential of CDK4 activation has been related to the deregulation of the G1 phase by increasing the hyperphosphorylation of retinoblastoma tumor suppressor protein helping to cancel its growth-inhibitory effects.

The CDK4 (12q13) specific DNA probe is optimized to detect copy numbers of the CDK4 gene region at 12q13. The chromosome 12 satellite enumeration probe (SE 12) at D12Z3 is included to facilitate chromosome identification.

Cat.# KBI-10725 CDK4 (12q13) / SE 12



CDK4 (12q13) / SE 12 probe hybridized to liposarcoma tissue showing multiple amplification involving the CDK4 gene region at 12q13 (3+R2G). Image kindly provided by Dr. Sapi, Hungary.

Literature:

Kuhnen et al, 2002, Virchows Arch 441 ; 299-302.
Shimada et al, 2006, Hum Path 37(9) ; 1123-1129.

| Ordering information | Color | Tests | Cat# |
|-------------------------|-----------|-------|-----------|
| ON CDK4 (12q13) / SE 12 | red/green | 10 | KBI-10725 |



Different Cancer types

ON ERCC1 (19q13) & ZNF443 (19p13)

Nucleotide excision repair (NER) is the primary DNA repair mechanism that removes platinum-DNA adducts from genomic DNA. Excision repair cross-complementing rodent repair deficiency, complementation group 1 (ERCC1) is a critical gene in the NER pathway. A growing list of reports links cisplatin, carboplatin, and oxaliplatin based chemotherapy resistance to ERCC1 expression levels in several tumors. This relationship has been suggested for patients with gastric, bladder, ovarian, colorectal and non-small-cell lung cancers (NSCLC). ERCC1 has been shown to be an important marker to predict responsiveness to cisplatin-based chemotherapy. Low ERCC1 gene expression correlates with prolonged survival after cisplatin-based chemotherapy.

The ERCC1 (19q13) specific DNA Probe has been optimized to detect copy numbers of the ERCC1 gene region at 19q13. The ZNF443 (19p13) probe is included to facilitate chromosome identification.

Cat.# KBI-10739 ERCC1 (19q13) & ZNF443 (19p13)



ERCC1 (19q13) & ZNF443 (19p13) probe hybridized to paraffin embedded tissue (2R2G).

Literature:

Olaussen et al, 2006, N. Engl. J. Med. 335; 983-991.
Ceppi et al, 2006, Ann. Oncol. 17; 1818-1825.

| Ordering information | Color | Tests | Cat# |
|-----------------------------------|-----------|-------|-----------|
| ON ERCC1 (19q13) / ZNF443 (19p13) | red/green | 10 | KBI-10739 |

ON AURKA (20q13) / 20q11



Aurora kinase A (AURKA) has the fundamental role of regulating proper centrosome function, important to maintain genomic stability during cell division and to ensure equal segregation of replicated chromosomes to daughter cells.

Deregulated duplication and distribution of centrosomes has been implicated in mechanisms leading to mitotic spindle aberrations, aneuploidy, and genomic instability that are seen in many different tumor types. Consistent with this, AURKA amplification has been detected in approximately 12% of primary breast tumors, as well as in breast, bladder, ovarian, colon, prostate, neuroblastoma and cervical cancer cell lines.

Recent investigations on new drugs developments have focused on the importance of aurora kinases for tumor suppression.

The AURKA (20q13) specific DNA probe is optimized to detect copy numbers of the AURKA gene region at region 20q13. The 20q11 specific DNA probe is included to facilitate chromosome identification.

Cat.# KBI-10721 AURKA (20q13) / 20q11



AURKA (20q13) / 20q11 probe hybridized to colorectal carcinoma material showing amplification of AURKA, gene region at 20q13.
Material kindly provided by Dr. Carvalho, Amsterdam.

Literature:

Uchida et al., 2010, Cancer Genet Cytogenet 203; 324-327.
Sen et al, 2002, J of Nat Canc Inst 94; 1320-1329.
Lassmann et al, 2007, Clin Cancer Res 13; 4083-4091.

| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON AURKA (20q13) / 20q11 | red/green | 10 | KBI-10721 |

ON AURKB (17p13) / SE 17



Aurora kinase B (AURKB) localizes to microtubules, and is a key regulator of the mitotic cell division and chromosome segregation processes. Gain of function of AURKB correlates with cell proliferation, induction of multinuclear cells, and chromosomal instability.

The significant interest of the gene in cancer diagnostics is related to the driving function of AURKB in tumor progression, histological differentiation, and metastasis. AURKB is predictive for the aggressive recurrence of many different types of tumors, including hepatocellular carcinoma and oral squamous cell carcinoma.

Recently new drugs have been under investigation for their capacity of interfering with the aurora kinases activity related to tumor-suppressor effects.

The AURKB (17p13) specific DNA probe is optimized to detect copy numbers of the AURKB gene region at region 17p13. The Chromosome 17 Satellite Enumeration (SE 17) probe at D17Z1 is included to facilitate chromosome identification.

Cat.# KBI-10722 AURKB (17p13) / SE 17



AURKB (17p13) / SE 17 probe hybridized to tumor tissue (2R2G).

Literature:
Smith et al, 2005, Br J Cancer, 93; 719-729.

| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON AURKB (17p13) / SE 17 | red/green | 10 | KBI-10722 |

ON CCND1 (11q13) / SE 11



CCND1 (also named Cyclin D1 or BCL1) is a key cell cycle regulator of the G1 to S phase progression. The binding of cyclin D1 to cyclin-dependent kinase (CDKs) leads to the phosphorylation of retinoblastoma protein (pRb), subsequently triggering the release of E2F transcription factors to allow G1 to S phase progression of the cell cycle.

Consistent with this function, overexpression of cyclin D1 results in a more rapid progression from the G1 to S phase transition and in a reduced serum dependency in fibroblast cells, characteristics typically seen in cancer cells.

Amplification of cyclin D1 plays pivotal roles in the development of a subset of human cancers including parathyroid adenoma, breast cancer, colon cancer, lymphoma, melanoma, and prostate cancer.

The CCND1 (11q13) specific DNA Probe is optimized to detect copy numbers of the CCND1 gene region at region 11q13. The Chromosome 11 Satellite Enumeration (SE 11) probe at D11Z1 is included to facilitate chromosome identification.

Cat.# KBI-10734 CCND1 (11q13) / SE 11



CCDN1 (11q13) / SE 11 probe hybridized to patient interphases/ metaphase showing CCDN1 (11q13) amplification with polyploidy for chromosome 11.

Literature:
Okami et al, 1999, Oncogene 18; 3541-3545.
Freier et al, 2003, Cancer Res; 1179-1182.



| Ordering information | Color | Tests | Cat# |
|--------------------------|-----------|-------|-----------|
| ON CCND1 (11q13) / SE 11 | red/green | 10 | KBI-10734 |

ON IGH (14q32) Break (tissue)

Chromosomal rearrangements involving the immunoglobulin heavy chain gene (IGH) at 14q32 are observed in 50% of patients with B-cell non-Hodgkin's lymphoma (NHL) and many other types of Lymphomas. More than 50 translocation partners with IGH have been described. In particular t(8;14), is associated with Burkitt's lymphoma, t(11;14) is associated with Mantle cell lymphoma, t(14;18) is observed in a high proportion of follicular lymphomas and t(3;14) is associated with Diffuse Large B-Cell Lymphoma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14q32 in a dual-color, split assay. Kreatech has developed this probe for the specific use on cell material (KBI-10601), or for the use on tissue (KBI-10729).

Cat.# KBI-10729 IGH (14q32) Break (tissue)



IGH (14q32) Break probe hybridized to patient material showing a partial deletion of 14q32 (1R1R).

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969.
Gozetti et al, 2002, Cancer Research, 62: 5523-5527.



| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON IGH (14q32) Break (tissue) | red/green | 10 | KBI-10729 |

ON BCL2 (18q21) Break (tissue)

Follicular lymphoma is a mature B-cell lymphoma characterized by the presence of the t(14;18) translocation that juxtaposes the BCL2 locus on chromosome 18q21 to the immunoglobulin H (IGH) locus on chromosome 14q32, resulting in the overexpression of the anti-apoptotic protein BCL2. Next to IGH, additional translocation partners to BCL2 have been identified (e.g. IGK at 2p11.2 and IGL at 22q11). A break or split assay is therefore best suited to detect rearrangements of the BCL2 gene region at 18q21.

The BCL2 (18q21) Break probe is optimized to detect translocations involving the BCL2 gene region at 18q21 in a dual-color, split assay on paraffin embedded tissue sections. Kreatech has developed this probe for the specific use on cell material (KBI-10612), or for the use on tissue (KBI-10745).

Cat# KBI-10745 BCL2 (18q21) Break (tissue)



ON BCL2 (18q21) Break hybridized to paraffin embedded tissue (2RG).

Literature:

Taniwaki M et al, 1995, Blood, 86; 1481-1486.
Poetsch M et al, 1996, J Clin Oncol, 14; 963- 969.
Einerson R et al, 2005, Am J Clin Pathol, 124; 421-429.



| Ordering information | Color | Tests | Cat# |
|--------------------------------|-----------|-------|-----------|
| ON BCL2 (18q21) Break (tissue) | red/green | 10 | KBI-10745 |

ON BCL6 (3q27) Break (tissue)

Chromosomal translocations involving band 3q27 with various different partner chromosomes represent a recurrent cytogenetic abnormality in B-cell non-Hodgkin's lymphoma. A FISH strategy using two differently labeled flanking BCL6 probes provides a robust, sensitive, and reproducible method for the detection of common and uncommon abnormalities of BCL6 gene in interphase nuclei. Kreatech has developed this probe for the specific use on cell material (KBI-10607), or for the use on tissue (KBI-10730).

Two different breakpoint regions have been identified; the major breakpoint region (**MBR**) is located within the 5' noncoding region of the BCL6 proto-oncogene, while the atypical breakpoint region (**ABR**) is located approximately 200 kb distal to the BCL6 gene. The BCL6 (3q27) Break probe is designed to flank both possible breakpoints, thereby providing clear split signals in either case.

Cat.# KBI-10730 BCL6 (3q27) Break (tissue)



BCL6 (3q27) Break probe hybridized to patient material showing both normal (2RG) and aberrant signals (1RG1R1G).

Image kindly provided by Prof Siebert, Kiel.

Literature:

Butler et al, 2002, Cancer Res, 62; 4089-4094.
Sanchez-Izquierdo, 2001, Leukemia, 15; 1475-1484.



| Ordering information | Color | Tests | Cat# |
|-------------------------------|-----------|-------|-----------|
| ON BCL6 (3q27) Break (tissue) | red/green | 10 | KBI-10730 |

ON MALT (18q21) Break (tissue)

Low grade malignant lymphomas arising from mucosa associated lymphoid tissue (MALT) represent a distinct clinicopathological entity. The three major translocations seen in MALT lymphomas are t(11;18)(q21;q21)/API2-MALT1, t(14;18)(q32;q21)/IGH-MALT1 and t(1;14)(p22;q32)/IGH-BCL10. A break or split probe for MALT (18q21) is best used to analyze translocation of the MALT gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Kreatech has developed this probe for the specific use on cell material (KBI-10608), or for the use on tissue (KBI-10731).

Cat.# KBI-10731 MALT (18q21) Break (tissue)



MALT (18q21) Break tissue probe hybridized to paraffin embedded material (2RG).

Literature:

Morgan et al, 1999, Cancer Res, 59; 6205-6213.
Dierlamm et al, 2000, Blood, 96; 2215-2218.



| Ordering information | Color | Tests | Cat# |
|--------------------------------|-----------|-------|-----------|
| ON MALT (18q21) Break (tissue) | red/green | 10 | KBI-10731 |

ON p53 (17p13) / SE 17 (tissue)

The p53 tumor suppressor gene at 17p13, has been shown to be implicated in the control of normal cellular proliferation, differentiation, and apoptosis. Allelic loss, usually by deletion, and inactivation of p53 have been reported in numerous tumor types and are associated with poor prognosis in CLL.

The p53 (17p13) specific DNA probe is optimized to detect copy numbers of the p53 gene region at 17p13. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification. Kreatech has developed this probe for the specific use on cell material (KBI-10112/KBI-12112), or for the use on tissue (KBI-10738).

Cat.# KBI-10738 p53 (17p13) / SE 17 (tissue)



P53 (17p13) / SE 17 (tissue) probe hybridized to paraffin embedded tissue (2R2G).

Literature:

Amiel A et al, 1997, Cancer Gener.Cytogenet., 97; 97-100.
Drach J et al, 1998, Blood, 92; 802-809.



| Ordering information | Color | Tests | Cat# |
|---------------------------------|-----------|-------|-----------|
| ON p53 (17p13) / SE 17 (tissue) | red/green | 10 | KBI-10738 |

ON TFE3 (Xp11) Break

Abnormalities of Xp11.2 region have often been observed in papillary renal cell carcinomas and are sometimes the sole cytogenetic abnormality present. The transcription factor binding to IGHM enhancer 3 (TFE3) gene, which encodes a member of the helix-loop-helix family of transcription factors, is located in this critical region and can be fused to various other chromosomal regions by translocation. Known fusion partners are NONO (Xq12), PRCC (1q21), SFPQ (1p34), CLTC (17q23) and ASPSCR1 (17q25).

The TFE3 (Xp11) Break probe is optimized to detect translocations involving the TFE3 gene region at Xp11.2 in a dual-color, split assay.

Cat.# KBI-10741 TFE3 (Xp11) Break



TFE3 (Xp11) Break probe hybridized to renal cell carcinoma showing a translocation at Xp11 (1R1G1G).
Image kindly provided by Dr. Desangles, Paris .

Literature:

Sidhar et al, 1996, Hum Mol Genet, 5; 1333-1338.
Weternan et al., 1996, Proc Natl, Acad Sci, 93; 15294-15298.


| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON TFE3 (Xp11) Break | red/green | 10 | KBI-10741 |

ON ALK (2p23) Break

Translocations of the ALK (anaplastic lymphoma kinase) gene at 2p23 have originally been associated with anaplastic lymphomas, B-cell lymphomas, neuroblastomas and myofibroblastic tumors. To date at least 21 translocation partners have been described, however 80% of the translocations involves the NPM1 gene (5q35). More recently ALK rearrangements have been described in non-small cell lung cancer (NSCLC) cases. Promising results have been obtained with specific anaplastic lymphoma kinase or ALK inhibitors like Crizotinib (Xalkori®) in patients carrying the fusion gene ALK-EML4.

The ALK (2p23) Break probe is optimized to detect translocations involving the ALK gene region at 2p23.

Cat# KBI-10747 ALK (2p23) Break



ALK (2p23) Break probe hybridized to lung adenocarcinoma tissue showing translocation involving the ALK region at 2p23 (1RG1R1G).
Image kindly provided by Prof. B. Terris, Dr. P.A. Just, Hôpital Cochin, Paris.

Literature:

Soda et al., Nature, 2007, 448, 561-566.
Kwak et al, J Clin Oncol., 27(26):4247-53.
Koivunen et al. Clin Cancer Res, 2008, 14, 4275-4283.

| Ordering information | Color | Tests | Cat# |
|----------------------|-----------|-------|-----------|
| ON ALK (2p23) Break | red/green | 10 | KBI-10747 |

ONCOLOGY – CHROMOGENIC IN SITU HYBRIDIZATION


Chromogenic *In Situ* Hybridization (CISH) is increasingly emerging as a viable alternative to FISH and often selected as the method of choice for molecular pathologists for visualizing over-expression of genes involved in tumor development. CISH, like FISH, directly visualizes the number of gene copies present in the nucleus, and it produces a permanent record of the slide that can be interpreted with a regular light microscope in the context of the tumor histopathology.

TwinStar

Our latest product range TwinStar is designed for dual-color CISH enabling the possibility to study ratio of genes in a light microscope, like for the Her-2 gene in relation to the centromeric region of chromosome 17.

Each kit includes the corresponding FISH probes and the specific TwinStar detection module converting both signals into chromogenic signals via a colorimetric assay system.

The TwinStar CISH Detection Kit is a module of all TwinStar assays including all reagents required to perform dual-color CISH. The detection part includes two proprietary substrates converting both fluorescent signals into distinct colorimetric signals, the red fluor is converted into a red signal and the green fluor into a grayish-green signal. TwinStar provides a universal solution capable of transforming all of Kreatech's dual-color FISH probes into chromogenic signals.




Amplified Her-2/Neu (ERBB2) on a breast cancer specimen visualized with the TwinStar dual color CISH kit.





UniStar

UniStar utilizes a novel assay based on KREATECH's well established POSEIDON™ DNA probes, enabling the clinician to consecutively perform FISH and CISH on the same sample.

In a first step, the sample is hybridized to specific fluorescent-labeled DNA probes enabling analysis under a fluorescence microscope. Subsequently, the fluorescent signal of the gene of interest can be converted into a chromogenic signal, which can then be analyzed with a regular bright field microscope. Alternatively, FISH can be omitted for direct CISH analysis.



Each of the UniStar kits include the corresponding POSEIDON™ DNA probes for copy number detection of the gene of interest, as well as a control probe for performing FISH in a dual-color assay. In addition, they include a specific detection module converting the signal of the critical probe into a chromogenic single-color signal via a colorimetric assay system.



Example of FISH and CISH performed consecutively on the same slide: breast cancer tissue hybridized with the Her-2/Neu POSEIDON™ fluorescent probe followed by conversion of the red Her-2/Neu signal into a brownish chromogenic signal using the UniStar CISH Detection Kit.

Glioblastoma specimen showing amplified EGFR. The slide was hybridized with the EGFR POSEIDON™ fluorescent probe followed by conversion of the red EGFR signal into a colorimetric signal using the UniStar CISH Detection Kit. Image kindly provided by Dr. K. Beiske, Oslo University Hospital, Norway.

Product and ordering information


| Product | Description | Tests | Cat# |
|-----------------------------|---|-------|-----------|
| UniStar CISH Detection Kit | UniStar CISH Detection Kit for the use with POSEIDON™ DNA Probes labeled in red | 10 | KBI-50001 |
| UniStar Her2/neu (17q12) | DNA probes specific for Her-2/Neu and SE 17, UniStar CISH Detection Kit | 10 | KBI-50701 |
| UniStar EGFR (7p11) | DNA probes specific for EGFR and and SE 7, UniStar CISH Detection Kit | 10 | KBI-50702 |
| UniStar C-MET (7q31) | DNA probes specific for C-MET and and SE 7, UniStar CISH Detection Kit | 10 | KBI-50719 |
| TwinStar CISH Detection Kit | For the use with POSEIDON™ DNA probes labeled in red and green | 10 | KBI-60010 |
| TwinStar Her2/neu (17q12) | DNA probes specific for Her-2/Neu and SE 17, TwinStar CISH Detection Kit | 10 | KBI-60701 |
| TwinStar EGFR (7p11) | DNA probes specific for EGFR and and SE 7, TwinStar CISH Detection Kit | 10 | KBI-60702 |
| TwinStar C-MET (7q31) | DNA probes specific for C-MET and and SE 7, TwinStar CISH Detection Kit | 10 | KBI-60719 |

PREIMPLANTATION GENETIC SCREENING

FISH is the current gold standard to determine the chromosomal constitution of an embryo. In contrast to karyotyping it can be used on interphase chromosomes, so that it can be applied on polar bodies, blastomeres and other single cell samples. FISH is therefore accepted as a routine method in preimplantation genetic screening (PGS) in determining chromosome aneuploidies prior to implanting an embryo and increases the success rate of an IVF-mediated pregnancy.

MultiStar 24 FISH

MultiStar 24 FISH consists of four DNA probe mixes each hybridizing to six different chromosomes using six different fluorochromes. It can be applied in lymphocytes, sperm and blastomeres.



Metaphase preparation from male lymphocyte cells visualizing all 24 chromosomes.
Image kindly provided by Prof. D. Griffin, University of Kent, United Kingdom.

The fully optimized protocol consists of three 15-30 minutes hybridizations followed by a 6 to 16 hours hybridization. In between the individual hybridizations, the preceding probes are washed off after imaging the results. The morphology of the cell types is retained despite repeated denaturation, hybridization







and post hybridization washes. Moreover the entire protocol can be completed within 24 hours, which fits the window for clinical PGS application. This novel method eliminates the bottleneck perceived when using FISH by omitting key chromosomes relevant to implantation failure not covered in a limited panel.

The probes required to provide information for 24 chromosomes comprise four panels each containing 6 chromosome-specific identifier sequences each labelled with a different fluorochrome. The first 3 panels use centromeric sequences (Panel 1: chromosomes 1,3,4,6,7,8; Panel 2: chromosomes 9,10,11,12,17,20 and Panel 3: chromosomes 2,15,16,18,X,Y) . In contrast, Panel 4 used for the final round of hybridization uses unique sequence probes for chromosomes 5,13,14,19,21,22 since centromeric sequences were not available for these chromosomes. A combination of separate probe denaturation for 10 minutes followed by a short co-denaturation between probe and sample is used. Panels 1-3 require short hybridization times (15 minutes) whereas Panel 4 requires a 6 to 16 hours hybridization.

Setup of the different layers

All DNA probes are labeled with PlatinumBright™ based on the Universal Linkage System (ULS™), KREATECH's proprietary non-enzymatic labeling technology capable of linking fluorescent labels or haptens to any nucleic acid of interest.


Chromosomes recognized by the different panels:

| Color |  |  |  |  |  |  |
|---------------------|---|---|---|--|---|---|
| Label | PlatinumBright 405 | PlatinumBright 415 | PlatinumBright 495 | PlatinumBright 547 | PlatinumBright 590 | PlatinumBright 647 |
| Color | Dark Blue | Blue | Green | Light Red | Dark Red | Far Red |
| Excitation/Emission | 410/455 | 429/470 | 495/517 | 547/565 | 587/612 | 647/665 |
| MultiStar Panel 1 | 7 | 1 | 6 | 8 | 3 | 4 |
| MultiStar Panel 2 | 11 | 9 | 20 | 12 | 10 | 17 |
| MultiStar Panel 3 | 18 | Y | X | 16 | 2 | 15 |
| MultiStar Panel 4 | 19 | 5 | 21 | 22 | 13 | 14 |

PreimpScreen PolB (13,16,18,21,22)

PreimpScreen PolB is designed for determining chromosome copy number in polar bodies.

The first polar body is removed from the unfertilized oocyte, and the second polar body from the zygote, shortly after fertilization. The main advantage of the use of polar bodies in preimplantation genetic screening (PGS) is that they are not necessary for successful fertilization or normal embryonic development, thus ensuring no deleterious effect for the embryo. In some countries, where the legislation bans the selection of preimplantation embryos, polar body analysis is the only possible method to perform PGS. The biopsy and analysis of the first and second polar bodies can be completed before syngamy, which is the moment from which the zygote is considered an embryo and becomes protected by the law.




Pseudo color image using PreimpScreen PolB (KBI-40050) on a metaphase spread from lymphocytes showing two signals each of chromosomes 13, 16, 18, 21, and 22, respectively.

PreimpScreen Blas (13,18,21,X,Y)

PreimpScreen Blas is designed for determination of chromosome copy number in blastomeres.

Cleavage-stage biopsy is generally performed the morning of day three post-fertilization, when normally developing embryos reach the eight-cell stage. The biopsy is usually performed on embryos with less than 50% of anucleated fragments and at an 8-cell or later stage of development. The main advantage of cleavage-stage biopsy over polar body (PB) analysis is that the genetic input of both parents can be studied, and therefore currently is the prevalent method when doing in situ hybridizations in preimplantation genetic screening.



Pseudo-color image on a healthy female blastomer using PreimpScreen Blas (13,18,21,X,Y) FISH panel, KBI-40051. Image kindly provided by Prof. D. Griffin, University of Kent, United Kingdom.

Literature:

Iaonnou D et al, 2012, Chromosome Res, 20:447-60.
Iaonnou D et al, 2011, Mol and Cel Probes, 25:199-205.

Product and ordering information

| Product | Description | Tests | Cat# |
|------------------------------------|--|-------|-----------|
| PreimpScreen PolB (13,16,18,21,22) | Five-color FISH-mix consisting of DNA probes specific for chromosomes 13, 16, 18, 21, and 22 | 20 | KBI-40050 |
| PreimpScreen Blas (13,18,21,X,Y) | Five-color FISH-mix consisting of DNA probes specific for chromosomes 13, 18, 21, X, and Y | 20 | KBI-40051 |
| MultiStar 24 FISH | FISH probe panel for visualizing all 24 chromosomes (including the four panels KBI-40061, KBI-40062, KBI-40063, and KBI-40064) | 10 | KBI-40060 |
| MultiStar FISH Panel 1 | FISH panel of centromeric probes for chromosomes 1, 3, 4, 6, 7, and 8 | 10 | KBI-40061 |
| MultiStar FISH Panel 2 | FISH panel of centromeric probes for chromosomes 9, 10, 11, 12, 17, and 20 | 10 | KBI-40062 |
| MultiStar FISH Panel 3 | FISH panel of centromeric probes for chromosomes 2, 15, 16, 18, X, and Y | 10 | KBI-40063 |
| MultiStar FISH Panel 4 | FISH panel of unique sequence probes for chromosomes 5, 13, 14, 19, 21, and 22 | 10 | KBI-40064 |



PRENATAL DNA PROBES

Prenatal cytogenetic analysis requires the isolation of metaphase chromosomes and takes 7-14 days for the final results. This waiting period tends to cause psychological distress for pregnant women and their families. Aneuploidies of 5 chromosomes (13, 18, 21, X, Y) account for 95% of the chromosomal aberrations that cause infants born with defects. Fluorescent labeled DNA probes of the 13, 18, 21, X, Y chromosomes can be used on uncultured cells obtained directly from amniotic fluid. The FISH rapid technique allows to reliably detect numerical aberrations for these chromosomes. While not all chromosome abnormalities can be identified simply by counting specific chromosomes within a cell, the majority of the most common abnormalities of chromosome number, including Down syndrome (trisomy 21), trisomy 18, trisomy 13, Klinefelter syndrome (47,XXY), triple-X syndrome (47,XXX), Turner syndrome (45,X) and 47,XYY can be reliably determined. The FISH analysis does not detect structural chromosome abnormalities, mosaicism, and other numerical chromosome abnormalities (excluding X, Y, 13, 18, and 21). In addition, false-positive or negative results, as well as maternal cell contamination, have been demonstrated in prenatal FISH analysis. It is recommended (e.g. American College of Medical Genetics) that irreversible therapeutic action should not be initiated on the basis of FISH results alone.

Trisomy 21 – Down Syndrome

Down syndrome is caused by an extra chromosome 21. It is the most common single cause of human birth defects, with an occurrence in 1 out of every 660 births. Congenital heart defects are frequently present in Down syndrome children. The normal life span mainly is shortened in Down syndrome by congenital heart disease and by increased incidence of acute leukemia. Mental retardation is variable, and usually moderate. Some adults live independently and are accomplished individuals. The chromosome 21 specific region probe is optimized to detect copy numbers of chromosome 21 at 21q22.1 on uncultured amniotic cells. In all PN combinations the 21q specific DNA probe is direct-labeled in red with PlatinumBright550.

Cat.# KBI-40002 PN 21 (21q22)





PN 21 (21q22) probe hybridized to a normal metaphase (2R).

Trisomy 13 – Patau Syndrome

Trisomy 13, also called Patau syndrome, is a chromosomal condition that is associated with severe mental retardation and certain physical abnormalities. Affected individuals rarely live past infancy because of the life-threatening medical problems associated with this condition. Trisomy 13 affects approximately 1 in 10,000 newborns. The risk of having a child with trisomy 13 increases as a woman gets older.

The chromosome 13 specific region probe is optimized to detect copy numbers of chromosome 13 at 13q14.2 on uncultured amniotic cells. In all PN combinations the 13q14 specific DNA probe is direct-labeled in green with PlatinumBright495.

Cat.# KBI-40001 PN 13 (13q14)




PN 13 (13q14) probe hybridized to a normal metaphase (2G).

Trisomy 18 – Edward Syndrome

Trisomy 18 is caused by an extra chromosome 18 and usually consists of mental retardation, small birth size, and many developmental anomalies, including severe microcephaly, prominent occiput, low-set malformed ears, and a characteristic pinched facial appearance. Trisomy 18 occurs 1 in 6000 live births, but spontaneous abortions are common. More than 95% of affected children have complete trisomy 18. The extra chromosome is almost always maternally derived, and advanced maternal age increases risk.

Cat.# KBI-20018-B SE 18 (D18Z1)




SE 18 (D18Z1) probe hybridized to a normal metaphase showing two blue signals (2B).

The chromosome 18 specific Satellite Enumeration (SE 18) probe (D18Z1) is optimized to detect copy numbers of chromosome 18 at 18p11-18q11 on uncultured amniotic cells. In all PN combinations the 18 SE centromeric DNA probe is offered direct-labeled in blue with PlatinumBright415.


Other combinations

Cat.# KBI-40003



PN 13 (13q14) / 21 (21q22) probe hybridized to a normal interphase (2R2G).

Cat.# KBI-40008



Ploidyscreen (21, X, Y) showing trisomy 21.

Material kindly provided by Prof. Wegner, Berlin.

Sex Chromosome Abnormalities

Chromosomal abnormalities involving the X and Y chromosome (sex chromosomes) are slightly less common than autosomal abnormalities and are usually much less severe in their effects. The high frequency of people with sex chromosome aberrations is partly due to the fact that they are rarely lethal conditions.

- **Turner syndrome** occurs when females inherit only one X chromosome – their genotype is XO.
- **Metafemales or triple-X females**, inherit three X chromosomes – their genotype is XXX or more rarely XXXX or XXXXX.
- **Klinefelter syndrome** males inherit one or more extra X chromosomes – their genotype is XXY or more rarely XXXY, XXXXY, or XY/XXY mosaic.
- **XYY syndrome** males inherit an extra Y chromosome – their genotype is XYY.

The chromosome X specific SE probe (DXZ1) is optimized to detect copy numbers of chromosome X at Xp11-Xq11 on uncultured amniotic cells. The chromosome Y specific SE probe (DYZ3) is optimized to detect copy numbers of chromosome Y at Yp11-Yq11 on uncultured amniotic cells. In all Prenatal Probes combinations the X SE centromeric DNA probe is offered direct-labeled in green with PlatinumBright495.

Technical information

In most Prenatal Probes combinations the Y SE centromeric DNA probe is offered direct-labeled in red with PlatinumBright550, except for the KBI-40008 where the Y SE is labeled in blue with PlatinumBright415. All prenatal probes are in a Ready-to-Use for more convenience. This format still allows adding of a SE 18 probe.

Product and ordering information

| Description | Color | Tests | Cat# |
|-----------------------------|----------------|-------|-------------|
| PN 13 (13q14) | green | 10 | KBI-40001 |
| PN 21 (21q22) | red | 10 | KBI-40002 |
| PN 13 (13q14) / 21 (21q22) | green/red | 10 | KBI-40003 |
| SE 18 (D18Z1) 5x conc | blue | 10 | KBI-20018-B |
| SE X (DXZ1) / SE Y (DYZ3) | green/red | 10 | KBI-20030 |
| SE 7 (D7Z1) / SE 8 (D8Z1) | red/green | 10 | KBI-20031 |
| SE (X,Y,18) | green/red/blue | 10 | KBI-20032 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 10 | KBI-40005 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 30 | KBI-40006 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 50 | KBI-40007 |
| PloidyScreen (21, X, Y) | red/green/blue | 20 | KBI-40008 |

MICRODELETION DNA PROBES

Microdeletion syndromes are usually caused by a chromosomal deletion spanning one or several genes that are too small to be detected under the microscope using conventional cytogenetic methods. Fluorescence *In Situ* Hybridization (FISH) can be employed to identify such deletions and therefore becomes the method of choice for diagnosing microdeletion syndromes.

The REPEAT-FREE™ POSEIDON™ Microdeletion Probes are direct labeled, Ready-to-Use in hybridization buffer and available in a 5 or a 10 test kit.

Table of Contents

| Description | Cat# | Page | Description | Cat# | Page |
|--|-----------|-------|--|-----------|------|
| MD Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) | KBI-40101 | 87 | MD Angelman UBE3A (15q11) / PML(15q24) | KBI-40110 | 86 |
| MD DiGeorge "N25" (22q11) / 22q13 (SHANK3) | KBI-40102 | 81,83 | MD Williams-Beuren ELN (7q11) / 7q22 | KBI-40111 | 86 |
| MD DiGeorge Tuple (22q11) / 22q13 (SHANK3) | KBI-40103 | 81,83 | MD Short Stature (Xp22) / SE X | KBI-40112 | 89 |
| MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) | KBI-40104 | 82,83 | MD NSD1 (5q35)/ hTERT (5p15) | KBI-40113 | 84 |
| MD DiGeorge II (10p14) / SE 10 | KBI-40105 | 83 | MD NF1 (17q11) / MPO (17q22) | KBI-40114 | 84 |
| MD Cri-Du-Chat CTNND (5p15) / 5q31 | KBI-40106 | 88 | MD STS (Xp22) / KAL (Xp22) / SE X TC | KBI-40115 | 90 |
| MD Wolf-Hirschhorn WHSC1 (4p16) / SE 4 | KBI-40107 | 88 | MD IGF1R (15q26) / 15q11 | KBI-40116 | 91 |
| MD X-Inactivation XIST (Xq13) / SE X | KBI-40108 | 89 | MD GATA4 (8p23) / SE 8 | KBI-40118 | V!! |
| MD Prader-Willi SNRPN (15q11) / PML(15q24) | KBI-40109 | 85 | | | |

DiGeorge / Velocardiofacial Syndrome (VCFS)

Microdeletion of chromosome 22 accounts for more than 90% of cases of DiGeorge anomaly which has an incidence of 1 in 4000 live births. Deletions of chromosome 22q11.2 are found in the vast majority of patients with DiGeorge anomaly and VCFS.

Most deletions are *de novo*, with 10% or less inherited from an affected parent. All probes that are currently in use to detect deletions in DiGeorge and VCFS are located within the described minimal critical region of 1.5 Mb.

MD DiGeorge "N25" (22q11) / 22q13 (SHANK3)

The DiGeorge "N25" probe was the first commercial microdeletion probe for chromosome 22q and detects the locus D22S75. This marker is located between DGCR2 and CLH22 (Clathrin). Both genes have been extensively investigated and their role in DiGeorge syndrome is well established.



The DiGeorge "N25" region probe covers the marker "N25" (D22S75) and adjacent region of CLH22 (Clathrin gene region) and DGCR2 (DiGeorge critical region gene 2). The SHANK3 probe at 22q13 is serving as internal control.

MD DiGeorge "TUPLE" (22q11) / 22q13 (SHANK3)

The DiGeorge "Tuple" probe targets a putative transcriptional regulator (TUPLE1 or HIRA, HIR histone cell cycle regulation defective homolog A) which also has been identified to lie within the commonly deleted region DiGeorge syndrome. This probe is located distally to the "N25" probe. The DiGeorge "Tuple" region probe is optimized to detect copy numbers of the Tuple (Hira) gene region at 22q11.2.



The SHANK3 probe at 22q13 is serving as internal control.

Cat.# KBI-40102 MD DiGeorge "N25" (22q11) / 22q13 (SHANK3)



MD DiGeorge "N25" (22q11) / 22q13 (SHANK3) probe hybridized to a normal metaphase (2R2G).

Cat.# KBI-40103 DiGeorge "Tuple" (22q11) / 22q13 (SHANK3)



MD DiGeorge "Tuple" (22q11) / 22q13 (SHANK3) probe hybridized to a normal metaphase (2R2G).

Literature:

Sirotkin et al, 1996, Hum Mol Genet, 5: 617-624.
Holmes et al, 1997, Hum Mol Genet, 6: 357-367.
Wilson, et al, 2003, J Med Genet 40; 575-584.
Luciani, et al, 2003, J Med Genet 40; 690-696.

Literature:

Lorain et al, 1996, Genome Res, 6: 43-50.

| Ordering information | Gene Region | Tests | Cat# |
|--|-------------|-------|-----------|
| MD DiGeorge "N25" (22q11) / 22q13 (SHANK3) | N25 | 10 | KBI-40102 |

| Ordering information | Gene Region | Tests | Cat# |
|--|-------------|-------|-----------|
| MD DiGeorge Tuple (22q11) / 22q13 (SHANK3) | TUPLE | 10 | KBI-40103 |

MICRODELETION DNA PROBES

MD DiGeorge T-Box1 (22q11) / 22q13 (SHANK3)



The 22q11 deletion in DiGeorge syndrome/VCFS is characterized by defects in the derivatives of the pharyngeal apparatus. TBX1, a member of the T-box transcription factor family, is required for normal development of the pharyngeal arch arteries. Haploinsufficiency of TBX1 has been demonstrated to be sufficient to generate at least one important component of the DiGeorge syndrome phenotype in mice. The TBX1 is also located within the minimal critical DiGeorge region in humans.

The DiGeorge TBX1 region probe is optimized to detect copy numbers of the TBX1 gene region at 22q11.2. The subtelomeric (ST) 22qter DNA probe is included as control probe. The SHANK3 probe at 22q13 is serving as internal control.

Phelan-McDermid Syndrome

The 22q13 deletion syndrome (or Phelan-McDermid syndrome) is characterized by moderate to profound mental retardation, delay/absence of expressive speech, hypotonia, normal to accelerated growth, and mild dysmorphic features. A terminal deletion including the SHANK3 gene region has been identified for this syndrome.

Cat.# KBI-40104 MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3)



MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) probe hybridized to DiGeorge patient material showing a deletion of the TBX1 gene region at 22q11 (1R2G).

Image kindly provided by Dr. F. Girard-Lemaire, Service de Cytogénétique (Dr. Flori), CHU Strasbourg.


Literature:

Lindsay et al, 2001, Nature, 410: 97-101.
Merscher et al, 2001, Cell, 104: 619-629.
Paylor et al, 2006, PNAS, 103: 7729-7734.

| Ordering information | Gene Region | Tests | Cat# |
|---|-------------|-------|-----------|
| MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) | TBX1 | 10 | KBI-40104 |

MD DiGeorge (22q11) / 22q13 (SHANK3)

The 22q13 DNA probe is optimized to detect copy numbers of the SHANK3 gene region at 22q13. The DiGeorge region probe at 22q11 is serving as internal control.



MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) probe hybridized to patient material showing a deletion of the SHANK3 region at 22q13 (2R1G).

Literature:

Wilson, et al, 2003, J Med Genet 40; 575-584.
Luciani, et al, 2003, J Med Genet 40; 690-696.


| Ordering information | Gene Region | Tests | Cat# |
|---|-------------|-------|-----------|
| MD DiGeorge "N25" (22q11) / 22q13 (SHANK3) | N25 | 10 | KBI-40102 |
| MD DiGeorge Tuple (22q11) / 22q13 (SHANK3) | TUPLE | 10 | KBI-40103 |
| MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) | TBX1 | 10 | KBI-40104 |

MD DiGeorge II (10p14) / SE 10

DiGeorge and VCFS present many clinical problems and are frequently associated with deletions within 22q11.2 (see previous probes), but a number of cases have no detectable molecular defect of this region. A number of single case reports with deletions of 10p suggest genetic heterogeneity of DiGeorge syndrome. FISH analysis demonstrates that these patients have overlapping deletions at the 10p13/10p14 boundary. The shortest region of deletion overlap (SRO) has been identified in a 1 cM interval including makers D10S547 and D10S585.

The DiGeorge II region probe is optimized to detect copy numbers of the DGSII at 10p14. The chromosome 10 satellite enumeration (SE 10) probe at D10Z1 is included to facilitate chromosome identification.

Cat.# KBI-40105 MD DiGeorge II (10p14) / SE 10



MD DiGeorge II(10p14) / SE 10 probe hybridized to DiGeorge II patient material showing a deletion of the DGSII region at 10p14 (1R2G).

Image kindly provided by Azzedine Aboura, Hôpital Robert Debré Paris.

Literature:

Monaco et al, 1991, Am J Med Genet, 39: 215-216.
Schuffenhauer et al, 1998, Eur J Hum Genet, 6: 213-225.

| Ordering information | Gene Region | Tests | Cat# |
|--------------------------------|-------------|-------|-----------|
| MD DiGeorge II (10p14) / SE 10 | 10p- | 10 | KBI-40105 |

MICRODELETION DNA PROBES



MD NSD1 (5q35) / hTERT (5p15)

NSD1 microdeletions (chromosome 5q35) are the major cause of Sotos syndrome, and occur in some cases of Weaver syndrome. Sotos is a childhood overgrowth characterized by distinctive craniofacial features, advanced bone age, and mental retardation. Weaver syndrome is characterized by the same criteria but has its own specific facial characteristics.


Sotos syndrome is inherited in an autosomal dominant manner. While 50% of Sotos patients in Asia are showing a chromosomal microdeletion, only 9% deletion cases are observed in the affected European population.

The NSD1 (5q35) region probe is optimized to detect copy numbers of the NSD1 gene region at 5q35.2. The hTERT region specific DNA probe at 5p15 is included as control probe.

Cat.# KBI-40113 NSD1 (5q35) / hTERT (5p15)



NSD1 (5q35) / hTERT (5p15) probe hybridized to a normal metaphase (2R2G).



NSD1 (5q35) / hTERT (5p15) probe hybridized to patient material showing a microdeletion of the NSD1 gene region at 5q35 (1R2G).

Literature:

Douglas et al, 2003, Am. J. Hum. Genet. 72; 132-143.
Rio et al, 2003, J. Med. Genet. 40; 436-440.

| Ordering information | Gene Region | Tests | Cat# |
|-------------------------------|-------------|-------|-----------|
| MD NSD1 (5q35) / hTERT (5p15) | NSD1 | 10 | KBI-40113 |



MD NF1 (17q11) / MPO (17q22)

NF1, or von Recklinghausen disease, is one of the most common hereditary neurocutaneous disorders in humans and one of the most common single gene syndromes. Clinically, NF1 is characterized by café-au-lait spots, freckling, skin neurofibroma, plexiform neurofibroma, bone defects, Lisch nodules and tumors of the central nervous system. The responsible gene, NF1 (neurofibromin), was identified on chromosome 17q11. Whole NF1 gene deletions occur in 4%-5% of individuals with NF1 and can be detected by FISH analysis.

The NF1 (17q11) region probe is optimized to detect copy numbers of the NF1 gene region at 17q11.2.

The MPO region specific DNA probe at 17q22 is included as control probe.

Cat.# KBI-40114 NF1 (17q11) / MPO (17q22)



NF1 (17q11) / MPO (17q22) probe hybridized to patient material showing a deletion of NF1 gene region at 17q11 (1R2G).

Literature:

Riva P et al, 2000, Am.J.Hum.Genet. 66; 100-109.
Dorschner et al, 2000, Hum.Mol.Genet. 9; 35-46.

| Ordering information | Gene Region | Tests | Cat# |
|------------------------------|-------------|-------|-----------|
| MD NF1 (17q11) / MPO (17q22) | NF1 | 10 | KBI-40114 |

Prader-Willi/Angelman Syndrome



Prader-Willi syndrome (PWS) and Angelman syndrome (AS) are clinically distinct complex disorders mapped to chromosome 15q11-q13. They both have characteristic neurologic, developmental, and behavioral phenotypes plus other structural and functional abnormalities. However, the cognitive and neurologic impairment is more severe in AS, including seizures and ataxia. The behavioral and endocrine disorders are more severe in PWS, including obsessive-compulsive symptoms and hypothalamic insufficiency. Both disorders can result from microdeletion, uniparental disomy, or an imprinting center defect in 15q11-q13.

MD Prader-Willi SNRPN (15q11) / PML (15q24)

Prader-Willi syndrome (PWS) is a clinically distinct disorder including diminished fetal activity, obesity, hypotonia, mental retardation, short stature, hypogonadotropic hypogonadism, strabismus, and small hands and feet.

Approximately 70% of cases of PWS arise from paternal deletion of the 15q11-q13 region including the gene SNRPN (small nuclear ribonucleoprotein polypeptide N). The PWS SNRPN region probe is optimized to detect copy numbers of the SNRPN gene region at 15q11. The PML (promyelocytic leukemia) gene specific DNA probe at 15q24 is included as control probe.

Cat.# KBI-40109 Prader-Willi SNRPN (15q11) / PML (15q24)



Prader-Willi SNRPN (15q11) / PML (15q24) probe hybridized to a normal interphase/metaphase (2R2G).

Literature:
Knoll et al, 1989, Am J Med Genet, 32: 285-290.
Ozcelik et al, 1992, Nat Genet, 2: 265-269.

| Ordering information | Gene Region | Tests | Cat# |
|--|-------------|-------|-----------|
| MD Prader-Willi SNRP (15q11) PML (15q24) | SNRPN | 10 | KBI-40109 |



MICRODELETION DNA PROBES

MD Angelman UBE3A (15q11) / PML (15q24)

Angelman syndrome (AS) is characterized by severe developmental delay or mental retardation, severe speech impairment, gait ataxia and/or tremulousness of the limbs, and an unique behavior with an inappropriate happy demeanor that includes frequent laughing, smiling, and excitability. In addition, microcephaly and seizures are common. AS is caused by absence of a maternal contribution to the imprinted region on chromosome 15q11-q13 including the UBE3A gene.

The AS UBE3A region probe is optimized to detect copy numbers of the UBE3A gene region at 15q11. The PML (promyelocytic leukemia) gene specific DNA probe at 15q24 is included as control probe.

Cat.# KBI-40110 Angelman UBE3A (15q11) / PML (15q24)



Angelman UBE3A (15q11) / PML (15q24) probe hybridized to a normal interphase/metaphase (2R2G).

Literature:

Matsuura et al, 1997, Nat Genet, 15: 74-77.
Burger et al, 2002, Am J Med Genet, 111: 233-237.



| Ordering information | Gene Region | Tests | Cat# |
|---|-------------|-------|-----------|
| MD Angelman UBE3A (15q11) / PML (15q24) | UBE3A | 10 | KBI-40110 |

MD Williams-Beuren ELN (7q11) / 7q22

Williams-Beuren syndrome (WS) is characterized by cardiovascular disease, distinctive facial features, connective tissue abnormalities, mental retardation and endocrine abnormalities. Over 99% of individuals with the clinical diagnosis of WS have this contiguous gene deletion, that encompasses the elastin (ELN) gene region including ELN, LIMK1, and the D7S613 locus.

The Williams-Beuren region probe is optimized to detect copy numbers of the ELN gene region at 7q11. The 7q22 region specific DNA probe at 7q22 is included as control probe.

Cat.# KBI-40111 Williams-Beuren ELN (7q11) / 7q22



Williams-Beuren ELN (7q11) / 7q22 probe hybridized to a normal metaphase (2RG).

Literature:

Ewart, et al, 1993, Nat Genet, 5: 11-16.
Botta et al, 1999, J Med Genet, 36: 478-480.



| Ordering information | Gene Region | Tests | Cat# |
|--------------------------------------|-------------|-------|-----------|
| MD Williams-Beuren ELN (7q11) / 7q22 | ELN | 10 | KBI-40111 |

MD Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11)

Miller-Dieker Syndrome (MDS) is characterized by classical lissencephaly and distinct facial features. The lissencephaly represents the severe end of the spectrum with generalized agyria or agyria with some frontal pachygyria. Submicroscopic deletions of 17p13.3 including the LIS1 (now called PAFAH1B1, platelet-activating factor acetylhydrolase) gene are found in almost 100% of patients.

The Miller-Dieker region probe is optimized to detect copy numbers of the PAFAH1B1 gene (LIS1) region at 17p13.3. The Smith-Magenis RAI1 region probe at 17p11.2 is serving as internal control.

Cat.# KBI-40101 Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11)




Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) probe hybridized to a normal metaphase (2R1G).

Literature:
Kuwnano et al, 1991, Am J Hum Genet, 49: 707-714.
Cardoso et al, 2003, Am J Hum Genet, 72: 918-930.

Smith-Magenis Syndrome (SMS) is characterized by distinctive facial features that progress with age, developmental delay, cognitive impairment, and behavioral abnormalities. Molecular cytogenetic analysis by FISH using a DNA probe specific for the SMS critical region is recommended in cases of submicroscopic deletions and/or to resolve equivocal cases. RAI1 is the only gene known to account for a majority of features in SMS. All 17p11.2 deletions associated with SMS include a deletion of RAI1.

The Smith-Magenis region probe is optimized to detect copy numbers of the RAI1 gene region involved in Smith-Magenis syndrome at 17p11.2. The Miller-Dieker LIS1 probe at 17p13.3 is serving as internal control.



Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) probe hybridized to Smith-Magenis patient material showing a deletion of the RAI1 gene region at 17p11 (2R1G).
Image kindly provided by Prof. Jauch, Heidelberg.

Literature:
Smith et al, 1986, Am J Med Genet, 24: 393-414.
Greenberg et al, 1991, Am J Med Genet, 49: 1207-1218.
Vlangos et al, 2005, Am J Med Genet, 132: 278-282.

| Ordering information | Gene Region | Tests | Cat# |
|--|-------------|-------|-----------|
| MD Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) | LIS1/RAI1 | 10 | KBI-40101 |



MICRODELETION DNA PROBES

MD Wolf-Hirschhorn WHSC1 (4p16)

Wolf-Hirschhorn syndrome (WHS) affected individuals have prenatal-onset growth deficiency followed by postnatal growth retardation and hypotonia with muscle under-development. Developmental delay/mental retardation of variable degree is present in all. FISH analysis using a WHSC1 specific probe for chromosomal locus 4p16.3 detects more than 95% of deletions in WHS.

The Wolf-Hirschhorn region probe is optimized to detect copy numbers of the Wolf-Hirschhorn critical region at 4p16. The chromosome 4 satellite enumeration (SE 4) probe at D4Z1 is included to facilitate chromosome identification.

Cat.# KBI-40107 Wolf-Hirschhorn WHSC1 (4p16) / SE 4



Wolf-Hirschhorn WHSC1 (4p16) / SE 4 probe hybridized to Wolf-Hirschhorn patient material showing a deletion of the WHSC1 gene region at 4p16 (1R2G). Image kindly provided by Prof. Zollino, Rome.

Literature:

Gandelman et al, 1992, Am J Hum Genet, 51: 571-578.
Wright et al, 1997, Hum Mol Genet, 6: 317-324.



| Ordering information | Gene Region | Tests | Cat# |
|--|-------------|-------|-----------|
| MD Wolf-Hirschhorn WHSC1 (4p16) / SE 4 | WHSC1 | 10 | KBI-40107 |

MD Cri-Du-Chat CTNND (5p15) / 5q31

Cri-Du-Chat syndrome is an autosomal deletion syndrome caused by a partial deletion of chromosome 5p. It is characterized by a distinctive, high-pitched, catlike cry in infancy with growth failure, microcephaly, facial abnormalities, and mental retardation throughout life. Loss of a small region in band 5p15.2 (Cri-Du-Chat critical region) correlates with all the clinical features of the syndrome with the exception of the catlike cry, which maps to band 5p15.3 (catlike cry critical region).

The Cri-Du-Chat region probe is optimized to detect copy numbers at the CTNND2 gene region in the Cri-Du-Chat critical region at 5p15.2. The 5q31 specific DNA probe is included as control probe.

Cat.# KBI-40106 Cri-Du-Chat CTNND (5p15) / 5q31



Cri-Du-Chat CTNND (5p15) / 5q31 probe hybridized to a normal metaphase (2R2G).

Literature:

Overhauser et al, 1994, Hum Mol Genet, 3: 247-252.
Gersh et al, 1997, Cytogenet Cell Genet, 77: 246-251.



| Ordering information | Gene Region | Tests | Cat# |
|------------------------------------|-------------|-------|-----------|
| MD Cri-Du-Chat CTNND (5p15) / 5q31 | 5p- | 10 | KBI-40106 |

MD X-Inactivation XIST (Xq13) / SE X

The XIST locus is expressed only from the inactive X chromosome, resides at the putative X inactivation centre, and is considered a prime player in the initiation of mammalian X dosage compensation. The severe phenotype of human females whose karyotype includes tiny ring X chromosomes has been attributed to the inability of the small ring X chromosome to inactivate. Many of the ring chromosomes lack the XIST locus, consistent with XIST being necessary for cis inactivation.

The XIST specific DNA probe is optimized to detect copy numbers of the XIST region at Xq13. The chromosome X satellite enumeration (SE X) probe at DXZ1 is added to facilitate chromosome identification.

Cat.# KBI-40108 X-Inactivation XIST (Xq13) / SE X



X-Inactivation XIST (Xq13) / SE X probe hybridized to a male metaphase (1R1G).

Literature:

Migeon et al, 1993, PNAS, 90: 12025-12029.
Jani et al, 1995, Genomics, 27: 182-188.



| Ordering information | Gene Region | Tests | Cat# |
|--------------------------------------|-------------|-------|-----------|
| MD X-Inactivation XIST (Xq13) / SE X | XIST | 10 | KBI-40108 |

MD Short Stature (Xp22) / SE X

Individuals with SHOX-related short stature have disproportionate short stature and/or wrist abnormalities consistent with those described in Madelung deformity. The SHOX genes located on the pseudoautosomal regions of the X and Y chromosomes are the only genes known to be associated with SHOX-related haploinsufficiency.

The SHOX region probe is optimized to detect copy numbers of the SHOX gene region at Xp22. The chromosome X satellite enumeration (SE X) probe at DXZ1 is added to facilitate chromosome identification.

Cat.# KBI-40112 Short stature (Xp22) / SE X



Short stature (Xp22) / SE X probe hybridized to a male metaphase (2R1G).

Literature:

Rao et al, 1997, Hum Genet, 100: 236-239.
Morizio et al, 2003, Am J Med Genet, 119: 293-296.

| Ordering information | Gene Region | Tests | Cat# |
|--------------------------------|-------------|-------|-----------|
| MD Short Stature (Xp22) / SE X | SHOX | 10 | KBI-40112 |



MICRODELETION DNA PROBES

MD STS (Xp22) / KAL (Xp22) / SE X TC

STS (Steroid Sulfatase) disease is a chromosome X-linked disorder associated with a microdeletion of the gene within the Xp22.3 region. Deletion of the steroid sulfatase gene has been detected in individuals with recessive X-linked ichthyosis, the disease been considered one of the most frequent human enzyme deficient disorders. KAL1 (Kallmann syndrome interval gene-1) maps to the Kallmann syndrome critical region on the distal short arm of the human X chromosome. Individuals with Kallmann syndrome suffers from hypogonadotropic hypogonadism and anosmia, with clinical features of variable phenotype. It affects approximately 1 in 8000 males and 1 in 40000 females.

The STS (Xp22) region probe is optimized to detect copy numbers of the STS gene region at Xp22. The KAL (Xp 22) region probe is optimized to detect copy numbers of the KAL gene region at Xp22. The Chromosome X Satellite Enumeration (SE X) probe at DXZ1 is included to facilitate chromosome identification.

Cat.# KBI-40115 STS (Xp22) / KAL (Xp22) / SE X TC



STS (Xp22) / KAL (Xp22) / SE X TC probe hybridized to male patient material showing a deletion of the STS gene region (1R1B).
Material kindly provided by Necker hospital, Paris.

Literature:

Alper in et al, 1997, J. Biol. Chem 272; 20756-20763.
Meroni et al, 1996, Hum. Mol. Genet. 5; 423-431.



| Ordering information | Gene Region | Tests | Cat# |
|--------------------------------------|-------------|-------|-----------|
| MD STS (Xp22) / KAL (Xp22) / SE X TC | STS/KAL | 10 | KBI-40115 |

MD IGF1R (15q26) / 15q11

Congenital diaphragmatic hernia (CDH) is a severe, life-threatening, congenital anomaly characterized by variable defects in the diaphragm, pulmonary hypoplasia, and postnatal pulmonary hypertension. Deletion of the IGF1R (insulin-like growth factor 1 receptor) gene region at 15q25 is the most frequent anomaly found in CDH. The type 1 IGF receptor at 15q26 is required for normal embryonic and postnatal growth. Deletions, but also gain of an approximately 5 Mb region including the IGF1R gene, have been found to have a profound effect on prenatal and early postnatal growth.

The IGF1R (15q26) specific probe is optimized to detect copy numbers of the IGF1R gene region at region 15q26. The 15q11 (SNRPN / UBE3A) specific region probe is included to facilitate chromosome identification.

Cat.# KBI-40116 IGF1R (15q26) / 15q11



IGF1R (15q26) / 15q11 probe hybridized to patient material showing a deletion of the IGF1R gene region at 15q26 (1R2G).

Literature:

Faivre et al, 2002, Eur J Hum Genet. 10 ; 699-706.


Okubo et al, 2003, J Clin Endocrinol. Metab 88 ; 5981-5988.

| Ordering information | Gene Region | Tests | Cat# |
|--------------------------|-------------|-------|-----------|
| MD IGF1R (15q26) / 15q11 | IGF1R | 10 | KBI-40116 |


SUB-TELOMERE DNA PROBES

Telomeres are specialized DNA-protein structures containing long stretches of (TTAGGG) n repeats at the end of all chromosomes. They protect chromosomes from degradation and end-to-end fusion with other chromosomes. The region adjacent to sequences containing telomeric repeats is called the sub-telomer which has been found to be relative gene-rich. Cytogenetic analysis of sub-telomeric regions is difficult due to low resolution using conventional banding techniques. Cryptic deletions and rearrangements have been associated with unexplained mental retardation and congenital abnormalities. FISH probes specific for the sub-telomeric regions of the terminal chromosome regions are essential to detect such subtle rearrangements.

REPEAT-FREE™ POSEIDON™ Sub-Telomeric (ST) DNA Probes are supplied in a 5 test format together with hybridization buffer, direct labeled in either red (PlatinumBright™550), green (PlatinumBright495) or blue (PlatinumBright415). The ST probes are provided in a 5x concentrated format to allow, mixing of up to 5 ST probes in a single hybridization assay.




Dual-Color



Sub-Telomere probes 6pter and 2pter showing trisomy of tel 6p and monosomy of tel 2p.
Image kindly provided by Institute Pasteur, Paris.

Triple-Color



Combination of Sub-Telomere probes 6qter (red), 8qter (green) and 9pter (blue).

Product and ordering information

| <i>Location</i> | <i>Marker</i> | <i>Probe Size (kb)</i> | <i>Distance from telomere (kb)</i> | <i>Color</i> | <i>Tests</i> | <i>Cat.# *</i> |
|---------------------|---------------|------------------------|------------------------------------|--------------------|--------------|----------------|
| Sub-Telomere 1pter | D1S2217 | 170 | 800 | red, green or blue | 5 | KBI-40201 |
| Sub-Telomere 1qter | D1S555 | 170 | 350 | red, green or blue | 5 | KBI-40202 |
| Sub-Telomere 2pter | D2S2147 | 210 | 300 | red, green or blue | 5 | KBI-40203 |
| Sub-Telomere 2qter | D2S2142 | 185 | 800 | red, green or blue | 5 | KBI-40204 |
| Sub-Telomere 3pter | D3S4558 | 175 | 450 | red, green or blue | 5 | KBI-40205 |
| Sub-Telomere 3qter | D3S4168 | 170 | 900 | red, green or blue | 5 | KBI-40206 |
| Sub-Telomere 4pter | D4S3360 | 180 | 100 | red, green or blue | 5 | KBI-40207 |
| Sub-Telomere 4qter | D4S2283 | 190 | 700 | red, green or blue | 5 | KBI-40208 |
| Sub-Telomere 5pter | D5S2488 | 175 | 180 | red, green or blue | 5 | KBI-40209 |
| Sub-Telomere 5qter | D5S2006 | 270 | 600 | red, green or blue | 5 | KBI-40210 |
| Sub-Telomere 6pter | RH40931 | 110 | 350 | red, green or blue | 5 | KBI-40211 |
| Sub-Telomere 6qter | D6S2523 | 165 | 250 | red, green or blue | 5 | KBI-40212 |
| Sub-Telomere 7pter | D7S2644 | 220 | 850 | red, green or blue | 5 | KBI-40213 |
| Sub-Telomere 7qter | D7S427 | 180 | 200 | red, green or blue | 5 | KBI-40214 |
| Sub-Telomere 8pter | RH65733 | 180 | 550 | red, green or blue | 5 | KBI-40215 |
| Sub-Telomere 8qter | D8S595 | 210 | 200 | red, green or blue | 5 | KBI-40216 |
| Sub-Telomere 9pter | D9S917 | 190 | 450 | red, green or blue | 5 | KBI-40217 |
| Sub-Telomere 9qter | D9S1838 | 185 | 500 | red, green or blue | 5 | KBI-40218 |
| Sub-Telomere 10pter | D10S2488 | 180 | 350 | red, green or blue | 5 | KBI-40219 |
| Sub-Telomere 10qter | D10S2290 | 230 | 350 | red, green or blue | 5 | KBI-40220 |
| Sub-Telomere 11pter | D11S1363 | 155 | 1050 | red, green or blue | 5 | KBI-40221 |
| Sub-Telomere 11qter | D11S4437 | 170 | 300 | red, green or blue | 5 | KBI-40222 |
| Sub-Telomere 12pter | D12S158 | 185 | 150 | red, green or blue | 5 | KBI-40223 |
| Sub-Telomere 12qter | D12S399 | 190 | 180 | red, green or blue | 5 | KBI-40224 |
| Sub-Telomere 13qter | D13S1160 | 190 | 90 | red, green or blue | 5 | KBI-40225 |
| Sub-Telomere 14qter | D14S1419 | 170 | 250 | red, green or blue | 5 | KBI-40226 |
| Sub-Telomere 15qter | RH54179 | 150 | 250 | red, green or blue | 5 | KBI-40227 |
| Sub-Telomere 16pter | D16S521 | 150 | 40 | red, green or blue | 5 | KBI-40228 |
| Sub-Telomere 16qter | RH25942 | 180 | 240 | red, green or blue | 5 | KBI-40229 |
| Sub-Telomere 17pter | D17S643 | 110 | 80 | red, green or blue | 5 | KBI-40230 |
| Sub-Telomere 17qter | D17S724 | 70 | 500 | red, green or blue | 5 | KBI-40231 |
| Sub-Telomere 18pter | D18S1244 | 180 | 200 | red, green or blue | 5 | KBI-40232 |
| Sub-Telomere 18qter | D18S1390 | 220 | 160 | red, green or blue | 5 | KBI-40233 |
| Sub-Telomere 19pter | D19S814 | 220 | 550 | red, green or blue | 5 | KBI-40234 |
| Sub-Telomere 19qter | D19S989 | 160 | 300 | red, green or blue | 5 | KBI-40235 |
| Sub-Telomere 20pter | D20S1156 | 240 | 180 | red, green or blue | 5 | KBI-40236 |
| Sub-Telomere 20qter | RH44234 | 170 | 350 | red, green or blue | 5 | KBI-40237 |
| Sub-Telomere 21qter | D21S1446 | 190 | 80 | red, green or blue | 5 | KBI-40238 |
| Sub-Telomere 22qter | D22S1056 | 200 | 850 | red, green or blue | 5 | KBI-40239 |
| Sub-Telomere XYpter | DXYS130 | 180 | 400 | red, green or blue | 5 | KBI-40240 |
| Sub-Telomere XYqter | DXYS224 | 160 | 70 | red, green or blue | 5 | KBI-40241 |

* Add -G for Green, -R for Red, -B for Blue (available on request)

SATELLITE ENUMERATION DNA PROBES

The primary constriction, called the centromere, is a common feature of chromosomes necessary for cell division. Presence of repetitive sequences in the centromeric regions have been proven to be essential. In humans, and many other species, specific repetitive sequences, called 'Satellites' are characteristic for the centromere in general. Most chromosomes also have repetitive sequences which are specific for individual chromosomes and can be used for precise identification and enumeration of human chromosomes in metaphase and interphase cells. The REPEAT-FREE™ POSEIDON™ Satellite Enumeration DNA probes allow rapid and specific chromosome analysis, marker chromosome identification and the detection of aneuploidy. Essentially, these probes can be used in all aspects of routine work in genetics and oncology/pathology. Due to the sharp and bright signals produced, the Satellite probes can easily be used on various sample types, such as cultured cells, touch preparations buccal smears, cytopins, frozen and paraffin-embedded tissue sections, sputum samples, sperm samples, and bladder washes.

POSEIDON Satellite Enumeration (SE) DNA Probes are supplied in a 10 test format together with hybridization buffer, direct labeled in either red (PlatinumBright™550), green (PlatinumBright495) or blue (PlatinumBright415). The SE probes are provided in a 5x concentrated format to allow, mixing of up to 5 SE probes in a single hybridization assay. The SE combination kits are supplied in a Ready-to-Use dual color format.


All Human Centromere (AHC) probe

This probe specifically hybridizes to the centromeric region of all human chromosomes, and is labeled with high fluorescence intensity in red or green color. POSEIDON™ AHC probe is ideal for studying numerical chromosome aberrations, studying aneuploidy, polyploidy, dicentrics, tracentrics, and other complex aberrations. It can also be used for general numerical chromosome analysis. POSEIDON AHC probe is for RUO and is not meant to be used for medical purposes or as a diagnostic tool.


Cat.# KBI-20000 All Human Centromere (AHC) probe



Cat.# KBI-20000



Cat# KBI-20018-B



Cat# KBI-20031

SE7/SE8 showing trisomy 8.


Material kindly provided by Dr. Balogh, Budapest.

Note: Satellite sequences share some degree of homology between the sequences from chromosome to chromosome. Therefore the recommended stringency conditions in hybridization and posthybridization washes must be followed to provide optimal results. Chromosomes 1/5/19, 13/21, and 14/22 share the same repetitive sequences and cannot be differentiated by chromosome specific repeats.

Acro-P-Arms NOR

The NOR (Nucleolar Organizer Region) is located on every p-arm of the human acrocentric chromosomes. Enlargement of the acrocentric p-arms can be caused by an unusual variant or a translocation event. NOR stain of the p-arms is useful to detect such a p-arm variant. In the classification of small supernumerary marker chromosomes (SMCs) the Acro-P-Arms NOR probe can detect the origin of DNA, in which about 80% will turn out to be derived from the acrocentric chromosomes.

The Acro-P-Arms NOR probe is optimized to detect the short (p) arm of all acrocentric human chromosomes. The probe is intended to be used on metaphase/interphase spreads.



Cat# KBI-20033-R

Image kindly provided by Dr. Reboul, Nîmes.

Literature:

Starke H et al., 2003, Hum Genet, 114; 51-67.

Starke H et al., 2005, J Histochem Cytochem, 53, 359-360.

Product and ordering information


| Description | Chromosome | Color | DNA Class | Tests | Cat.# * |
|---------------------------------|----------------------|--------------------|---------------------|-------|-----------|
| SE 1 (1qh) | 1, 1qh | red, green or blue | Satellite III | 10 | KBI-20001 |
| SE 2 (D2Z) | 2 | red, green or blue | α -satellite | 10 | KBI-20002 |
| SE 3 (D3Z1) | 3 | red, green or blue | α -satellite | 10 | KBI-20003 |
| SE 4 (D4Z1) | 4 | red, green or blue | α -satellite | 10 | KBI-20004 |
| SE 6 (D6Z1) | 6 | red, green or blue | α -satellite | 10 | KBI-20006 |
| SE 7 (D7Z1) | 7 | red, green or blue | α -satellite | 10 | KBI-20007 |
| SE 8 (D8Z1) | 8 | red, green or blue | α -satellite | 10 | KBI-20008 |
| SE 9 (classical) | 9 | red, green or blue | α -satellite | 10 | KBI-20009 |
| SE 10 (D10Z1) | 10 | red, green or blue | α -satellite | 10 | KBI-20010 |
| SE 11 (D11Z1) | 11 | red, green or blue | α -satellite | 10 | KBI-20011 |
| SE 12 (D12Z3) | 12 | red, green or blue | α -satellite | 10 | KBI-20012 |
| SE 15 (D15Z) | 15 | red, green or blue | α -satellite | 10 | KBI-20015 |
| SE 16 (D16Z2) | 16 | red, green or blue | α -satellite | 10 | KBI-20016 |
| SE 17 (D17Z1) | 17 | red, green or blue | α -satellite | 10 | KBI-20017 |
| SE 18 (D18Z1) | 18 | red, green or blue | α -satellite | 10 | KBI-20018 |
| SE 20 (D20Z1) | 20 | red, green or blue | α -satellite | 10 | KBI-20020 |
| SE X (DXZ1) | X | red, green or blue | α -satellite | 10 | KBI-20023 |
| SE Y (DYZ3) | Y, centromeric | red, green or blue | α -satellite | 10 | KBI-20024 |
| SE Y classical (DYZ1) | Y, Yqh | red, green or blue | Satellite III | 10 | KBI-20025 |
| SE 1/5/19 (D1Z7) (D5Z2) (D19Z3) | 1, 5, 19 | red, green or blue | α -satellite | 10 | KBI-20026 |
| SE 13/21 (D13Z1) (D21Z1) | 13 and 21 | red, green or blue | α -satellite | 10 | KBI-20027 |
| SE 14/22 (D14Z1) (D22Z1) | 14 and 22 | red, green or blue | α -satellite | 10 | KBI-20028 |
| SE combinations | | | | | |
| SE X (DXZ1) / SE Y (DYZ3) RtU | X and Y | green/red | α -satellite | 10 | KBI-20030 |
| SE 7 (D7Z1) / SE 8 (D8Z1) RtU | 7 and 8 | red/green | α -satellite | 10 | KBI-20031 |
| AHC probe | All Human Centromere | red or green | | 10 | KBI-20000 |
| Acro-P-Arms NOR | 13,14,15,21,22 | red, green or blue | | 10 | KBI-20033 |

* Add -G for Green, -R for Red, -B for Blue (available on request)


WHOLE CHROMOSOME DNA PROBES

The Whole Chromosome probes are used for identifying whole human chromosomes, analysis of translocation events, chromosome rearrangement studies and determining the origin of marker chromosomes. Mutagenesis analysis, radiation and sensitivity testing, and identification of human chromosomes on hybrid cells are investigated using these probes. The REPEAT-FREE™ POSEIDON™ Whole Chromosome DNA probes hybridize to unique sequences spanning the entire length of the target chromosome. These probes are derived from flow-sorted or microdissected chromosomes and provide accurate coverage with excellent signal specificity and high fluorescent intensity. Some minor cross-hybridization may occur at the short arm of acrocentric chromosomes (13, 14, 15, 21 and 22) and in the pseudo-autosomal regions of chromosome X and Y.


REPEAT-FREE POSEIDON Whole Chromosome (WC) DNA Probes are supplied in a 5 test format together with hybridization buffer, direct labeled in either red (PlatinumBright™550), green (PlatinumBright495) or blue (PlatinumBright415). The WC probes are provided in a 5x concentrated format to allow, mixing of up to 5 WC probes in a single hybridization assay.




KBI-30013-R



WC Dual-Color



KBI-30023-G



WC Triple-Color

Note: Due to decondensation of chromosomal DNA signals on interphase, cells can be very diffuse. We therefore recommend not to use WC probes for interphase FISH.

Product and ordering information

| <i>Description</i> | <i>Chromosome</i> | <i>Color</i> | <i>Tests</i> | <i>Cat.# *</i> |
|---------------------|---------------------|--------------------|--------------|----------------|
| Whole Chromosome 1 | Chromosome 1 Paint | Red, green or blue | 5 | KBI-30001 |
| Whole Chromosome 2 | Chromosome 2 Paint | Red, green or blue | 5 | KBI-30002 |
| Whole Chromosome 3 | Chromosome 3 Paint | Red, green or blue | 5 | KBI-30003 |
| Whole Chromosome 4 | Chromosome 4 Paint | Red, green or blue | 5 | KBI-30004 |
| Whole Chromosome 5 | Chromosome 5 Paint | Red, green or blue | 5 | KBI-30005 |
| Whole Chromosome 6 | Chromosome 6 Paint | Red, green or blue | 5 | KBI-30006 |
| Whole Chromosome 7 | Chromosome 7 Paint | Red, green or blue | 5 | KBI-30007 |
| Whole Chromosome 8 | Chromosome 8 Paint | Red, green or blue | 5 | KBI-30008 |
| Whole Chromosome 9 | Chromosome 9 Paint | Red, green or blue | 5 | KBI-30009 |
| Whole Chromosome 10 | Chromosome 10 Paint | Red, green or blue | 5 | KBI-30010 |
| Whole Chromosome 11 | Chromosome 11 Paint | Red, green or blue | 5 | KBI-30011 |
| Whole Chromosome 12 | Chromosome 12 Paint | Red, green or blue | 5 | KBI-30012 |
| Whole Chromosome 13 | Chromosome 13 Paint | Red, green or blue | 5 | KBI-30013 |
| Whole Chromosome 14 | Chromosome 14 Paint | Red, green or blue | 5 | KBI-30014 |
| Whole Chromosome 15 | Chromosome 15 Paint | Red, green or blue | 5 | KBI-30015 |
| Whole Chromosome 16 | Chromosome 16 Paint | Red, green or blue | 5 | KBI-30016 |
| Whole Chromosome 17 | Chromosome 17 Paint | Red, green or blue | 5 | KBI-30017 |
| Whole Chromosome 18 | Chromosome 18 Paint | Red, green or blue | 5 | KBI-30018 |
| Whole Chromosome 19 | Chromosome 19 Paint | Red, green or blue | 5 | KBI-30019 |
| Whole Chromosome 20 | Chromosome 20 Paint | Red, green or blue | 5 | KBI-30020 |
| Whole Chromosome 21 | Chromosome 21 Paint | Red, green or blue | 5 | KBI-30021 |
| Whole Chromosome 22 | Chromosome 22 Paint | Red, green or blue | 5 | KBI-30022 |
| Whole Chromosome X | Chromosome X Paint | Red, green or blue | 5 | KBI-30023 |
| Whole Chromosome Y | Chromosome Y Paint | Red, green or blue | 5 | KBI-30024 |

* Add **-G** for Green, **-R** for Red, **-B** for Blue (available on request)

ARM SPECIFIC / BAND SPECIFIC DNA PROBES

Arm Specific Probes

In addition to our Whole Chromosome probes, we also provide the entire series of Arm Specific Probes (ASP). These probes hybridize to unique sequences comprising either p- or q-arms of all human chromosomes (except the p-arm of the acrocentric chromosomes), and they span the entire length of the respective chromosome arm. They are derived from flow-sorted or microdissected chromosomes, to be highly specific for each chromosome.

ASP applications permit the detection of chromosomal aberrations at the resolution of chromosome arms. This allows the analysis of chromosome partners involved in translocations, the identification of the chromosome of origin of marker chromosomes, analyses of complex chromosomal rearrangements in neoplastic cells and studying the inborn supernumerary marker chromosome as well as confirmation of results obtained from M-FISH and SKY testing. ASP may be of particular interest to those studying mutagenesis of human chromosomes, for instance as a result of exposure to genotoxic agents.


ASP are developed for Research Use Only (RUO) and are not meant to be used for medical purposes or as a diagnostics tool. Our Arm Specific Probes are supplied in a ready-to-use format, and are available in two colors of choice: green and red.

"Please note that only the Q arm is available for acrocentric chromosomes 13, 14, 15, 21 and 22. Heterochromatic areas (e.g. 1qh and 9qh) will not or only partially be covered by this type of probe. Due to the possibility of diffuse signals these ARM Specific Probes are not recommended for Interphase cell analysis."

Band Specific Probes

Rearrangements affecting regions smaller than an average G-band can be visualized using band-specific FISH probes. These particular probes enhance the resolution typically obtained with whole chromosome probes when identifying chromosomal abnormalities.

Band-Specific probes are capable of detecting small chromosomal segments, such as those involved in subtle translocations with breakpoints localized in distinct bands. They are amplified from microdissected chromosome material and fluorescently labeled to allow detection of these subchromosomal regions. For an updated list of Band-Specific probes please visit our website www.kreatech.com



Arm Specific Probes 2p and 6p showing a translocation of chromosome 2p on chromosome 6p.

Image kindly provided by Dr. Chantal Hamon, Paris.


Product and ordering information

| <i>Description</i> | <i>Color</i> | <i>Tests</i> | <i>Cat#</i> |
|--------------------------------|--------------|--------------|-------------|
| Arm Specific Probe 1p | green or red | 5 | KBI-30100 |
| Arm Specific Probe 1q | green or red | 5 | KBI-30101 |
| Arm Specific Probe 2p | green or red | 5 | KBI-30102 |
| Arm Specific Probe 2q | green or red | 5 | KBI-30103 |
| Arm Specific Probe 3p | green or red | 5 | KBI-30104 |
| Arm Specific Probe 3q | green or red | 5 | KBI-30105 |
| Arm Specific Probe 4p | green or red | 5 | KBI-30106 |
| Arm Specific Probe 4q | green or red | 5 | KBI-30107 |
| Arm Specific Probe 5p | green or red | 5 | KBI-30108 |
| Arm Specific Probe 5q | green or red | 5 | KBI-30109 |
| Arm Specific Probe 6p | green or red | 5 | KBI-30110 |
| Arm Specific Probe 6q | green or red | 5 | KBI-30111 |
| Arm Specific Probe 7p | green or red | 5 | KBI-30112 |
| Arm Specific Probe 7q | green or red | 5 | KBI-30113 |
| Arm Specific Probe 8p | green or red | 5 | KBI-30114 |
| Arm Specific Probe 8q | green or red | 5 | KBI-30115 |
| Arm Specific Probe 9p | green or red | 5 | KBI-30116 |
| Arm Specific Probe 9q | green or red | 5 | KBI-30117 |
| Arm Specific Probe 10p | green or red | 5 | KBI-30118 |
| Arm Specific Probe 10q | green or red | 5 | KBI-30119 |
| Arm Specific Probe 11p | green or red | 5 | KBI-30120 |
| Arm Specific Probe 11q | green or red | 5 | KBI-30121 |
| Arm Specific Probe 12p | green or red | 5 | KBI-30122 |
| Arm Specific Probe 12q | green or red | 5 | KBI-30123 |
| Arm Specific Probe 13q | green or red | 5 | KBI-30124 |
| Arm Specific Probe 14q | green or red | 5 | KBI-30125 |
| Arm Specific Probe 15q | green or red | 5 | KBI-30126 |
| Arm Specific Probe 16p | green or red | 5 | KBI-30127 |
| Arm Specific Probe 16q | green or red | 5 | KBI-30128 |
| Arm Specific Probe 17p | green or red | 5 | KBI-30129 |
| Arm Specific Probe 17q | green or red | 5 | KBI-30130 |
| Arm Specific Probe 18p | green or red | 5 | KBI-30131 |
| Arm Specific Probe 18q | green or red | 5 | KBI-30132 |
| Arm Specific Probe 19p | green or red | 5 | KBI-30133 |
| Arm Specific Probe 19q | green or red | 5 | KBI-30134 |
| Arm Specific Probe 20p | green or red | 5 | KBI-30135 |
| Arm Specific Probe 20q | green or red | 5 | KBI-30136 |
| Arm Specific Probe 21q | green or red | 5 | KBI-30137 |
| Arm Specific Probe 22q | green or red | 5 | KBI-30138 |
| Arm Specific Probe Xp | green or red | 5 | KBI-30139 |
| Arm Specific Probe Xq | green or red | 5 | KBI-30140 |
| Arm Specific Probe Yq | green or red | 5 | KBI-30141 |
| | | | |
| Band Specific Probes - inquire | green or red | 20 | KBI-302xx |


*Add -G for Green, -R for Red

MOUSE DNA PROBES


New applications for mouse molecular cytogenetics are becoming apparent. Such applications include the definition of transgene integration sites in epigenetics studies, the characterization of the mouse genome as a result of increasingly sophisticated techniques for its engineering and the screening for cytogenetic abnormalities in cell lines, such as embryonic stem cells. FISH mapping in mouse is complicated by the relative difficulty of mouse chromosome identification (karyotyping) by laboratories that are not accustomed to mouse chromosome banding techniques. Fluorescent karyotyping can be made easier by the use of counterstains or multiple reference probes but each approach needs specialized equipment and experience.



Cat.# KBI-30501




Cat.# KBI-30502



Cat.# KBI-30503

Kreatech has developed a couple of mouse region-specific probes prepared from defined BAC clones which are labeled with Platinum*Bright* dyes. These mouse probes are commonly used when developing model systems for studies in the field of genetics, mutagenesis, developmental cancer biology and also as human surrogates for studying the effects of genotoxic agents.



Cat.# KBI-30500

Hybridization of the All Mouse Centromere probe on mouse metaphase chromosomes.

Product and ordering information

| Description | Color | Tests | Cat# |
|-----------------------------|--------------|-------|-----------|
| All Mouse Centromere (AMC) | red or green | 10 | KBI-30500 |
| TK (11qE1) / AurKa (2qH3) | red/green | 10 | KBI-30501 |
| TK (11qE1) / WC Y | red/green | 10 | KBI-30502 |
| RAB9B (XqF1) / DSCR (16qC4) | red/green | 10 | KBI-30503 |
| RAB9B (XqF1) / WC Y | red/green | 10 | KBI-30505 |

Further detailed information and the availability of these probes with other labels or other combinations are available on request.

Literature:

Sabhnani et al, 2011, Repr. BioMed. Onl, 22: 621-631.
Torchia et al, 2012, J of Inv. Derm, 133: 78-86.

CELL CULTURE MEDIA / ACCESSORIES

KREAVITAL CYTOGENETIC MEDIA

Our product line KREAvital and related products have been carefully selected to optimally address cytogenetic applications. They have been further optimized for best performance and undergo strict quality procedures ensuring consistent performance and superior results.

KREAvital Prenatal Medium (Complete)

The *in vitro* cultivation of amniotic fluid cells and chorionic villi is an essential part of every diagnostic cytogenetics laboratory, since the preparation of metaphase chromosome spreads is dependent upon obtaining cells in division. Amniocentesis and chorionic villi sampling are the major invasive diagnostic procedures used for the detection of fetal chromosomal abnormalities. KREAvital Prenatal Medium is specifically optimized for the primary culture of human amniotic fluid cells and chorionic villi samples used in prenatal diagnostic testing. The medium contains serum, glutamine and antibiotics, and greatly reduces karyotyping time compared to conventional media.

In addition to the complete ready-to-use medium, we also provide KREAvital Prenatal as basal medium and supplement as separate components.

KREAvital Prenatal Medium PLUS (Complete)


Growing cells from amniotic fluids yield a mixture of epithelial cells, fibroblasts and amniocytes. KREAvital Prenatal Medium PLUS is an optimized formulation specifically enriching for fibroblasts and amniocytes, which are the cells best suited for genetic analysis. This circumstance ensures optimal chromosome morphology and metaphase structure for microscopic observation, while reducing the relatively high background generated from having too much epithelial cells.

KREAvital Prenatal Medium PLUS provides

- Less epithelial cells
- Clearer chromosome morphology for optimal banding analysis
- Enhanced buffering capacity in closed systems
- Prolonged stability when stored at 4°C

KREAvital Lymphocyte Karyotyping Medium

KREAvital Lymphocyte Karyotyping Medium is intended for use in short-term cultivation of peripheral blood lymphocytes for chromosome evaluation. It is based on RPMI-1640 basal medium supplemented with L-Glutamine, fetal bovine serum and antibiotics. We provide this medium with and without the addition of phytohaemagglutinin (PHA), respectively.




Human Amniocytes

KREAvital Bone Marrow Karyotyping Medium

An increasing number of cytogenetic analyses are carried out using bone marrow aspirates for studying chromosomal abnormalities in hematology. KREAvital Bone Marrow Karyotyping Medium is intended for use in short-term cultivation of primary bone marrow cells for chromosome evaluation and has been optimized for providing a high mitotic index. It is based on RPMI-1640 basal medium supplemented with L-Glutamine, fetal bovine serum, and antibiotics. The medium does not contain any mitogens or conditioned medium.

KREAvital Myeloid Cell Medium

Fresh cells or cells grown in short-term cultures often yield an insufficient number of mitotic cells and repeated aspirations are required. KREAvital Myeloid Cell Medium was developed to stimulate the proliferation of human hematopoietic cells from bone marrow as well as peripheral blood. This medium is particularly effective for karyotyping of acute non-lymphocytic leukemias and various stages of chronic myelogenous leukemia as well as other hematological disorders such as myelodysplastic syndrome and polycythemia vera. KREAvital Myeloid Cell Medium is based on MEM-Alpha basal medium supplemented with L-Glutamine, fetal bovine serum, antibiotics and conditioned medium.



Lymphocytes

All KREAvital media include L-Glutamine and antibiotics. To secure rigorous quality assurance, each KREAvital batch is tested for cell growth in a leading clinical cytogenetics laboratory.

Colchicine

Colchicine is a secondary metabolite originally extracted from plants of the genus *Colchicum*. In cell biology it is traditionally used as a mitosis inhibitor to arrest cells in metaphase and allowing cell harvest and karyotyping to be performed.

Colcemid

Colcemid, also known as demecolcine, is related to colchicine, but less toxic, and therefore increasingly popular as an equivalent mitosis inhibitor in cytogenetics.

Phytohaemagglutinin

Phytohaemagglutinin (PHA) is a lectin found in plants, especially beans. It has a number of physiological effects and is used in medical research to trigger cell division in T-lymphocytes. In this function, PHA is the most commonly used agent to induce mitosis in nondividing cells, such as lymphocytes and mature cells. Our PHA is provided as the mucoprotein form (PHA-M). It is supplied sterile, in a lyophilized form for constituting 5 ml of solution.



Product and ordering information

| Description | contents | Cat# |
|--|----------|-----------|
| KREAvital Prenatal Medium (Basal) | 90ml | KBI-90010 |
| KREAvital Prenatal Medium (Basal) | 450ml | KBI-92010 |
| KREAvital Prenatal Medium (Supplement) | 10ml | KBI-90011 |
| KREAvital Prenatal Medium (Supplement) | 50ml | KBI-92011 |
| KREAvital Prenatal Medium (Complete) | 100ml | KBI-90012 |
| KREAvital Prenatal Medium (Complete) | 500ml | KBI-92012 |
| KREAvital Prenatal Medium PLUS (Complete) | 100ml | KBI-90013 |
| KREAvital Prenatal Medium PLUS (Complete) | 500ml | KBI-92013 |
| KREAvital Lymphocyte Karyotyping Medium (without PHA) | 100ml | KBI-90020 |
| KREAvital Lymphocyte Karyotyping Medium (without PHA) | 500ml | KBI-92020 |
| KREAvital Lymphocyte Karyotyping Medium (including PHA) | 100ml | KBI-90021 |
| KREAvital Lymphocyte Karyotyping Medium (including PHA) | 500ml | KBI-92021 |
| KREAvital Bone Marrow Karyotyping Medium | 100ml | KBI-90030 |
| KREAvital Bone Marrow Karyotyping Medium | 500ml | KBI-92030 |
| KREAvital Myeloid Cell Medium | 100ml | KBI-90031 |
| KREAvital Myeloid Cell Medium | 500ml | KBI-92031 |
| Accessories | | |
| Colchicine Solution (10µg/ml, in PBS) | 25ml | KBI-90050 |
| Colcemid Solution (10µg/ml, in PBS) | 10ml | KBI-90051 |
| Potassium Chloride (0.075M) | 100ml | KBI-90052 |
| Phytohaemagglutinin M-Form | 5ml | KBI-90053 |
| Sodium Citrate Solution (0.8%) | 500ml | KBI-90054 |
| Trypsin EDTA 10X (EDTA 0.2%, Trypsin 0.5%, in saline solution) | 20ml | KBI-90055 |
| Trypsin EDTA 10X (EDTA 0.2%, Trypsin 0.5%, in saline solution) | 100ml | KBI-92055 |

PRETREATMENT KITS / REAGENTS

The ready-to-use pretreatment kits are recommended to be used with POSEIDON™ FISH DNA probes to produce high quality results. The kits contain all necessary ready-to-use reagents used for slide pretreatment and washing steps for FISH. The reagents increase the permeabilization of the cell membranes to facilitate penetration of the POSEIDON FISH DNA probes. The pretreatment kits will allow to process up to 25 slides in batches of 5 slides per experiment.

Pretreatment Kits

FISH Reagent Kit (KBI-60005)

FISH Reagent Kit contains all necessary ready-to-use reagents to be used for basic pretreatment of freshly prepared cytological samples.

FISH Digestion Kit (KBI-60006)

FISH Digestion Kit consists of ready-to-use reagents designed to obtain optimal results with older/difficult cytological samples or samples with cytoplasmic background which have been fixed in alcohol based fixatives (e.g. Carnoy's). Reagents provided allow to perform mild digestion on difficult cytogenetic samples, such as uncultured amniocytes, direct blood smears, buccal scrapings, urine, touch preps and others.

Tissue Digestion Kit I (KBI-60007)

Tissue Digestion Kit I contains all necessary ready-to-use reagents to be used for pretreatment of conventional paraffin-embedded tissues.

Tissue Digestion Kit II (KBI-60004)

Tissue Digestion Kit II contains all necessary ready-to-use reagents to be used for pretreatment of heavily cross-linked paraffin-embedded tissue. The Tissue Digestion Kit II provides a more intense pretreatment for optimal performance.

The selection table can guide you which pretreatment kit is the most optimal for your sample type:

| Sample | KBI-60005 | KBI-60006 | KBI-60007 | KBI-60004 |
|---|-----------|-----------|-----------|-----------|
| Freshly prepared cytological samples | + | (+) | | |
| Older / difficult / cytological samples | (+) | + | | |
| Standard tissue | | | + | (+) |
| Heavily cross-linked tissue | | | (+) | + |

Reagents

FISH Hybridization Buffer

FISH Hybridization Buffer (FHB) is a ready-to-use hybridization solution used for ULS™ labeled FISH DNA probes. FHB contains formamide, SSC, and Dextran Sulfate. Qualified for use for standard cytological and all kind of direct samples (uncultured amniocytes, blood, bone marrow, buccal scrapin, urine etc.). FHB is provided with all POSEIDON Satellite Enumeration and Sub-Telomeric probes to obtain a ready-to-use format.

Paraffin Tissue Buffer

Paraffin Tissue Buffer (PTB) is a ready-to-use hybridization solution for ULS labeled FISH DNA probes. PTB contains formamide, SSC, and Dextran Sulfate. Qualified for use on Paraffin tissues.

Whole Chromosome Buffer

Whole Chromosome Buffer (WCB) is a ready- to- use hybridization solution for ULS labeled Whole Chromosome probes. WCB contains formamide, SSC, and Dextran Sulfate. Qualified for use on standard cytological samples. WCB is provided with all POSEIDON Whole Chromosome probes.

Rubber Cement, Fixogum

For use in the formation of an air tight seal around the perimeter of the glass coverslip during probe hybridization for FISH.

Product and ordering information

| Description | Contents | Cat# |
|---------------------------------------|------------|-----------|
| POSEIDON™ Tissue Digestion Kit II | 5x5 slides | KBI-60004 |
| POSEIDON FISH Reagent Kit | 5x5 slides | KBI-60005 |
| POSEIDON FISH Digestion Kit | 5x5 slides | KBI-60006 |
| POSEIDON Tissue Digestion Kit I | 5x5 slides | KBI-60007 |
| Rubber Cement, Fixogum | 125 ml | LK-071A |
| DAPI Counterstain (0.1µg/ml) | 1 ml | LK-095A |
| DAPI Counterstain (1µg/ml) | 1 ml | LK-096A |
| Counterstain Diluent | 1 ml | LK-097A |
| Pepsin Solution | 2.5 ml | LK-101A |
| Wash Buffer I (0.4 x SSC/0.3% Igepal) | 100 ml | LK-102A |
| Wash Buffer II (2 x SSC/0.1% Igepal) | 100 ml | LK-103A |
| FISH Hybridization Buffer (FHB) | 100 µl | KBI-FHB |
| Paraffin Tissue Buffer (PTB) | 100 µl | KBI-PTB |
| Whole Chromosome Buffer (WCB) | 50 µl | KBI-WCB |

EQUIPMENT

ThermoBrite™

This programmable system automates the denaturation and hybridization steps in slide-based in situ hybridization procedures, and provides walk-away convenience for clinical and research personnel. The ThermoBrite accepts a wide range of sample types, is easy to use, and reduces hands-on time by more than 50% while ensuring overall precision and accuracy in all FISH assays.



Thermobrite
Slide Denaturation
Hybridization System

The ThermoBrite accommodates up to 12 slides and maintains uniform temperature across all slide positions. The lid seals tightly when closed providing optimal chamber humidity. The numeric keypad allows for easy programming with 40 user programmable settings and 3 modes of operation; denaturation/hybridization, hybridization, and fixed temperature.

ThermoBrite and CytoFuge 2 are trademarks of StatSpin, a division of IRIS Sample Processing.

CytoFuge® 2

The Cytofuge 2 is a low cost personal cytocentrifuge. It has simply understandable controls and new snap-seal Filter Concentrators that make operation easy. Samples are processed quickly, silently, and conveniently on any bench or in a safety cabinet. Results are consistent: easy to scan monolayer cell presentations of excellent morphologic detail.



The CytoFuge 2
is ideal for a complete
range of body fluids:

- Pleural fluids
- Sputum
- Urine
- Synovial fluids
- Fine needle aspirates
- Bronchial washings
- Cerebrospinal fluids

Specifications

| | | |
|-------------------------|--|--|
| Safety features: | Cover interlock system Leak-resistant rotor with transparent lid Low voltage drive system (24V) Designed to meet requirements of IEC-1010-2-020 | |
| Rotor capacity: | One to four slides | |
| Speed: | 600 RPM (20 x g) 700 RPM (27x g) 850 RPM (40 x g) 1000 RPM (55 x g) | 1300 RPM (93 x g) 1600 RPM (140 x g) 2200 RPM (265 x g) 4400 RPM (1060 x g) |
| Dimensions: | 15,2 cm x 16,8 cm x 21,8 cm | |
| Weight: | 2,5 kg | |

Product and ordering information


| Description | | Cat# |
|-------------------------------------|-----|--------|
| ThermoBrite (120V, 50-60Hz) | | TS-01 |
| ThermoBrite (240V, 50-60Hz) | | TS-02 |
| Humidity Control Cards | 10 | HC-10 |
| Cytofuge2 (100 - 240V, 50 / 60 Hz)* | | CF-02 |
| Reusable Filter Concentrators | 20 | FFR1 |
| Filter Concentrators (disposable) | 192 | FF01-B |

*) not available in all countries. Please inquire


THE ULS LABELING TECHNOLOGY

The proprietary ULS™ (Universal Linkage System) technology provides the basis for KREATECH's broad range of labeling applications. ULS labeling is based on the stable coordinative binding properties of platinum to nucleic acids and proteins. The ULS molecule consists of a platinum complex, a detectable molecule and a leaving group which is displaced upon reaction with the target. This reaction results in a coordinative bond, firmly coupling the ULS to the target. ULS labels DNA and RNA by binding to the N7 position of guanine. In proteins, ULS binds to sulfur and nitrogen containing side chains of methionine, cysteine and histidine. ULS is coupled to a variety of fluorophores and haptens.

Principle of the ULS technology



Optimal template amount for using Kreatech's ULS labeling kits



Targets labeled with ULS

| Nucleic Acids | Proteins |
|---|---|
| Guanine | Methionine, Cysteine, Histidine |
| ULS labels DNA and RNA by forming a coordinative bond on the N7 position of guanine | ULS labels proteins by forming a coordinative bond on the sulfur atoms of methionine, cysteine and the nitrogen atom of histidine |

Unique features of the ULS Labeling Technology

- Time saving: Entire labeling procedure in only 30 minutes for DNA and 15 minutes for RNA
- Extremely easy to use
- Robust and reproducible – no enzymes involved
- Superior labeling technology for FFPE samples
- One technology for labeling DNA, RNA and proteins
- Compatible with all types of microarray platforms.

Table of Contents

| <i>Description</i> | <i>Cat#</i> | <i>Page</i> |
|--|-------------|-------------|
| ArrayCGH | | |
| ULS arrayCGH Labeling Kit (with Cy3 and Cy5) | EA-005 | 107 |
| ULS arrayCGH Labeling Kit (with Cy3) | EA-005A | 107 |
| ULS arrayCGH Labeling Kit (with Cy5) | EA-005B | 107 |
| DNA products | | |
| ArrayGrade KREAcot DNA | EA-020 | 108 |
| ArrayGrade KREAcot DNA | EA-035 | 108 |
| Megapool Reference DNA (male) | EA-100M | 108 |
| Megapool Reference DNA (female) | EA-100F | 108 |
| Gene expression | | |
| ULS aRNA Labeling Kit (with Cy3 and Cy5) | EA-006 | 109 |
| ULS aRNA Labeling Kit (with Biotin for Affymetrix Genechips) | EA-010 | 109 |
| ULS aRNA Labeling Kit (with Biotin) | EA-018 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3 and Cy5) | EA-021 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy5) | EA-022 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3) | EA-023 | 109 |
| ULS Labeling Kit for CombiMatrix arrays (with Cy5) | EA-025 | 109 |
| ULS Labeling Kit for CombiMatrix arrays (with Biotin) | EA-027 | 109 |
| MicroRNA | | |
| ULS™ microRNA Labeling Kit (with Cy3 and Cy5) | EA-036 | 110 |
| ULS microRNA Labeling Kit (with Cy3) | EA-037 | 110 |
| ULS microRNA Labeling Kit (with Cy5) | EA-038 | 110 |
| General nucleic acid labeling | | |
| PlatinumBright™ Nucleic Acid Labeling Kit (495 Green) | GLK-001 | 112 |
| PlatinumBright Nucleic Acid Labeling Kit (547 Light Red) | GLK-002 | 112 |
| PlatinumBright Nucleic Acid Labeling Kit (647 Far Red) | GLK-003 | 112 |
| PlatinumBright Nucleic Acid Labeling Kit (550 Red) | GLK-004 | 112 |
| PlatinumBright Nucleic Acid Labeling Kit (415 Blue) | GLK-006 | 112 |
| PlatinumBright Nucleic Acid Labeling Kit (Biotin) | GLK-007 | 112 |
| FISHBright™ Labeling Kit (415 Blue) | FLK-001 | 112 |
| FISHBright Labeling Kit (495 Green) | FLK-002 | 112 |
| FISHBright Labeling Kit (505 Green) | FLK-003 | 112 |
| FISHBright Labeling Kit (550 Red) | FLK-004 | 112 |
| FISHBright Labeling Kit (547 Light Red) | FLK-005 | 112 |
| FISHBright Labeling Kit (647 Far Red) | FLK-006 | 112 |
| FISHBright Labeling Kit (Biotin) | FLK-007 | 112 |
| FISH Grade CoT | KB-COT | 112 |
| Protein labeling / Antibody labeling | | |
| PlatinumLink Protein Labeling Kit (BIO) | PLK-007 | 113 |
| PlatinumLink Protein Labeling Kit (FLU) | PLK-009 | 113 |
| PlatinumLink Protein Labeling Kit (RHO) | PLK-010 | 113 |
| Supporting products | | |
| KREApure™ columns | KP-020 | 114 |
| KREApure columns | KP-050 | 114 |
| KREApure 96* | KP-096 | 114 |

ARRAYCGH


ULS™ arrayCGH Labeling Kit

The ULS arrayCGH Labeling Kit offers a novel procedure that allows direct (non-enzymatic) labeling of both intact genomic DNA as well as fragmented genomic DNA isolated from formalin-fixed paraffin-embedded (FFPE) samples. The ULS arrayCGH Labeling Kit yields highly reproducible fluorescent labeled DNA within minutes. The ULS technology has been validated on a wide range of oligo and BAC arrayCGH platforms. All ULS kits for arrayCGH analysis include the KREApure purification technology and KREAblock blocking reagents to ensure efficient purification of labeled DNA samples and low background levels during hybridization.

Unique features of ULS arrayCGH Labeling Kit


- 30-minute labeling procedure
- Labeling independent of fragment size - ideal for FFPE samples
- ULS labeling is not affected by cross links present in genomic DNA from FFPE samples
- No enzymatic bias
- KREApure column purification ideally suited for fragmented DNA
- No bias when amplified DNA need to be labeled.

ULS Labeling vs Random Prime Labeling



Ratio plots of 3 independent array CGH hybridizations using genomic DNA from healthy liver which is the reference vs. genomic DNA isolated from a liver tumor FFPE sample. (A-C) ULS labeled samples. (D-F) Random Prime labeled samples.

Time of ULS labeling vs. enzymatic labeling procedure



Product and ordering information

| Description | Contents | Cat# |
|--|----------------------------|---------|
| arrayCGH Labeling | | |
| ULS arrayCGH Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 20 µg DNA | EA-005 |
| ULS arrayCGH Labeling Kit (with Cy3) | for labeling 40 µg DNA | EA-005A |
| ULS arrayCGH Labeling Kit (with Cy5) | for labeling 40 µg DNA | EA-005B |

DNA PRODUCTS


ArrayGrade KREAcot DNA

ArrayGrade KREAcot DNA is extracted from human placental DNA and subsequently fragmented, denatured, and re-annealed under conditions that enrich for repetitive DNA sequences (1,2). ArrayGrade KREAcot DNA can be used to suppress cross-hybridization to human repetitive DNA sequences. ArrayGrade KREAcot DNA has specifically been optimized to be used in array CGH applications. Amounts of ArrayGrade KREAcot DNA needed in a microarray experiment should be determined empirically, but will be in the order of 12.5-50 times excess to the amount of labeled genomic DNA. In addition, ArrayGrade KREAcot is suitable for blocking repetitive sequences during the target enrichment procedure for high throughput sequencing purposes. Please visit the Kreatech website (www.kreatech.com) for a detailed protocol for the analysis of genomic DNA using arrayCGH.

Unique features of ArrayGrade KREAcot

- Specifically optimized for use in arrayCGH applications
- Quality controlled with arrayCGH analysis.

Comparison of performance of ArrayGrade KREAcot DNA and human C₀t-1 DNA from competitor X.



ArrayGrade KREAcot

²log ratio for X chromosome = -0.65
Mean absolute deviation for chrom 1-22 = 0.05

Human C₀t-1 DNA competitor X

²log ratio for X chromosome = -0.66
Mean absolute deviation for chrom 1-22 = 0.11

Shown are the chromosome plots obtained from an arrayCGH hybridization comparing healthy male Cy5-ULS labeled genomic DNA and healthy female Cy3-ULS labeled genomic DNA on BAC microarrays. Although the spread observed for the X-chromosome is comparable for both blockers, the use of ArrayGrade KREAcot minimizes the variation observed for chromosomes 1-22.


Megapool Reference DNA

This product is a homogeneous DNA pool from male or female human genomic DNA which has been isolated from 100 different anonymous healthy individuals. The genomic DNA is of high quality and the DNA from each of these 100 individuals contributes evenly to the DNA pool. This product is specifically developed as a reference for genomic microarray-based comparative genomic hybridization experiments (array CGH).

Unique features of the MegaPool Reference DNA


- The product will guarantee perfectly comparable profiles for male and female individuals as well as complete absence of pure homozygous deletions
- Sex-matched arrayCGH experiments will be possible and enable high-quality genotyping data for both chromosomes X and Y as well as unbiased CNV profiles between male and female test samples.

Megapool Reference DNA (Female vs Male)



Left: Female Megapool Reference DNA is hybridized against the Male Megapool Reference DNA on an Agilent's 4 x 180 k human aCGH array. The X and Y chromosomes can be clearly identified. The Derivative Log-Ratio Spread (DLRS) value is **0.084** which is considered as excellent.

Male patient vs. Male Megapool Reference DNA



Right: Male Megapool Reference DNA is hybridized against a male mental retardation patient on an Agilent's 4 x 180k human aCGH array. Variations are clearly visible at Chromosomes 5, 15 and 16. The Derivative Log-Ratio Spread (DLRS) value is **0.1189**. Data are kindly provided by Bauke Ylstra and colleagues VU University Medical Center, Amsterdam.

Product and ordering information

| Description | Contents | Cat# |
|---------------------------------|----------|---------|
| ArrayGrade KREAcot DNA | 500 µg | EA-020 |
| ArrayGrade KREAcot DNA | 10 mg | EA-035 |
| Megapool Reference DNA (male) | 200 µg | EA-100M |
| Megapool Reference DNA (female) | 200 µg | EA-100F |

GENE EXPRESSION

Amplification and aRNA Labeling Kits


aRNA Labeling Kits

The Universal Linkage System (ULS™) allows direct labeling of unmodified, amplified RNA (aRNA). Specific aRNA labeling kits are both available for the use with commercial platforms from Agilent, Affymetrix® (Genechips®) and CombiMatrix, and for the use with self-spotted DNA arrays.

Unique features of aRNA Labeling Kits


- Higher yields of unmodified aRNA compared to modified aRNA
- Storage of unmodified aRNA for subsequent use possible
- No waste of labeled material; label only the amount aRNA needed
- Reduced labeling cost per hybridization when using multiple arrays per slide.

Higher yields of unmodified aRNA as compared to modified aRNA



Unmodified antisense RNA (aRNA) vs. modified aRNA.

Lower cost if multiple arrays per slide are used



Costs for ULS labeling of unmodified aRNA vs. costs for enzymatic labeling generating modified aRNA.

Product and ordering information

| Description | Contents | Cat# |
|--|-----------------------------|--------|
| aRNA Labeling | | |
| ULS aRNA Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 50 µg aRNA | EA-006 |
| ULS aRNA Labeling Kit (with Biotin for Affymetrix Genechips) | for labeling 500 µg aRNA | EA-010 |
| ULS aRNA Labeling Kit (with Biotin) | for labeling 250 µg aRNA | EA-018 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3 and Cy5) | for labeling 2 x 50 µg aRNA | EA-021 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy5) | for labeling 50 µg aRNA | EA-022 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3) | for labeling 50 µg aRNA | EA-023 |
| ULS Labeling Kit for CombiMatrix arrays (with Cy5) | for labeling 125 µg aRNA | EA-025 |
| ULS Labeling Kit for CombiMatrix arrays (with Biotin) | for labeling 125 µg aRNA | EA-027 |

MICRORNA



microRNA Labeling Kit

The new ULS microRNA Labeling Kit provides a very simple and fast method for miRNA detection on microarray platforms. The Universal Linkage System (ULS) labeling technology is used to directly label the naturally occurring miRNAs which are present in the total RNA mixture without the need to enrich for small RNA. ULS can be used to label total RNA from all kind of organisms.

Unique features of the ULS microRNA Labeling Kit

- 15-minute labeling step
- Easy two-step procedure to label total RNA
- Compatible with most existing miRNA platforms
- Not size discriminative.

Time of ULS labeling vs. competitive labeling technologies.



Product and ordering information

| Description | Contents | Cat# |
|--|----------------------------|--------|
| microRNA Labeling | | |
| ULS microRNA Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 25 µg RNA | EA-036 |
| ULS microRNA Labeling Kit (with Cy3) | for labeling 50 µg RNA | EA-037 |
| ULS microRNA Labeling Kit (with Cy5) | for labeling 50 µg RNA | EA-038 |

GENERAL NUCLEIC ACID LABELING


PLATINUMBright™ Labeling Kits

The PlatinumBright Nucleic Acid Labeling Kits are based on KREATECH's patented Universal Linkage System (ULS). The ULS molecule is a proprietary platinum compound linked to a hapten or a fluorophore. The PlatinumBright kits offer ULS coupled to a variety of fluorophores or biotin providing a non-enzymatic, easy-to-use alternative to conventional techniques such as nick translation, random priming, end labeling, etc.


The labeled nucleic acid can be used for the downstream application of your choice. In addition to the PlatinumBright Nucleic Acid Labeling Kits, KREATECH also offers the FISHBright kits which, in addition to reagents for labeling, include all reagents needed for FISH hybridization.

Unique features of the PlatinumBright Labeling Kits

- Only a two-step, very fast labeling procedure
- Easy to perform
- A variety of labels to choose from
- Not dependent on the performance of enzymes.



Centromere 11 labeled with PlatinumBright 415.



4-color FISH.

FISHBright™ labeling kits

The FISHBright labeling kits are based on KREATECH's Universal Linkage System (ULS™). The ULS molecule is a proprietary platinum compound linked to a hapten or a fluorophore. The FISHBright kits offer ULS molecules coupled to a variety of fluorophores or haptens providing a non-enzymatic, easy-to-use alternative to conventional techniques such as nick translation, random priming, end labeling, etc.

The FISHBright labeling kits are specially developed for FISH applications: labeled DNA probes can be used on all type of samples including metaphase spreads, direct interphase cell preparation (e.g. blood smears, bone marrow smears, urine samples, paraffin-embedded or frozen tissue sections).

KREATECH's REPEAT-FREE™ POSEIDON™ DNA probes are labeled with the same technology and all POSEIDON pretreatment kits can be used in combination with DNA probes which are labeled with the FISHBright labeling kits.

Besides the ULS technology, the kits also contain Cell specific Hybridization Buffer (CHB), and Tissue Hybridization Buffer (THB). In addition, the KREAboost solution, which is developed to generate a high signal-to-noise ratio especially on paraffin-embedded tissue sections, is included in these kits.

The FISHBright labeling kits will offer you all advantages of the ULS technology:

- Easy to use
- 30-minute labeling procedure
- Use of KREApure columns: a unique purification column which guarantees a very high removal of non-reacted ULS (> 99%) and high sample recoveries (> 95%).

Content:

- PlatinumBright component (fluorochromes or haptens)
- 10 x labeling solution
- KREApure purification columns
- Hybridization buffer for cytological preps (CHB)
- Hybridization buffer for paraffin preps (THB+KREAboost).

The FISHBright kit contains sufficient solution to label up to 10 µg of BAC DNA, which generally allows for 100 or more FISH assays using 10 µl labeled probe in hybridization buffer (amount of FISH assays is dependent on the final concentration of labeled DNA. Usually 50 – 100 ng/µl of labeled BAC DNA is necessary per assay).

FISHgrade C₀t

FISHgrade C₀t DNA is extracted from human placental DNA treated under conditions that enrich for repetitive DNA sequences (1,2). These conditions have been optimized by Kretech Diagnostics for the specific use in FISH experiments.

FISHgrade C₀t DNA suppresses cross-hybridization to human repetitive DNA. Amounts of FISHgrade C₀t DNA needed in a FISH experiment should be determined empirically, but will be in the order of 12.5 – 50 times excess to the amount of labeled genomic DNA.


Unique features of FISHgrade C₀t

- High quality human C₀t-1 DNA optimized for use in FISH application
- Specific quality control with FISH analysis
- Works perfectly for probes labeled with FISHBright ULS-Kits.

References












1. Weiner, A.M., et al. (1986) Ann. Rev. Biochem. 55, 631.
2. Britten, R.J., et al. (1974) Methods Enzymol. 29, 363.

Cat.# KB-COT FISHgrade Cot



FISH image on metaphase spread performed using a classical FISH probe labeled in red/green colors using a FISHBright™ labeling kit. FISHgrade C₀t DNA has been used to block all aspecific DNA sequences. The probe specifically recognizes deletions at 20q chromosome that have been related to myelodysplastic syndromes.

Product and ordering information

| Description | Color | Excitation/Emission | Contents | Cat# |
|--|---|---------------------|-----------------------------|---------|
| PlatinumBright™ Nucleic Acid Labeling Kit (495 Green) |  | 495/517 | for labeling 20 µg template | GLK-001 |
| PlatinumBright Nucleic Acid Labeling Kit (547 Light Red) |  | 547/565 | for labeling 20 µg template | GLK-002 |
| PlatinumBright Nucleic Acid Labeling Kit (647 Far Red) |  | 647/665 | for labeling 20 µg template | GLK-003 |
| PlatinumBright Nucleic Acid Labeling Kit (550 Red) |  | 550/580 | for labeling 20 µg template | GLK-004 |
| PlatinumBright Nucleic Acid Labeling Kit (415 Blue) |  | 429/470 | for labeling 20 µg template | GLK-006 |
| PlatinumBright Nucleic Acid Labeling Kit (Biotin) | | | for labeling 20 µg template | GLK-007 |
| FISHBright™ Labeling Kit (415 Blue) |  | 429/470 | for labeling 10 µg DNA | FLK-001 |
| FISHBright Labeling Kit (495 Green) |  | 495/517 | for labeling 10 µg DNA | FLK-002 |
| FISHBright Labeling Kit (505 Green) |  | 500/528 | for labeling 10 µg DNA | FLK-003 |
| FISHBright Labeling Kit (550 Red) |  | 550/580 | for labeling 10 µg DNA | FLK-004 |
| FISHBright Labeling Kit (547 Light Red) |  | 547/565 | for labeling 10 µg DNA | FLK-005 |
| FISHBright Labeling Kit (647 Far Red) |  | 647/665 | for labeling 10 µg DNA | FLK-006 |
| FISHBright Labeling Kit (Biotin) | | | for labeling 10 µg DNA | FLK-007 |
| FISH Grade CoT | | | 500 µg | KB-COT |

PROTEIN LABELING / ANTIBODY LABELING

PlatinumLink Antibody Labeling Kits

The PlatinumLink Kits provide high coverage of the proteome through specific labeling of methionine, cysteine and histidine residues. Unlike lysine, which is targeted by N-hydroxysuccinimidyl (NHS) esters, the amino acids labeled with ULS are less likely to be involved in protein - protein interactions. Therefore, labeling of proteins and antibodies with ULS will reduce negative effects on interaction domains and epitope recognition.


The PlatinumLink Antibody Labeling Kits have been optimized to label 100 µg of purified antibody (or recombinant protein) per reaction with fluorescent dyes or haptens.

Unique features of the PlatinumLink Antibody Labeling Kits

- Easy and robust labeling of antibodies and recombinant proteins
- ULS-Trap removes unreacted label, no size exclusion
- Labeling of methionines, cysteines and histidines in antibodies and recombinant proteins
- Labeling reaction is independent of pH
- Compatible with commonly-used detergents including Tris-HCl.

Schematic overview of the PlatinumLink labeling procedure.

Step 1: Protein labeling Step 2: Removal of unreacted ULS-label



Recovery comparison of different proteins.




Fig.18 Protein recovery using ULS-Trap columns. 100 µg of protein was labeled and purified using PlatinumLink protein labeling kits. Protein concentration was determined using Biorad's Rc Dc protein assay. Recovery is presented as percentage protein recovered as compared to unpurified sample.

Product and ordering information


| Description | Reactions | Cat# |
|---|--------------------|---------|
| PlatinumLink Protein Labeling Kit (BIO) | 4 single labelings | PLK-007 |
| PlatinumLink Protein Labeling Kit (FLU) | 4 single labelings | PLK-009 |
| PlatinumLink Protein Labeling Kit (RHO) | 4 single labelings | PLK-010 |

SUPPORTING PRODUCTS


KREApure™ purification columns and plates

KREApure purification columns and plates are based on a proprietary purification matrix optimized for purifying ULS-labeled target nucleic acids. The procedure quantitatively removes non-attached label and provides optimal recovery of labeled nucleic acids.

KREApure spin column purification.



Yield comparison of different purification methods.




Cy3 – ULS-labeled DNA has been purified using five different purification systems. Experiments were performed in triplicate.

Unique features of the KREApure purification

- Maximal removal of unreacted ULS-label
- > 95% recovery of labeled nucleic acid
- Efficient recovery of small fragments

KREApure 96-well procedure.



Product and ordering information

| Description | Reactions | Cat# |
|-------------------|-----------|--------|
| KREApure™ columns | 20 pcs | KP-020 |
| KREApure columns | 50 pcs | KP-050 |
| KREApure 96* | 1 plate | KP-096 |

* specific high-throughput labeling modules available upon request

PRODUCT INDEX

CATALOGUE 2013-2014 - NEW PRODUCTS

ON = Oncology probes

| Description | Color | Tests | Cat# | Page |
|---|----------------|-------|-----------|------|
| ONCOLOGY - HEMATOLOGY DNA PROBES | | | | |
| CML ON Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | green/red | 10 | KBI-10013 | II |
| AML ON MECOM / RUNX1 t(3;21) Fusion | red/green | 10 | KBI-10310 | II |
| ON NUP98 (11p15) Break | red/green | 10 | KBI-10311 | III |
| ALL ON MLL/AFF1 t(4;11) Fusion | red/green | 10 | KBI-10404 | III |
| MM ON MAFB/IGH@ t(14;20) Fusion | red/green | 10 | KBI-10510 | IV |
| ONCOLOGY – SOLID TUMOR DNA PROBES | | | | |
| Lung Cancer ON ROS1 (6q22) Break | green/red | 10 | KBI-10752 | V |
| ON RET (10q11) Break | green/red | 10 | KBI-10753 | V |
| ON FGFR1 (8p11) / SE 8 (D8Z1) | red/green | 20 | KBI-12754 | VI |
| ON FGFR1 (8p11) / SE 8 (D8Z1) | red/green | 50 | KBI-14754 | VI |
| Sarcoma ON COL1A1/PDGFB t(17;22) DC, S-Fusion | red/green | 10 | KBI-10742 | VI |
| ON EWSR1/NFATC2 t(20;22) DC, S-Fusion | red/green | 10 | KBI-10751 | VII |
| Lymphoma (tissue) ON MYC (8q24) Break, TC (tissue) | red/green/blue | 10 | KBI-10749 | VII |
| ON BCL2/IGH@ t(14;18) Fusion (tissue) | red/green | 10 | KBI-10755 | VIII |
| MICRODELETION DNA PROBES | | | | |
| MD GATA4 (8p23) / SE 8 | red/green | 10 | KBI-40118 | VIII |
| MD GATA4 (8p23) / SE 8 | red/green | 5 | KBI-45118 | VIII |

ONCOLOGY – HEMATOLOGY DNA PROBES

| Description | Color | Tests | Cat# | Page |
|---|----------------|-------|-----------|-------|
| CML ON FIP1L1-CHIC2-PDGFRα (4q12) Del, Break | red/green | 10 | KBI-10003 | 14 |
| ON PDGFRB (5q33) Break | red/green | 10 | KBI-10004 | 15 |
| ON BCR/ABL t(9;22) Fusion | red/green | 10 | KBI-10005 | 11,35 |
| ON BCR/ABL t(9;22) Fusion | red/green | 20 | KBI-12005 | 11,35 |
| ON BCR/ABL t(9;22) TC, D-Fusion | red/green/blue | 10 | KBI-10006 | 11,35 |
| ON FIP1L1-CHIC2-PDGFRα (4q12) Del, Break, TC | red/green/blue | 10 | KBI-10007 | 15 |
| ON BCR/ABL t(9;22), DC, S-Fusion, ES | red/green | 10 | KBI-10008 | 12,35 |
| ON BCR/ABL t(9;22) DC, S-Fusion | red/green | 10 | KBI-10009 | 12,35 |
| ON p53 (17p13) / MPO (17q22) "ISO 17q" | green/red | 10 | KBI-10011 | 13 |
| ON JAK2 (9p24) Break | green/red | 10 | KBI-10012 | 13 |
| CLL ON DLEU (13q14) / 13qter | red/green | 10 | KBI-10102 | 17 |
| ON ATM (11q22) / SE 11 | red/green | 10 | KBI-10103 | 18 |
| ON GLI (12q13) / SE 12 | red/green | 10 | KBI-10104 | 20 |
| ON 6q21 / SE 6 | red/green | 10 | KBI-10105 | 18 |
| ON C-MYC (8q24) / SE 8 | red/green | 10 | KBI-10106 | 19 |
| ON ATM (11q22) / GLI (12q13) | red/green | 10 | KBI-10108 | 21 |
| ON 6q21 / MYC (8q24) | red/green | 10 | KBI-10117 | 21 |
| ON hTERC (3q26) / 3q11 | red/green | 10 | KBI-10110 | 19 |
| ON p53 (17p13) / SE 17 | red/green | 10 | KBI-10112 | 17 |
| ON p53 (17p13) / SE 17 | red/green | 20 | KBI-12112 | 17 |
| ON DLEU (13q14) / p53 (17p13) | red/green | 10 | KBI-10113 | 20 |
| ON p53 (17p13) / ATM (11q22) | green/red | 10 | KBI-10114 | 22 |

| | Description | Color | Tests | Cat# | Page |
|-------------------------|---|----------------|-------|-----------|----------|
| MDS | ON MDS 7q- (7q22; 7q36) | green/red | 10 | KBI-10202 | 25 |
| | ON MDS 20q- (PTPRT 20q12) / 20q11 | red/green | 10 | KBI-10203 | 26 |
| | ON EVI t(3;3); inv(3) (3q26) Break | red/green | 10 | KBI-10204 | 26 |
| | ON EVI t(3;3); inv(3) (3q26) Break, TC | red/green/blue | 10 | KBI-10205 | 27 |
| | ON MDS 7q- (7q22; 7q36) / SE 7 TC | green/red/blue | 10 | KBI-10207 | 25 |
| | ON hTERT (5p15) / 5q31 | red/green | 10 | KBI-10208 | 23 |
| | ON MDS 5q- (5q31; 5q33) | green/red | 10 | KBI-10209 | 24 |
| | ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC | green/red/blue | 10 | KBI-10210 | 24 |
| AML | ON AML/ETO t(8;21) Fusion | green/red | 10 | KBI-10301 | 28 |
| | ON PML/RARA t(15;17) Fusion | green/red | 10 | KBI-10302 | 29 |
| | ON PML/RARA t(15;17) Fusion | green/red | 20 | KBI-12302 | 29 |
| | ON MLL (11q23) Break | red/green | 10 | KBI-10303 | 29,35 |
| | ON CBFB t(16;16), inv(16) Break | red/green | 10 | KBI-10304 | 32 |
| | ON RARA (17q21) Break | red/green | 10 | KBI-10305 | 31 |
| | ON DEK / NUP214 t(6;9) Fusion | red/green | 10 | KBI-10306 | 32 |
| | ON MLL/MLLT1 t(11;19) Fusion | red/green | 10 | KBI-10307 | 30 |
| | ON MLL/MLLT3 t(9;11) Fusion | red/green | 10 | KBI-10308 | 30 |
| | ON MLL/MLLT4 t(6;11) Fusion | red/green | 10 | KBI-10309 | 31 |
| ALL | ON TEL/AML t(12;21) Fusion | red/green | 10 | KBI-10401 | 33 |
| | ON p16 (9p21) / 9q21 | red/green | 10 | KBI-10402 | 34 |
| | ON ETV6 (TEL) (12p13) Break | red/green | 10 | KBI-10403 | 34 |
| Multiple Myeloma | ON MM 11q23 / DLEU (13q14) | green/red | 10 | KBI-10502 | 38 |
| | ON MM 1q21 / 8p21 | green/red | 10 | KBI-10503 | 39 |
| | ON MM 15q22 / 6q21 | green/red | 10 | KBI-10504 | 40 |
| | ON MM 1q21 / SRD (1p36) | green/red | 10 | KBI-10507 | 41 |
| | ON MM 15q22 / 9q34 | green/red | 10 | KBI-10508 | 40 |
| | ON MM 19q13 / p53 (17p13) | green/red | 10 | KBI-10509 | 39 |
| Lymphoma related probes | ON IGH (14q32) Break | red/green | 10 | KBI-10601 | 22,41,44 |
| | ON FGFR3/IGH t(4;14) Fusion | green/red | 10 | KBI-10602 | 38,47 |
| | ON MYC/IGH t(8;14) Fusion | green/red | 10 | KBI-10603 | 42 |
| | ON BCL1/IGH t(11;14) Fusion | green/red | 10 | KBI-10604 | 43 |
| | ON MYEOV/IGH t(11;14) Fusion | green/red | 10 | KBI-10605 | 36,47 |
| | ON BCL2/IGH t(14;18) Fusion | green/red | 10 | KBI-10606 | 43 |
| | ON BCL6 (3q27) Break | red/green | 10 | KBI-10607 | 45 |
| | ON MALT (18q21) Break | red/green | 10 | KBI-10608 | 45 |
| | ON CCND1 (BCL1;11q13) Break | red/green | 10 | KBI-10609 | 46 |
| | ON MAF/IGH t(14;16) Fusion | green/red | 10 | KBI-10610 | 37 |
| | ON MYC (8q24) Break, TC | red/green/blue | 10 | KBI-10611 | 47 |
| | ON BCL2 (18q21) Break | red/green | 10 | KBI-10612 | 44 |

ONCOLOGY – SOLID TUMOR DNA PROBES

ON = Oncology probes

| Description | Color | Tests | Cat# | Page |
|---|----------------|-------|-----------|------|
| ON ERBB2, Her-2/Neu (17q12) / SE 17 | red/green | 10 | KBI-10701 | 49 |
| ON ERBB2, Her-2/Neu (17q12) / SE 17 | red/green | 50 | KBI-14701 | 49 |
| ON EGFR, Her-1 (7p11) / SE 7 | red/green | 10 | KBI-10702 | 54 |
| ON CC hTERT (3q26) C-MYC (8q24) / SE 7 TC | red/green/blue | 10 | KBI-10704 | 52 |
| ON MYCN (2p24) / LAF (2q11) | red/green | 10 | KBI-10706 | 59 |
| ON PPARγ (3p25) Break | red/green | 10 | KBI-10707 | 58 |
| ON hTERT (5p15) / 5q31 (tissue) | red/green | 10 | KBI-10709 | 55 |
| ON p16 (9p21) / 9q21 (tissue) | red/green | 10 | KBI-10710 | 52 |
| ON MLL (11q23) / SE 11 | red/green | 10 | KBI-10711 | 60 |
| ON SRD (1p36) / SE 1(1qh) | red/green | 10 | KBI-10712 | 60 |
| ON SYT (18q11) Break | red/green | 10 | KBI-10713 | 63 |
| ON CHOP (12q13) Break | red/green | 10 | KBI-10714 | 63 |
| ON FUS (16p11) Break | red/green | 10 | KBI-10715 | 64 |
| ON FKHR (13q14) Break | red/green | 10 | KBI-10716 | 64 |
| ON MDM2 (12q15) / SE 12 | red/green | 10 | KBI-10717 | 65 |

| <i>Description</i> | <i>Color</i> | <i>Tests</i> | <i>Cat#</i> | <i>Page</i> |
|---|----------------|--------------|-------------|-------------|
| ON PTEN (10q23) / SE 10 | red/green | 10 | KBI-10718 | 57 |
| ON C-MET (7q31) / SE 7 | red/green | 10 | KBI-10719 | 55 |
| ON AR (Xq12) / SE X | red/green | 10 | KBI-10720 | 56 |
| ON AURKA (20q13) / 20q11 | red/green | 10 | KBI-10721 | 67 |
| ON AURKB (17p13) / SE 17 | red/green | 10 | KBI-10722 | 68 |
| ON TOP2A (17q21) / SE 17 | red/green | 10 | KBI-10724 | 49 |
| ON CDK4 (12q13) / SE 12 | red/green | 10 | KBI-10725 | 66 |
| ON TMPRSS2-ERG (21q22) Del, Break, TC | red/green/blue | 10 | KBI-10726 | 57 |
| ON IGH (14q32) Break (tissue) | red/green | 10 | KBI-10729 | 70 |
| ON BCL6 (3q27) Break (tissue) | red/green | 10 | KBI-10730 | 71 |
| ON MALT (18q21) Break (tissue) | red/green | 10 | KBI-10731 | 71 |
| ON ZNF217 (20q13) / 20q11 | red/green | 10 | KBI-10733 | 50 |
| ON CCND1 (11q13) / SE 11 | red/green | 10 | KBI-10734 | 69 |
| ON TOP2A (17q21) / Her-2/neu (17q12) / SE 17 TC | red/green/blue | 10 | KBI-10735 | 50 |
| ON MDM4 (1q32) /SE 1 | red/green | 10 | KBI-10736 | 61 |
| ON FGFR1 (8p12) Break | red/green | 10 | KBI-10737 | 16,51 |
| ON p53 (17p13) / SE 17 (tissue) | red/green | 10 | KBI-10738 | 72 |
| ON ERCC1 (19q13) / ZNF443 (19p13) | red/green | 10 | KBI-10739 | 66 |
| ON TFE3 (Xp11) Break | red/green | 10 | KBI-10741 | 72 |
| ON BCL2 (18q21) Break (tissue) | red/green | 10 | KBI-10745 | 70 |
| ON ALK/EML4 t(2;2); inv(2) Fusion | red/green | 10 | KBI-10746 | 54 |
| ON ALK (2p23) Break | red/green | 10 | KBI-10747 | 53,73 |
| ON EWSR1 (22q12) Break | red/green | 10 | KBI-10750 | 62 |

ONCOLOGY - CHROMOGENIC IN SITU HYBRIDIZATION

| <i>Description</i> | <i>Tests</i> | <i>Cat#</i> | <i>Page</i> |
|-----------------------------|--------------|-------------|-------------|
| UniStar CISH Detection Kit | 10 | KBI-50001 | 74 |
| UniStar Her2/neu (17q12) | 10 | KBI-50701 | 74 |
| UniStar EGFR (7p11) | 10 | KBI-50702 | 74 |
| UniStar C-MET (7q31) | 10 | KBI-50719 | 74 |
| TwinStar CISH Detection Kit | 10 | KBI-60010 | 74 |
| TwinStar Her2/neu (17q12) | 10 | KBI-60701 | 74 |
| TwinStar EGFR (7p11) | 10 | KBI-60702 | 74 |
| TwinStar C-MET (7q31) | 10 | KBI-60719 | 74 |

PREIMPLANTATION GENETIC SCREENING

| <i>Description</i> | <i>Color</i> | <i>Tests</i> | <i>Cat#</i> | <i>Page</i> |
|------------------------------------|---------------------|--------------|-------------|-------------|
| PreimpScreen PolB (13,16,18,21,22) | five color | 20 | KBI-40050 | 77 |
| PreimpScreen Blas (13,18,21,X,Y) | five color | 20 | KBI-40051 | 77 |
| MultiStar 24 FISH | six color, 4 panels | 10 | KBI-40060 | 77 |
| MultiStar FISH Panel 1 | six color | 10 | KBI-40061 | 77 |
| MultiStar FISH Panel 2 | six color | 10 | KBI-40062 | 77 |
| MultiStar FISH Panel 3 | six color | 10 | KBI-40063 | 77 |
| MultiStar FISH Panel 4 | six color | 10 | KBI-40064 | 77 |

PRENATAL DNA PROBES

PN = Prenatal probes

| Description | Color | Tests | Cat# | Page |
|-----------------------------|----------------|-------|-------------|------|
| PN 13 (13q14) | green | 10 | KBI-40001 | 78 |
| PN 21 (21q22) | red | 10 | KBI-40002 | 78 |
| PN 13 (13q14) / 21 (21q22) | green/red | 10 | KBI-40003 | 79 |
| SE 18 (D18Z1) 5x conc | blue | 10 | KBI-20018-B | 79 |
| SE X (DXZ1) / SE Y (DYZ3) | green/red | 10 | KBI-20030 | 79 |
| SE 7 (D7Z1) / SE 8 (D8Z1) | red/green | 10 | KBI-20031 | 79 |
| SE (X,Y,18) | green/red/blue | 10 | KBI-20032 | 79 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 10 | KBI-40005 | 79 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 30 | KBI-40006 | 79 |
| PrenatScreen (13/21, X/Y18) | green/red/blue | 50 | KBI-40007 | 79 |
| PloidyScreen (21, X, Y) | red/green/blue | 20 | KBI-40008 | 79 |

MICRODELETION DNA PROBES

MD = Microdeletion probes

| Description | Color | Tests | Cat# | Page |
|--|----------------|-------|-----------|-------|
| MD Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) | red/green | 10 | KBI-40101 | 87 |
| MD DiGeorge "N25" (22q11) / 22q13 (SHANK3) | red/green | 10 | KBI-40102 | 81,83 |
| MD DiGeorge Tuple (22q11) / 22q13 (SHANK3) | red/green | 10 | KBI-40103 | 81,83 |
| MD DiGeorge T-box1 (22q11) / 22q13 (SHANK3) | red/green | 10 | KBI-40104 | 82,83 |
| MD DiGeorge II (10p14) / SE 10 | red/green | 10 | KBI-40105 | 83 |
| MD Cri-Du-Chat CTNND (5p15) / 5q31 | red/green | 10 | KBI-40106 | 88 |
| MD Wolf-Hirschhorn WHSC1 (4p16) / SE 4 | red/green | 10 | KBI-40107 | 88 |
| MD X-Inactivation XIST (Xq13) / SE X | red/green | 10 | KBI-40108 | 89 |
| MD Prader-Willi SNRPN (15q11) / PML(15q24) | red/green | 10 | KBI-40109 | 85 |
| MD Angelman UBE3A (15q11) / PML(15q24) | red/green | 10 | KBI-40110 | 86 |
| MD Williams-Beuren ELN (7q11) / 7q22 | red/green | 10 | KBI-40111 | 86 |
| MD Short Stature (Xp22) / SE X | red/green | 10 | KBI-40112 | 89 |
| MD NSD1 (5q35)/ hTERT (5p15) | red/green | 10 | KBI-40113 | 84 |
| MD NF1 (17q11) / MPO (17q22) | red/green | 10 | KBI-40114 | 84 |
| MD STS (Xp22) / KAL (Xp22) / SE X TC | red/green/blue | 10 | KBI-40115 | 90 |
| MD IGF1R (15q26) / 15q11 | red/green | 10 | KBI-40116 | 91 |

SATELLITE ENUMERATION DNA PROBES

SE = Satellite Enumeration probes
5x conc format

| Description | Color | Tests | Cat# | Page |
|--------------------------------|---------------------|-------|-----------|-------|
| All Human Centromer (AHC), RTU | green or red | 10 | KBI-20000 | 94 |
| SE 1 (1qh) | green, red, or blue | 10 | KBI-20001 | 94 |
| SE 2 (D2Z) | green, red, or blue | 10 | KBI-20002 | 94 |
| SE 3 (D3Z1) | green, red, or blue | 10 | KBI-20003 | 94 |
| SE 4 (D4Z1) | green, red, or blue | 10 | KBI-20004 | 94 |
| SE 6 (D6Z1) | green, red, or blue | 10 | KBI-20006 | 94 |
| SE 7 (D7Z1) | green, red, or blue | 10 | KBI-20007 | 94 |
| SE 8 (D8Z1) | green, red, or blue | 10 | KBI-20008 | 14,94 |
| SE 9 (classical) | green, red, or blue | 10 | KBI-20009 | 94 |
| SE 10 (D10Z1) | green, red, or blue | 10 | KBI-20010 | 94 |
| SE 11 (D11Z1) | green, red, or blue | 10 | KBI-20011 | 94 |
| SE 12 (D12Z3) | green, red, or blue | 10 | KBI-20012 | 22,94 |
| SE 15 (D15Z) | green, red, or blue | 10 | KBI-20015 | 94 |
| SE 16 (D16Z2) | green, red, or blue | 10 | KBI-20016 | 94 |
| SE 17 (D17Z1) | green, red, or blue | 10 | KBI-20017 | 94 |
| SE 18 (D18Z1) | green, red, or blue | 10 | KBI-20018 | 94 |
| SE 20 (D20Z1) | green, red, or blue | 10 | KBI-20020 | 94 |
| SE X (DXZ1) | green, red, or blue | 10 | KBI-20023 | 94 |
| SE Y (DYZ3) | green, red, or blue | 10 | KBI-20024 | 94 |
| SE Y classical q arm | green, red, or blue | 10 | KBI-20025 | 94 |
| SE 1/5/19 | green, red, or blue | 10 | KBI-20026 | 94 |
| SE 13/21 | green, red, or blue | 10 | KBI-20027 | 94 |
| SE 14/22 | green, red, or blue | 10 | KBI-20028 | 94 |
| SE X (DXZ1) / SE Y (DYZ3) | green/red | 10 | KBI-20030 | 94 |
| SE 7 (D7Z1) / SE 8 (D8Z1) | red/green | 10 | KBI-20031 | 14,94 |
| SE (X,Y,18) | green/red/blue | 10 | KBI-20032 | 94 |
| Acro-P-Arms NOR | green, red, or blue | 10 | KBI-20033 | 95 |

* Add -G for Green, -R for Red, -B for Blue (available on request)

SUB-TELOMERE DNA PROBES

5x conc format

| Description | Color | Tests | Cat# | Page |
|---------------------|---------------------|-------|-----------|------|
| Sub Telomere 1pter | green, red, or blue | 5 | KBI-40201 | 92 |
| Sub Telomere 1qter | green, red, or blue | 5 | KBI-40202 | 92 |
| Sub Telomere 2pter | green, red, or blue | 5 | KBI-40203 | 92 |
| Sub Telomere 2qter | green, red, or blue | 5 | KBI-40204 | 92 |
| Sub Telomere 3pter | green, red, or blue | 5 | KBI-40205 | 92 |
| Sub Telomere 3qter | green, red, or blue | 5 | KBI-40206 | 92 |
| Sub Telomere 4pter | green, red, or blue | 5 | KBI-40207 | 92 |
| Sub Telomere 4qter | green, red, or blue | 5 | KBI-40208 | 92 |
| Sub Telomere 5pter | green, red, or blue | 5 | KBI-40209 | 92 |
| Sub Telomere 5qter | green, red, or blue | 5 | KBI-40210 | 92 |
| Sub Telomere 6pter | green, red, or blue | 5 | KBI-40211 | 92 |
| Sub Telomere 6qter | green, red, or blue | 5 | KBI-40212 | 92 |
| Sub Telomere 7pter | green, red, or blue | 5 | KBI-40213 | 92 |
| Sub Telomere 7qter | green, red, or blue | 5 | KBI-40214 | 92 |
| Sub Telomere 8pter | green, red, or blue | 5 | KBI-40215 | 92 |
| Sub Telomere 8qter | green, red, or blue | 5 | KBI-40216 | 92 |
| Sub Telomere 9pter | green, red, or blue | 5 | KBI-40217 | 92 |
| Sub Telomere 9qter | green, red, or blue | 5 | KBI-40218 | 92 |
| Sub Telomere 10pter | green, red, or blue | 5 | KBI-40219 | 92 |
| Sub Telomere 10qter | green, red, or blue | 5 | KBI-40220 | 92 |
| Sub Telomere 11pter | green, red, or blue | 5 | KBI-40221 | 92 |
| Sub Telomere 11qter | green, red, or blue | 5 | KBI-40222 | 92 |
| Sub Telomere 12pter | green, red, or blue | 5 | KBI-40223 | 92 |
| Sub Telomere 12qter | green, red, or blue | 5 | KBI-40224 | 92 |
| Sub Telomere 13qter | green, red, or blue | 5 | KBI-40225 | 92 |
| Sub Telomere 14qter | green, red, or blue | 5 | KBI-40226 | 92 |
| Sub Telomere 15qter | green, red, or blue | 5 | KBI-40227 | 92 |
| Sub Telomere 16pter | green, red, or blue | 5 | KBI-40228 | 92 |
| Sub Telomere 16qter | green, red, or blue | 5 | KBI-40229 | 92 |
| Sub Telomere 17pter | green, red, or blue | 5 | KBI-40230 | 92 |
| Sub Telomere 17qter | green, red, or blue | 5 | KBI-40231 | 92 |
| Sub Telomere 18pter | green, red, or blue | 5 | KBI-40232 | 92 |
| Sub Telomere 18qter | green, red, or blue | 5 | KBI-40233 | 92 |
| Sub Telomere 19pter | green, red, or blue | 5 | KBI-40234 | 92 |
| Sub Telomere 19qter | green, red, or blue | 5 | KBI-40235 | 92 |
| Sub Telomere 20pter | green, red, or blue | 5 | KBI-40236 | 92 |
| Sub Telomere 20qter | green, red, or blue | 5 | KBI-40237 | 92 |
| Sub Telomere 21qter | green, red, or blue | 5 | KBI-40238 | 92 |
| Sub Telomere 22qter | green, red, or blue | 5 | KBI-40239 | 92 |
| Sub Telomere XYpter | green, red, or blue | 5 | KBI-40240 | 92 |
| Sub Telomere XYqter | green, red, or blue | 5 | KBI-40241 | 92 |

* Add -G for Green, -R for Red, -B for Blue (available on request)

ARM SPECIFIC / BAND SPECIFIC DNA PROBES

ready-to-use

| <i>Description</i> | <i>Color</i> | <i>Tests</i> | <i>Cat#</i> | <i>Page</i> |
|--------------------------------|--------------|--------------|-------------|-------------|
| Arm Specific Probe 1p | green or red | 5 | KBI-30100 | 99 |
| Arm Specific Probe 1q | green or red | 5 | KBI-30101 | 99 |
| Arm Specific Probe 2p | green or red | 5 | KBI-30102 | 99 |
| Arm Specific Probe 2q | green or red | 5 | KBI-30103 | 99 |
| Arm Specific Probe 3p | green or red | 5 | KBI-30104 | 99 |
| Arm Specific Probe 3q | green or red | 5 | KBI-30105 | 99 |
| Arm Specific Probe 4p | green or red | 5 | KBI-30106 | 99 |
| Arm Specific Probe 4q | green or red | 5 | KBI-30107 | 99 |
| Arm Specific Probe 5p | green or red | 5 | KBI-30108 | 99 |
| Arm Specific Probe 5q | green or red | 5 | KBI-30109 | 99 |
| Arm Specific Probe 6p | green or red | 5 | KBI-30110 | 99 |
| Arm Specific Probe 6q | green or red | 5 | KBI-30111 | 99 |
| Arm Specific Probe 7p | green or red | 5 | KBI-30112 | 99 |
| Arm Specific Probe 7q | green or red | 5 | KBI-30113 | 99 |
| Arm Specific Probe 8p | green or red | 5 | KBI-30114 | 99 |
| Arm Specific Probe 8q | green or red | 5 | KBI-30115 | 99 |
| Arm Specific Probe 9p | green or red | 5 | KBI-30116 | 99 |
| Arm Specific Probe 9q | green or red | 5 | KBI-30117 | 99 |
| Arm Specific Probe 10p | green or red | 5 | KBI-30118 | 99 |
| Arm Specific Probe 10q | green or red | 5 | KBI-30119 | 99 |
| Arm Specific Probe 11p | green or red | 5 | KBI-30120 | 99 |
| Arm Specific Probe 11q | green or red | 5 | KBI-30121 | 99 |
| Arm Specific Probe 12p | green or red | 5 | KBI-30122 | 99 |
| Arm Specific Probe 12q | green or red | 5 | KBI-30123 | 99 |
| Arm Specific Probe 13q | green or red | 5 | KBI-30124 | 99 |
| Arm Specific Probe 14q | green or red | 5 | KBI-30125 | 99 |
| Arm Specific Probe 15q | green or red | 5 | KBI-30126 | 99 |
| Arm Specific Probe 16p | green or red | 5 | KBI-30127 | 99 |
| Arm Specific Probe 16q | green or red | 5 | KBI-30128 | 99 |
| Arm Specific Probe 17p | green or red | 5 | KBI-30129 | 99 |
| Arm Specific Probe 17q | green or red | 5 | KBI-30130 | 99 |
| Arm Specific Probe 18p | green or red | 5 | KBI-30131 | 99 |
| Arm Specific Probe 18q | green or red | 5 | KBI-30132 | 99 |
| Arm Specific Probe 19p | green or red | 5 | KBI-30133 | 99 |
| Arm Specific Probe 19q | green or red | 5 | KBI-30134 | 99 |
| Arm Specific Probe 20p | green or red | 5 | KBI-30135 | 99 |
| Arm Specific Probe 20q | green or red | 5 | KBI-30136 | 99 |
| Arm Specific Probe 21q | green or red | 5 | KBI-30137 | 99 |
| Arm Specific Probe 22q | green or red | 5 | KBI-30138 | 99 |
| Arm Specific Probe Xp | green or red | 5 | KBI-30139 | 99 |
| Arm Specific Probe Xq | green or red | 5 | KBI-30140 | 99 |
| Arm Specific Probe Yq | green or red | 5 | KBI-30141 | 99 |
| Band Specific Probes - inquire | green or red | 20 | KBI-302xx | 99 |

*Add -G for Green, -R for Red

WHOLE CHROMOSOME DNA PROBES

5x conc format

| Description | Color | Tests | Cat# | Page |
|---------------------|---------------------|-------|-----------|------|
| Whole Chromosome 1 | green, red, or blue | 5 | KBI-30001 | 97 |
| Whole Chromosome 2 | green, red, or blue | 5 | KBI-30002 | 97 |
| Whole Chromosome 3 | green, red, or blue | 5 | KBI-30003 | 97 |
| Whole Chromosome 4 | green, red, or blue | 5 | KBI-30004 | 97 |
| Whole Chromosome 5 | green, red, or blue | 5 | KBI-30005 | 97 |
| Whole Chromosome 6 | green, red, or blue | 5 | KBI-30006 | 97 |
| Whole Chromosome 7 | green, red, or blue | 5 | KBI-30007 | 97 |
| Whole Chromosome 8 | green, red, or blue | 5 | KBI-30008 | 97 |
| Whole Chromosome 9 | green, red, or blue | 5 | KBI-30009 | 97 |
| Whole Chromosome 10 | green, red, or blue | 5 | KBI-30010 | 97 |
| Whole Chromosome 11 | green, red, or blue | 5 | KBI-30011 | 97 |
| Whole Chromosome 12 | green, red, or blue | 5 | KBI-30012 | 97 |
| Whole Chromosome 13 | green, red, or blue | 5 | KBI-30013 | 97 |
| Whole Chromosome 14 | green, red, or blue | 5 | KBI-30014 | 97 |
| Whole Chromosome 15 | green, red, or blue | 5 | KBI-30015 | 97 |
| Whole Chromosome 16 | green, red, or blue | 5 | KBI-30016 | 97 |
| Whole Chromosome 17 | green, red, or blue | 5 | KBI-30017 | 97 |
| Whole Chromosome 18 | green, red, or blue | 5 | KBI-30018 | 97 |
| Whole Chromosome 19 | green, red, or blue | 5 | KBI-30019 | 97 |
| Whole Chromosome 20 | green, red, or blue | 5 | KBI-30020 | 97 |
| Whole Chromosome 21 | green, red, or blue | 5 | KBI-30021 | 97 |
| Whole Chromosome 22 | green, red, or blue | 5 | KBI-30022 | 97 |
| Whole Chromosome X | green, red, or blue | 5 | KBI-30023 | 97 |
| Whole Chromosome Y | green, red, or blue | 5 | KBI-30024 | 97 |

* Add -G for Green, -R for Red, -B for Blue (available on request)

MOUSE DNA PROBES

| Description | Color | Tests | Cat# | Page |
|-----------------------------|--------------|-------|-----------|------|
| All Mouse Centromere (AMC) | red or green | 10 | KBI-30500 | 100 |
| TK (11qE1) / AurKa (2qH3) | red/green | 10 | KBI-30501 | 100 |
| TK (11qE1) / WC Y | red/green | 10 | KBI-30502 | 100 |
| RAB9B (XqF1) / DSCR (16qC4) | red/green | 10 | KBI-30503 | 100 |
| RAB9B (XqF1) / WC Y | red/green | 10 | KBI-30505 | 100 |

* Add -G for Green, -R for Red

CELL CULTURE MEDIA / ACCESSORIES

| <i>Description</i> | <i>Contents</i> | <i>Cat#</i> | <i>Page</i> |
|--|-----------------|-------------|-------------|
| KREAvital Prenatal Medium (Basal) | 90ml | KBI-90010 | 102 |
| KREAvital Prenatal Medium (Basal) | 450ml | KBI-92010 | 102 |
| KREAvital Prenatal Medium (Supplement) | 10ml | KBI-90011 | 102 |
| KREAvital Prenatal Medium (Supplement) | 50ml | KBI-92011 | 102 |
| KREAvital Prenatal Medium (Complete) | 100ml | KBI-90012 | 102 |
| KREAvital Prenatal Medium (Complete) | 500ml | KBI-92012 | 102 |
| KREAvital Prenatal Medium PLUS (Complete) | 100ml | KBI-90013 | 102 |
| KREAvital Prenatal Medium PLUS (Complete) | 500ml | KBI-92013 | 102 |
| KREAvital Lymphocyte Karyotyping Medium (without PHA) | 100ml | KBI-90020 | 102 |
| KREAvital Lymphocyte Karyotyping Medium (without PHA) | 500ml | KBI-92020 | 102 |
| KREAvital Lymphocyte Karyotyping Medium (including PHA) | 100ml | KBI-90021 | 102 |
| KREAvital Lymphocyte Karyotyping Medium (including PHA) | 500ml | KBI-92021 | 102 |
| KREAvital Bone Marrow Karyotyping Medium | 100ml | KBI-90030 | 102 |
| KREAvital Bone Marrow Karyotyping Medium | 500ml | KBI-92030 | 102 |
| KREAvital Myeloid Cell Medium | 100ml | KBI-90031 | 102 |
| KREAvital Myeloid Cell Medium | 500ml | KBI-92031 | 102 |
| Colchicine Solution (10µg/ml, in PBS) | 25ml | KBI-90050 | 102 |
| Colcemid Solution (10µg/ml, in PBS) | 10ml | KBI-90051 | 102 |
| Potassium Chloride (0.075M) | 100ml | KBI-90052 | 102 |
| Phytohaemagglutinin M-Form | 5ml | KBI-90053 | 102 |
| Sodium Citrate Solution (0.8%) | 500ml | KBI-90054 | 102 |
| Trypsin EDTA 10X (EDTA 0.2%, Trypsin 0.5%, in saline solution) | 20ml | KBI-90055 | 102 |
| Trypsin EDTA 10X (EDTA 0.2%, Trypsin 0.5%, in saline solution) | 100ml | KBI-92055 | 102 |

PRETREATMENT KITS / REAGENTS

| <i>Description</i> | <i>Contents</i> | <i>Cat#</i> | <i>Page</i> |
|---------------------------------------|-----------------|-------------|-------------|
| POSEIDON™ Tissue Digestion Kit II | 5x5 slides | KBI-60004 | 103 |
| POSEIDON FISH Reagent Kit | 5x5 slides | KBI-60005 | 103 |
| POSEIDON FISH Digestion Kit | 5x5 slides | KBI-60006 | 103 |
| POSEIDON Tissue Digestion Kit I | 5x5 slides | KBI-60007 | 103 |
| Rubber Cement, Fixogum | 125 ml | LK-071A | 103 |
| DAPI Counterstain (0.1µg/ml) | 1 ml | LK-095A | 103 |
| DAPI Counterstain (1µg/ml) | 1 ml | LK-096A | 103 |
| Counterstain Diluent | 1 ml | LK-097A | 103 |
| Pepsin Solution | 2.5 ml | LK-101A | 103 |
| Wash Buffer I (0.4 x SSC/0.3% Igepal) | 100 ml | LK-102A | 103 |
| Wash Buffer II (2 x SSC/0.1% Igepal) | 100 ml | LK-103A | 103 |
| FISH Hybridization Buffer (FHB) | 100 µl | KBI-FHB | 103 |
| Paraffin Tissue Buffer (PTB) | 100 µl | KBI-PTB | 103 |
| Whole Chromosome Buffer (WCB) | 50 µl | KBI-WCB | 103 |

EQUIPMENT

| <i>Description</i> | <i>Contents</i> | <i>Cat#</i> | <i>Page</i> |
|-------------------------------------|-----------------|-------------|-------------|
| ThermoBrite™ (240V, 50 - 60Hz) | | TS-02 | 104 |
| Humidity Control Cards | 10 | HC-10 | 104 |
| Cytofuge®2 (100 - 240V, 50 / 60 Hz) | | CF-02 | 104 |

ARRAYCGH

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|---|----------------------------|-------------|-------------|
| ULS™ arrayCGH Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 20 µg DNA | EA-005 | 107 |
| ULS arrayCGH Labeling Kit (with Cy3) | for labeling 40 µg DNA | EA-005A | 107 |
| ULS arrayCGH Labeling Kit (with Cy5) | for labeling 40 µg DNA | EA-005B | 107 |

DNA PRODUCTS

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|--------------------------------|------------------|-------------|-------------|
| Array-Grade KREAcot DNA | 500 µg | EA-020 | 108 |
| Array-Grade KREAcot DNA | 10 mg | EA-035 | 108 |
| Megapool Reference DNA, male | 200 µg | EA-100M | 108 |
| Megapool Reference DNA, female | 200 µg | EA-100F | 108 |

GENE EXPRESSION

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|--|-----------------------------|-------------|-------------|
| ULS aRNA Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 50 µg aRNA | EA-006 | 109 |
| ULS aRNA Labeling Kit (with Biotin for Affymetrix® Genechips®) | for labeling 500 µg aRNA | EA-010 | 109 |
| ULS aRNA Labeling Kit (with Biotin) | for labeling 250 µg aRNA | EA-018 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3 and Cy5) | for labeling 2 x 50 µg aRNA | EA-021 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy5) | for labeling 50 µg aRNA | EA-022 | 109 |
| ULS Fluorescent Labeling Kit for Agilent arrays (with Cy3) | for labeling 50 µg aRNA | EA-023 | 109 |
| ULS Labeling Kit for CombiMatrix arrays (with Cy5) | for labeling 125 µg aRNA | EA-025 | 109 |
| ULS Labeling Kit for CombiMatrix arrays (with Biotin) | for labeling 125 µg aRNA | EA-027 | 109 |
| (20 amplifications each) | | | |

MICRO RNA

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|---|----------------------------|-------------|-------------|
| ULS™ microRNA Labeling Kit (with Cy3 and Cy5) | for labeling 2 x 25 µg RNA | EA-036 | 110 |
| ULS microRNA Labeling Kit (with Cy3) | for labeling 50 µg RNA | EA-037 | 110 |
| ULS microRNA Labeling Kit (with Cy5) | for labeling 50 µg RNA | EA-038 | 110 |

GENERAL NUCLEIC ACID LABELING

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|---|-----------------------------|-------------|-------------|
| Platinum <i>Bright</i> [™] Nucleic Acid Labeling Kit (495 Green) | for labeling 20 µg template | GLK-001 | 111 |
| Platinum <i>Bright</i> Nucleic Acid Labeling Kit (547 Light Red) | for labeling 20 µg template | GLK-002 | 111 |
| Platinum <i>Bright</i> Nucleic Acid Labeling Kit (647 Far Red) | for labeling 20 µg template | GLK-003 | 111 |
| Platinum <i>Bright</i> Nucleic Acid Labeling Kit (550 Red) | for labeling 20 µg template | GLK-004 | 111 |
| Platinum <i>Bright</i> Nucleic Acid Labeling Kit (415 Blue) | for labeling 20 µg template | GLK-006 | 111 |
| Platinum <i>Bright</i> Nucleic Acid Labeling Kit (Biotin) | for labeling 20 µg template | GLK-007 | 111 |
| | | | |
| FISH <i>Bright</i> [™] 415 Blue | for labeling 10 µg DNA | FLK-001 | 111 |
| FISH <i>Bright</i> 495 Green | for labeling 10 µg DNA | FLK-002 | 111 |
| FISH <i>Bright</i> 505 Green | for labeling 10 µg DNA | FLK-003 | 111 |
| FISH <i>Bright</i> 550 Red | for labeling 10 µg DNA | FLK-004 | 111 |
| FISH <i>Bright</i> 547 Light Red | for labeling 10 µg DNA | FLK-005 | 111 |
| FISH <i>Bright</i> 647 Far Red | for labeling 10 µg DNA | FLK-006 | 111 |
| FISH <i>Bright</i> Biotin | for labeling 10 µg DNA | FLK-007 | 111 |
| | | | |
| FISH Grade CoT | 500 µg | KB-COT | 112 |

PROTEIN LABELING / ANTIBODY LABELING

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|---|--------------------|-------------|-------------|
| Platinum <i>Link</i> Protein Labeling Kit (BIO) | 4 single labelings | PLK-007 | 113 |
| Platinum <i>Link</i> Protein Labeling Kit (FLU) | 4 single labelings | PLK-009 | 113 |
| Platinum <i>Link</i> Protein Labeling Kit (RHO) | 4 single labelings | PLK-010 | 113 |

SUPPORTING PRODUCTS

| <i>Description</i> | <i>Reactions</i> | <i>Cat#</i> | <i>Page</i> |
|-------------------------------|------------------|-------------|-------------|
| KREApure [™] columns | 20 pcs | KP-020 | 114 |
| KREApure columns | 50 pcs | KP-050 | 114 |
| KREApure 96 | 1 plate | KP-096 | 114 |

CE MARKING & ISO CERTIFICATES

All probes from KREATECH are manufactured and marketed in accordance with the European IVD Directive 98/79/EC. Apart from those products designed for research use (Kreatech's ASP, BSP, and mouse FISH probes, as well as the ULS™ Labeling product line), all KREATECH products are CE marked. In addition, KREATECH has obtained certificates for CE-marking of conformity for our prenatal and preimplantation genetic screening probes.



All CE-marked FISH probes in this catalog are manufactured at KREATECH Diagnostics facilities in Amsterdam, The Netherlands, and are labeled with ULS (Universal Linkage System).

Manufacturing of all CE-marked products by KREATECH is carried out according to international standards ISO 9001:2008 and ISO 13485:2003 for which KREATECH holds certificates.

Compliance to the European IVD Directive 98/79/EC, ISO 9001:2008 and ISO 13485:2003 means that we are providing high quality in vitro diagnostic medical devices to be used for diagnosis, prevention, monitoring, treatment, alleviation or investigation of disease, congenital disorders, physiological and/or pathological state of the human body without being used in or on the human body.

We strive to continuously improve the quality of these products based on the latest developments and techniques.



TERMS & CONDITIONS

Ordering and Delivery

Orders for Kreatech products can be placed via e-mail, fax or over the telephone with one of our Sales and/or Customer Service representatives on the basis of a purchase order number. Purchase orders for regularly stocked items received by Kreatech by or before 14:00 CET Monday through Thursday are shipped the same working day. Purchase orders received after 14:00 CET are shipped the following working day. In the event that an order contains a non-stocked item, or if, for any reason, Kreatech is unable to completely fulfill an order, a partial shipment will be made unless instructed otherwise. You will be notified immediately of any delay. Delivery is made EXW Kreatech, Amsterdam (Incoterms 2010).

Prices, Payment and shipping terms

Terms of sales are net 30 days from the date shown on the invoice. Prices are subject to change without notice. All orders under € 400 will be charged with the shipping costs. There is no minimum order quantity; however Kreatech reserves the right to charge an additional processing fee of € 15 for orders received under € 200.

All orders are subject to review by Kreatech and additional information may be requested. Kreatech reserves the right to refuse shipment of material to anyone.

Standing Order / Volume Discounts

For your convenience, Kreatech can schedule routine shipments of any single purchase order. Discounted prices may be available when purchasing large volumes. Inquiries should be directed to your local Sales Representative.

Methods of payment including Credit Card orders

To better serve your needs with convenience and efficiency, Kreatech also accepts MasterCard and VISA as payment for orders placed via fax. Please provide the name on your credit card, account number, expiration date and security code when placing credit card orders.

The Kreatech Customer Satisfaction Policy

If you are not completely satisfied with the performance of any product, Kreatech will, at your option, either credit your account or send a replacement product to you at our expense. Requirements for this are that the product has not expired, it has been stored and processed according to Kreatech product literature and that the user followed all appropriate Kreatech protocols when using the product.

Returned Goods Policy

No returns will be approved without the authorization of Kreatech. Unauthorized returns will not be accepted and will be returned at the sender's expense. Complete documentation must accompany all authorized returned goods. If errors in ordering goods are made by Kreatech, Kreatech will exchange goods at Kreatech's expense.

Product Use

All products manufactured and/or distributed by Kreatech should be used in accordance with the products labeled intended use.

Kreatech's Poseidon™ Repeat-Free™ FISH products are supplied for in vitro diagnostics (or for research use only if not CE marked). They are NOT intended to be used as drugs, food, additives, cosmetics or household chemicals. Kreatech warrants that, at the time of shipment, the products sold are free from defects in materials and conform to the specifications described in Kreatech's literature. Kreatech makes no other warranty, expressed or implied, with respect to the products. This includes any warranty of merchantability or suitability for any particular purpose.

Liability

Kreatech accepts liability for damages suffered by its customers, up to a maximum amount of € 5.000,- if these damages are directly and exclusively due to gross negligence and/or willful misconduct by Kreatech.

Kreatech does not accept liability for (i) damages resulting from the use/stocking of Kreatech products by its customers contrary to Kreatech's instructions, (ii) damages resulting from off label use of the Kreatech products, (iii) the results of any diagnosis or research made on the basis of or with the Kreatech products or (iv) any consequential or indirect damages suffered by its customers.

Patents and Trademarks

These products or the use of these products is subject to intellectual property rights. FISH and CISH probes are produced under an exclusive license from Veridex LLC for its REPEAT-FREE™ technology. FISH probes are labeled with the Kreatech's proprietary Universal Linkage system (ULS™). The fluorophores used in the PlatinumBright™ Nucleic Acid Labeling Kits (-415; -547; -590; -647) are subject to patents and know-how, owned or controlled, and manufactured by DYOMICS GmbH. The fluorophores used in the ULS™Cy5 and ULS™Cy3 labeling compounds are subject to patents, owned or controlled, and manufactured by GE Healthcare. The ULS™ technology and products are covered by one or more of the following US patents and corresponding patent claims outside the US owned by KREATECH: US 5,580, 990; 5,714,327; 5,985,566; 6,133,038; US RE 40,557E, 6,797,818; 7,217,813.

KREATECH is a trade name of Kreatech Biotechnology B.V.; ULS™, PlatinumBright™, KreaPure™, KreaBlock™ and POSEIDON™ are trademarks of KREATECH Biotechnology B.V.

REPEAT-FREE™ is a trademark of Veridex LLC. ThermoBrite™ and CytoFuge™ 2 are trademarks of StatSpin, a division of IRIS Sample Processing.

Nothing disclosed herein is to be construed as a recommendation to use any of KREATECH's products in violation of any patents. Furthermore, KREATECH does not warrant that the product, either alone, or in combination with other products, is immune from charges of patent infringement. KREATECH will not be held responsible for patent infringements or other violations that may occur from the use of any of its products.

For full General Sale and Supply Terms, see www.kreatech.com

INTERNATIONAL DISTRIBUTORS

Algeria
Eurl Syphax Biotechnologie
N° 126 Quartier Benamara
Abdelkader
Haoche Chaouche El Achour, Alger
Phone: 021.60.70.91
Fax: 021.60.70.91
E-mail: syphaxbiotech@yahoo.fr

Australia
Diagnostic Technology Pty Ltd.
Suite 45, 7 Narabang Way
Belrose, NSW 2085
Phone: + 61 2 99862011
Fax: + 61 2 99862022
E-mail: info@diagnostictchnology.com.au

Brazil
Uniscience do Brasil
Av. Cândido Portinari, 933/937
05114-001 São Paulo - SP
Phone: + 55 (0)11 3622 2320
Fax: + 55 (0)11 3622 23 23
E-mail: info@uniscience.com

Bulgaria:
Elta 90m O O D
19, Dunav Str. Ent.A Fl.1
1000 Sofia
Phone: + 359 2 983 9649
Fax: + 359 2 983 2211
E-mail: elta90@dir.bg

Chile
Prolab
Vergara 24
Oficina 908 - Casilla 3645
Santiago
Fono / Fax 56-2- 698 7215
56-2- 698 9617
56-2- 695 8404
E-mail: ecerna@prolab.cl

China
You Ning Biotech Co. Ltd.
7C, Run Yue Building Block 2. No.0,
Tian He Road
Guangzhou
Phone: +86 20 37604588
Fax: +86 20 37604618
E-mail: info@youning.com

China
Will-Tek Electro-Optical Ltd.
2C, No. 1, Alley 889,
Wuzhong Road
Shanghai
Phone: 021-6406-4668
Fax: 021-5422-5188
E-mail: will-tek@will-tek.com.cn

Colombia
Andina de Tecnologías Ltda.
Calle 34 No. 81-59
Medellín
Phone: +57 4 4164040
Fax: +57 4 4164040
E-mail: biotecn@anditecnica.com

Croatia
Jasika d.o.o.
Remetinečka cesta 115,
10000 Zagreb
Phone: + 385 1 65 52 664
Fax: + 385 1 65 35 018
E-mail: ivana.tomicic@jasika.hr

Czech Republic
Intimex s.r.o.
Kavrina-Nove Mesto
Vrchlickeho sady 541/6
735 06
Tel: + 420 596 311612
Fax: + 421 596 311018
E-mail: intimex@post.cz

Denmark
AH diagnostics as
Runefoften 18, 8210 Aarhus
Phone: + 45 8745 9010
Fax: + 45 8745 1292
E-mail: ahdiag@ahdiag.dk

Ecuador
GenLife cia.Ltda
Urb.La Pampa 2, calle K, N5-314
Quito, EC170177
Phone: +593 22352992
Fax: +593 22354500
E-mail: info@gen-life.net

Egypt
Cornell Lab
3/1 Degla 5th Sector
El Maadi
Cairo
Phone: + 202 25160732
Fax: + 202 25200270
E-mail: corlab@cornell-lab.com.eg

Estonia
Quantum Eesti AS
Saekoja 36a
50107 Tartu
Phone: + 372 7301321
Fax: + 372 7304310
E-mail: quantum@quantum.ee

Finland
AH diagnostics Oy
Koetilantie 1 B 4
00790 Helsinki
Phone: + 358 9 3509100
Fax: + 358 9 35091022
E-mail: ahdiagnostics@ahdiagnostics.fi

Greece
Bioanalytica S.A.
Biotechnology, SafeBlood BioAnalytica
and Analytical Systems
3-5 Iliadion Str.
11528 Athens
Phone: +30 210 6400318
Fax: + 30 210 6462748
E-mail: l.katabelisi@bioanalytica.gr
and orders@bioanalytica.gr

Hungary
BioMarker Ltd.
2100 Gödöllő
Szölő utca 27 - 29.
Tel: +36 28 419 986
Fax +36 28 422 319
E-mail: biomarker@biomarker.hu

India
Applied Resonance Technologies
505, PUNIT TOWER-II, SECTOR-11,
CBD BELAPUR,
NAVI MUMBAI-400614
Mr. Rajeev Jaiswal
Phone: +91-98-2102 4375
E-mail: rajeev@appliedresonancetech.com

Iran
Sakan Teb
No 28, Akhavan St., Shariate Ave.
Passed Junction of Ayatollah Sadr
19316 Tehran
Phone: + 982 1 2223 3133
Fax: +982 1 2221 4773
E-mail: SakanTeb@hotmail.com

Ireland
Serosep Limited
Annacotty Business Park
Annacotty
Co. Limerick
Phone: + 353 61 35 8190
Fax: + 353 61 35 8191
E-mail: info@serosep.com

Israel
Biological Industries Ltd
Contact: Dr. Gary Bennet
Kibbutz Beit Haemek
M.P. Ashrat
ISRAEL, 25115
Phone: + 972 4 996 05 95
Fax: + 972 4 996 99 96
E-mail: gary@bioind.com

Italy
Resnova S.r.l.
Via Cadore, 14
00045 Genzano di Roma RM
Phone: + 39 (0)6 93955058
Fax: + 39 (0)6 93955059
E-mail: servizioclienti@resnovaweb.it

Italy
Teltec srl
Via Lecco 4,
20864 Agrate Brianza (MB)
Phone: + 39 039 6892171
Fax: + 39 039 633098
E-mail: info@teltec.com

Japan
Funakoshi Co. Ltd.
9-7 Hongo 2-Chome
Bunkyo-Ku
Tokyo 113-0033
Phone: + 81 3 5684 1620
Fax: + 81 3 5684 1775
E-mail: reagent@funakoshi.co.jp

Jordan
LeanGene
P.O.Box 697 Tlaa Alali
Amman
Phone: + 962 6 46 55 401
Fax: + 962 6 46 55 402
E-mail: support@leangene-jo.com

Kuwait
Tareq Co. W.L.L.
Area 1, Block 103, Ardiya,
P.O. Box 20506,
Safat 13066
Phone: + 965 243 106 15
Fax: + 965 223 130 85
E-mail: henry@tareqco.com

Latvia
Interlux, SIA
Lubānas iela 78
Rīga, LV-1073
Phone: +371 67795240
Fax: +371 67795241
E-mail: info@interlux.lv

Lithuania
Lithuania Interlux Ltd.
Avieciu 16, Vilnius 2021
Phone: + 370 5 278 68 50
Fax: + 370 5 279 67 28
E-mail: spirit@interlux.lt

Malaysia
Nano Life Quest Sdn. Bhd.
No.33, Jalan Puteri 5/10
Bandar Puteri, 47100 Puchong
Selangor D.E.
Phone: + 603 8063 2688
Fax: + 603 8063 2699
E-mail: info@nanolifequest.com

México
Uniparts S.A.
Galileo 92, Polanco
México, D.F. 11550
Phone: + 52 55 52 82 47 18
Fax: + 52 55 52 81 47 22
E-mail: benitor@uniparts.com.mx

Norway
AH diagnostics as
Fjellgata 1
N-0566 Oslo
Phone: +47 2323 3260
Fax: + 47 2323 3270
E-mail: ahdiagnostics@ahdiagnostics.no

Poland
Intimex sp. z o.o.
Krzysztof Stankiewicz
ul. Spacerowa 2
05-119 Legionowo
Phone: + 48 22 6686465
Fax: +48 22 8221147
E-mail: biuro@intimex.com.pl

Portugal
Bioportugal
Rua do Campo Alegre 1306
2º Sala 208
4150-174 Porto
Phone: + 351 22 6004800
Fax: + 351 22 6004801
E-mail: bioportugal@bioportugal.pt

Romania
Dexter Com
Str. Popa Rusu, nr. 9A, apt 6
Bucuresti
Phone: +4-021-212.23.69
Fax: +4-021-212-23.70
E-mail: dexter@itcnet.ro

Russia
Medico-Diagnostic Laboratory
Gamalai, 16
Moscow, 123098
Phone: + 499 190 28 68
Fax: + 499 190 30 47
E-mail: meddialab@mtu-net.ru

Slovenia
Probo d.o.o.
Prekorje 48
SI-3211 Škofja vas
Phone: + 386 (0)3 4928741
Fax: + 386 (0)3 4928741
E-mail: gabrijel.popovic@probo.si

Singapore
Sciencewerke Pte Ltd
67 Ayer Rajah Crescent #04-21
Singapore 139950
Phone: +65 6777 1045
Fax: +65 6777 3054
E-mail: customerservice@sciencewerke.com

South Africa
Iepsa
279 Brooklyn Rd, Brooklyn
Pretoria
Phone: + 27 12 3622029
Fax: + 27 12 3622030
E-mail: iepamed@iepsamed.co.za

South Korea
MIRAX Co, Ltd
334-1 Cheoncheon-dong Jangan-gu
Suwon-si,
Gyeonggi-do, Korea (ZIP:440-710)
Tel : 82-31-268-0383
Fax : 82-31-268-0384
E-mail: info@mirax.co.kr

Spain
Vitro S.A.
Via de los Poblados, no 175a Planta
Naves 13 y 2
Edificio Indubuilding
28033 Madrid
Phone: + 34 913 82 16 20
Fax: + 34 917 63 66 68
E-mail: admvm@vitroweb.com

Sweden
Micro System
Ålvkvarnsvägen 165
S-163 51 SPÅNGA
Phone: + 46 8 754 08 40
Fax: + 46 8 754 0840
E-mail: info@microsystem.se

Syria
LeanGene
P.O.Box 697 Tlaa Alali
Amman, Jordan
Phone: + 962 6 46 55 401
Fax: + 962 6 46 55 402
E-mail: support@leangene-jo.com

Taiwan
Major Instruments Co., Ltd
9F: 69-3, Chung-Cheng E. Road, Sec. 2,
Tanshui, Taipei
Phone: +866-2-2808-1452
Fax: +866-2-2808-2354
E-mail: major@major.com.tw

Thailand
Dscience Company Ltd
301/354 Ramkamhaeng 68
Ramkamhaeng Rd, Huamark
Bangkapi
Bangkok 10240
Phone: +66 (2) 729 0239
Fax: +66 (2) 729 0239 ext 0
E-mail: info@dscience.co.th

Tunisia
Biomax Tunisie
63 Avenue Abou Kacem Chebbi,
Residence Omnia
Escalier C, N°6
2080 Ariana
Phone: +216 98349693
Fax: +216 71710590
E-mail: biomax.tunisie@planet.tn

Turkey
Tokra Medical Tic. Ltd. Sti.
682. Sokak No: 19
06370 Ostim
Ankara
Phone: + 90 312 3956009
Fax: + 90 312 3953961
E-mail: import@tokra.com.tr

Vietnam
Vinh Tien Technology Solutions Co., Ltd.
118/34/18/1 Lien Khu 5-6 Str.
Binh Hung Hoa B ward
Binh Tan District
Ho Chi Minh City
Phone: + 84 9 3327 8868
Fax: +84 (8) 6266 9209
E-mail: info@vitechins.com

HUGO GENE SYMBOLS

The table below gives the approved HUGO Gene Symbol and synonyms used for the genes in our REPEAT-FREE™ POSEIDON™ probes according to the HUGO Gene Nomenclature Committee (HGNC).

| <i>Chromosome position</i> | <i>Genes involved</i> | <i>HUGO Gene Symbols</i> | <i>Synonyms</i> | <i>Cat.# KBI</i> | <i>Page</i> |
|----------------------------|-----------------------|--------------------------|--|------------------|-------------|
| 1p36 | CHD5 | CHD5 | | 10507 | 41 |
| | | | | 10712 | 60 |
| 1q21 | S100A10 | S100A10 | 42C, "annexin II tetramer (Alt) p11 subunit", CLP11, P11 | 10503 | 39 |
| | | | | 10507 | 41 |
| 1q32 | MDM4 | MDM4 | MDMX | 10736 | 61 |
| 2p21 | EML4 | EML4 | ELP120, ROPP120 | 10746 | 54 |
| 2p23 | ALK | ALK | CD246 | 10746 | 54 |
| | | | | 10747 | 53,73 |
| 2p24 | MYCN | MYCN | bHLHe37, N-myc | 10706 | 59 |
| 2q11 | LAF | AFF3 | MLLT2-like | 10706 | 59 |
| 3p25 | PPARγ | PPARG | PPARG1, PPARG2, NR1C3, PPARGgamma | 10707 | 58 |
| 3q26 | EVI | MECOM | MDS1-EVI1, PRDM3 | 10204 | 26 |
| | | | | 10205 | 27 |
| 3q26 | hTERC | TERC | hTR, SCARNA19, "small Cajal body-specific RNA 19", TR, TRC3 | 10110 | 19 |
| | | | | 10204 | 26 |
| | | | | 10205 | 27 |
| | | | | 10704 | 52 |
| 3q27 | BCL6 | BCL6 | BCL5, BCL6A, LAZ3, ZBTB27 | 10607 | 45 |
| | | | | 10730 | 71 |
| 4p16 | FGFR3 | FGFR3 | CD333, CEK2, JTK4 | 10602 | 38,47 |
| 4p16 | WHSC1 | WHSC1 | MMSET, NSD2 | 40107 | 88 |
| 4q12 | CHIC2 | CHIC2 | BTL | 10003 | 14 |
| | | | | 10007 | 15 |
| 4q12 | FIP1L1 | FIP1L1 | DKFZp586K0717 | 10003 | 14 |
| | | | | 10007 | 15 |
| 4q12 | PDGFRA | PDGFRA | CD140a, PDGFR2 | 10003 | 14 |
| | | | | 10007 | 15 |
| 4q21-22 | AFF1 | AFF1 | AF-4, AF4 | 10404 | III |
| 5p15 | CTNND | CTNND2 | GT24, "neural plakophilin-related arm-repeat protein", NPRAP | 40106 | 88 |
| 5p15 | hTERT | TERT | EST2, hEST2, TCS1, TP2, TRT | 10208 | 23 |
| | | | | 10210 | 24 |
| | | | | 10709 | 55 |
| | | | | 40113 | 84 |
| 5q31 | CDC25C | CDC25C | PPP1R60, "protein phosphatase 1, regulatory subunit 60" | 10208 | 23 |
| | | | | 10209 | 24 |
| | | | | 10210 | 24 |
| | | | | 10709 | 55 |
| 5q31 | EGR1 | EGR1 | TIS8, G0S30, NGFI-A, KROX-24, ZIF-268, AT225, ZNF225 | 10208 | 23 |
| | | | | 10209 | 24 |
| | | | | 10210 | 24 |
| | | | | 10709 | 55 |
| 5q33 | CSF1R | CSF1R | C-FMS, CD115, CSFR | 10209 | 24 |
| | | | | 10210 | 24 |
| 5q33 | PDGFRB | PDGFRB | CD140b, JTK12, PDGFR1 | 10004 | 15 |
| | | | | 40064 | 77 |
| 5q33 | RPS14 | RPS14 | "40S ribosomal protein S14", "emetine resistance", EMTB, S14 | 10209 | 24 |
| | | | | 10210 | 24 |
| 5q35 | NSD1 | NSD1 | ARA267, FLJ22263, KMT3B | 40113 | 84 |
| 6p22 | DEK | DEK | D6S231E | 10306 | 32 |
| 6q21 | SEC63 | SEC63 | DNAJC23, ERdj2, PRO2507, SEC63L | 10105 | 18 |
| | | | | 10117 | 21 |
| | | | | 10504 | 40 |
| 6q22 | ROS1 | ROS1 | c-ros-1, MCF3, ROS | 10752 | V |
| 6q27 | MLLT4 | MLLT4 | AF-6, AF6 | 10309 | 31 |
| 7p11 | EGFR, Her1 | EGFR | ERBB1 | 10702 | 54 |
| 7q11 | ELN | ELN | "supravalvular aortic stenosis", SVAS, "tropoelastin", WBS, "Williams-Beuren syndrome", WS | 40111 | 86 |
| 7q11 | LIMK1 | LIMK1 | LIMK | 40111 | 86 |
| 7q22 | CUTL1 | CUX1 | CASP, CDP, CDP/Cut, CDP/Cux, CDP1, Clox, CUT, CUX, Cux/CDP, "golgi integral membrane protein 6", GOLIM6" | 10202 | 25 |
| | | | | 10207 | 25 |
| 7q31 | C-MET | MET | HGFR, RCCP2 | 10719 | 55 |
| 8p11 | FGFR1 | FGFR1 | BFGFR, BFGFR, CD331, CEK, FLG, H2, H3, H4, H5, N-SAM, "Pfeiffer syndrome" | 12754 | VI |
| | | | | 14754 | VI |
| | | | | 10737 | 16,51 |
| 8p21 | PNOC | PNOC | "nocistatin", PPNOC | 10503 | 39 |
| 8p23 | GATA4 | GATA4 | | 40118 | VII |

| <i>Chromosome position</i> | <i>Genes involved</i> | <i>HUGO Gene Symbols</i> | <i>Synonyms</i> | <i>Cat.# KBI</i> | <i>Page</i> |
|----------------------------|-----------------------|--------------------------|--|------------------|-------------|
| 8p21 | ETO | RUNX1T | CDR, ETO, MTG8, ZMYND2 | 10301 | 28 |
| 8q24 | C-MYC | MYC | bHLHe39, c-Myc | 10106 | 19 |
| | | | | 10117 | 21 |
| | | | | 10603 | 42 |
| | | | | 10611 | 47 |
| | | | | 10704 | 52 |
| 9p21 | p16 | CDKN2A | CDK4I, p16, INK4a, MTS1, CMM2, ARF, p19, p14, INK4, p16INK4a, p19Arf | 10402 | 34 |
| | | | | 10710 | 52 |
| 9p21 | MLLT3 | MLLT3 | AF-9, AF9, YEATS3 | 10308 | 30 |
| 9p24 | JAK2 | JAK2 | JTK10 | 10012 | 13 |
| 9q34 | ABL | ABL1 | JTK7, c-ABL, p150 | 10005 | 11,35 |
| | | | | 10006 | 11,35 |
| | | | | 10008 | 12,35 |
| | | | | 10009 | 12,35 |
| | | | | 10013 | II |
| 9q34 | ASS | ASS1 | CTLN1 | 10005 | 11,35 |
| | | | | 10006 | 11,35 |
| | | | | 10008 | 12,35 |
| | | | | 10013 | II |
| | | | | 10508 | 40 |
| 9q34 | NUP214 | NUP214 | CAIN, CAN, "CAN protein, putative oncogene", D9S46E, N214, "nuclear pore complex protein Nup214" | 10306 | 32 |
| 10q11 | RET | RET | "cadherin-related family member 16", CDHF12, CDHR16, PTC, RET51 | 10753 | V |
| 10q23 | PTEN | PTEN | MMAC1, "mutated in multiple advanced cancers 1", PTEN1, TEP1 | 10718 | 57 |
| 11p15 | NUP98 | NUP98 | NUP96 | 10311 | III |
| 11q13 | BCL1 | CCND1 | "B-cell CLL/lymphoma 1", "G1/S-specific cyclin D1", "parathyroid adenomatosis 1", U21B31 | 10604 | 43 |
| | | | | 10605 | 36,57 |
| | | | | 10609 | 46 |
| | | | | 10734 | 69 |
| 11q13 | MYEOV | MYEOV | OCIM | 10605 | 36,47 |
| 11q22 | ATM | ATM | TEL1, "TEL1, telomere maintenance 1, homolog (S. cerevisiae)", TEL01 | 10103 | 18 |
| | | | | 10108 | 21 |
| | | | | 10114 | 22 |
| 11q23 | MLL | MLL | ALL-1, CXXC7, HRX, HTRX1, KMT2A, MLL1A, TRX1 | 10303 | 29,35 |
| | | | | 10307 | 30 |
| | | | | 10308 | 30 |
| | | | | 10309 | 31 |
| | | | | 10404 | III |
| | | | | 10711 | 60 |
| 11q23 | PLZF | ZBTB16 | PLZF | 10502 | 38 |
| 12p13 | TEL | ETV6 | TEL, "TEL oncogene" | 10401 | 33 |
| | | | | 10403 | 34 |
| 12q13 | CDK4 | CDK4 | PSK-J3 | 10725 | 66 |
| 12q13 | CHOP | DDIT3 | "C/EBP zeta", CHOP, CHOP10, GADD153" | 10714 | 63 |
| 12q13 | GLI | GLI1 | | 10104 | 20 |
| | | | | 10108 | 21 |
| 12q15 | MDM2 | MDM2 | HDM2, HDMX, MGC5370 | 10717 | 65 |
| 13q14 | DLEU | DLEU1 | "B-cell neoplasia-associated gene with multiple splicing", LEU1, LINC00021, "long intergenic non-protein coding RNA 21", NCRNA00021, "non-protein coding RNA 21", XTP6 | 10102 | 17 |
| | | | | 10113 | 20 |
| | | | | 10502 | 38 |
| 13q14 | FKHR | FOXO1 | FKH1 | 10716 | 64 |
| 13q14 | RB7 | RB1 | RB | 40001 | 78 |
| | | | | 40003 | 79 |
| | | | | 40005 | 79 |
| | | | | 40050 | 77 |
| | | | | 40051 | 77 |
| | | | | 40064 | 77 |
| 14q32 | IGH | IGH@ | IGH | 10510 | IV |
| | | | | 10601 | 22,41,44 |
| | | | | 10602 | 38,47 |
| | | | | 10603 | 42 |
| | | | | 10604 | 43 |
| | | | | 10605 | 36,47 |
| | | | | 10606 | 43 |
| | | | | 10610 | 37 |
| | | | | 10729 | 70 |
| | | | | 10755 | VIII |


| <i>Chromosome position</i> | <i>Genes involved</i> | <i>HUGO Gene Symbols</i> | <i>Synonyms</i> | <i>Cat.# KBI</i> | <i>Page</i> |
|----------------------------|-----------------------|--------------------------|---|------------------|-------------|
| 15q11 | SNRPN | SNRPN | HCERN3, RT-LI, "SM protein N", SM-D, "small nuclear ribonucleoprotein N", SMN, SNRNP-N, SNURF-SNRPN, "tissue-specific splicing protein" | 40109 | 85 |
| 15q11 | UBE3A | UBE3A | hUba3, "UBA3, ubiquitin-activating enzyme E1 homolog (yeast)" | 40110 | 86 |
| | | | | 40116 | 91 |
| 15q22 | SMAD6 | SMAD6 | HsT17432 | 10504 | 40 |
| | | | | 10508 | 40 |
| 15q24 | PML | PML | MYL, RNF71, TRIM19" | 10302 | 29 |
| | | | | 40109 | 85 |
| | | | | 40110 | 86 |
| 15q26 | IGF1R | IGF1R | CD221, IGFIR, IGFR, JTK13, MGC18216 | 40116 | 91 |
| 16p11 | FUS | FUS | FUS1, "heterogeneous nuclear ribonucleoprotein P2", hnRNP-P2, TLS, "translocated in liposarcoma" | 10715 | 64 |
| 16q22 | CBFB | CBFB | PEBP2B | 10304 | 32,35 |
| 16q23 | MAF | MAF | c-MAF | 10610 | 37 |
| 17p11 | RAI1 | RAI1 | DKFZP434A139, KIAA1820, MGC12824, SMS | 40101 | 87 |
| 17p13 | AURKB | AURKB | Aik2, AIM-1, ARK2, AurB, "aurora-1", "aurora-B", IPL1, PPP1R48, "protein phosphatase 1, regulatory subunit 48", STK5 | 10722 | 68 |
| 17p13 | LIS | PAFAH1B1 | LIS1, PAFAH | 40101 | 87 |
| 17p13 | p53 | TP53 | LFS1, "Li-Fraumeni syndrome", p53 | 10011 | 13 |
| | | | | 10112 | 17 |
| | | | | 10113 | 20 |
| | | | | 10114 | 22 |
| | | | | 10509 | 39 |
| | | | | 10738 | 72 |
| 17q11 | NF1 | NF1 | "neurofibromatosis", "von Recklinghausen disease", "Watson disease" | 40114 | 84 |
| 17q12 | ERBB2 | ERBB2 | CD340, HER-2, HER2, NEU | 10701 | 49 |
| | | | | 10735 | 50 |
| 17q21 | COL1A1 | COL1A1 | OI4 | 10742 | VI |
| 17q21 | RARA | RARA | NR1B1, RAR | 10302 | 29 |
| | | | | 10305 | 31 |
| 17q21 | TOP2A | TOP2A | | 10724 | 49 |
| | | | | 10735 | 50 |
| 17q22 | MPO | MPO | | 10011 | 13 |
| | | | | 40114 | 84 |
| 18q11 | SYT | SS18 | SYT | 10713 | 63 |
| 18q21 | BCL2 | BCL2 | Bcl-2, PPP1R50, "protein phosphatase 1, regulatory subunit 50" | 10606 | 43 |
| | | | | 10612 | 44 |
| | | | | 10745 | 70 |
| | | | | 10755 | VIII |
| 18q21 | MALT | MALT1 | "paracaspase" | 10608 | 45 |
| | | | | 10731 | 71 |
| 19p13 | MLLT1 | MLLT1 | ENL, LTG19, YEATS1 | 10307 | 30 |
| 19p13 | ZNF443 | ZNF443 | ZK1 | 10739 | 66 |
| 19q13 | CD37 | CD37 | TSPAN26 | 10509 | 39 |
| | | | | 40064 | 77 |
| 19q13 | ERCC1 | ERCC1 | RAD10 | 10739 | 66 |
| 20q11 | MAPRE1 | MAPRE1 | "adenomatous polyposis coli-binding protein EB1", EB1 | 10203 | 26 |
| | | | | 10721 | 67 |
| | | | | 10733 | 50 |
| 20q12 | MAFB | MAFB | "Kreiser (mouse) maf-related leucine zipper homolog", KR" | 10510 | IV |
| 20q12 | PTPRT | PTPRT | KIAA0283, RPTPrho | 10203 | 26 |
| 20q13 | AURKA | AURKA | AIK, ARK1, AurA, "Aurora-A kinase", BTAk, PPP1R47, "protein phosphatase 1, regulatory subunit 47", STK7 | 10721 | 67 |
| 20q13 | NFATC2 | NFATC2 | NF-ATP, NFAT1, NFATp | 10751 | VII |
| 20q13 | ZNF217 | ZNF217 | ZABC1 | 10733 | 50 |
| 21q22 | AML | RUNX1 | "aml1 oncogene", AMLCR1, PEBP2A2 | 10301 | 28 |
| 21q22 | DSCR1 | RCAN1 | | 40002 | 79 |
| | | | | 40003 | 79 |
| | | | | 40005 | 79 |
| | | | | 40006 | 79 |
| | | | | 40007 | 79 |
| | | | | 40008 | 79 |
| | | | | 40050 | 77 |
| | | | | 40051 | 77 |
| 21q22 | ERG | ERG | erg-3, p55, "TMPRSS2-ERG prostate cancer specific", "transcriptional regulator ERG (transforming protein ERG)", "v-ets avian erythroblastosis virus E26 oncogene related", "v-ets erythroblastosis virus E26 oncogene like" | 10726 | 57 |

| <i>Chromosome position</i> | <i>Genes involved</i> | <i>HUGO Gene Symbols</i> | <i>Synonyms</i> | <i>Cat.# KBI</i> | <i>Page</i> |
|----------------------------|-----------------------|--------------------------|--|------------------|-------------|
| 21q22 | TMPRSS2 | TMPRSS2 | PRSS10 | 10726 | 57 |
| 22q11 | BCR | BCR | ALL, CML, D22S662, PHL | 10005 | 11,35 |
| | | | | 10006 | 11,35 |
| | | | | 10008 | 12,35 |
| | | | | 10009 | 12,35 |
| | | | | 10013 | II |
| 22q11 | CLH22 | CLTCL1 | CHC22, CLH22, CLTD | 40102 | 81,83 |
| 22q11 | DGCR2 | DGCR2 | DGS-C, IDD, "integral membrane protein DGCR2", KIAA0163, LAN, SEZ-12 | 40102 | 81,83 |
| 22q11 | TUPLE | HIRA | DGCR1, "DiGeorge critical region gene 1", TUP1 | 40103 | 81,83 |
| 22q11 | TBX1 | TBX1 | | 40104 | 81,83 |
| 22q12 | EWSR1 | EWSR1 | EWS | 10750 | 62 |
| | | | | 10751 | VII |
| 22q13 | ARSA | ARSA | "metachromatic leucodystrophy" | 40102 | 81,83 |
| 22q13 | PDGFB | PDGFB | "becaplermin", "oncogene SIS", SSV | 10742 | VI |
| 22q13 | SHANK3 | SHANK3 | KIAA1650, "proline rich synapse associated protein 2", prosap2, PSAP2, "shank postsynaptic density protein", SPANK-2 | 40102 | 81,83 |
| | | | | 40103 | 81,83 |
| | | | | 40104 | 81,83 |
| Xp11 | TFE3 | TFE3 | bHLHe33, member A, TFEA, transcription factor E family | 10741 | 72 |
| | | | | 40051 | 77 |
| Xp22 | KAL | KAL1 | "anosmin-1", KALIG-1, "WAP four-disulfide core domain 19", WFDC19 | 40115 | 90 |
| Xp22 | SHOX | SHOX | GCFX, PHOG, SHOXY, SS | 40112 | 89 |
| Xp22 | STS | STS | ARSC, "arylsulfatase C" | 40115 | 90 |
| Xq12 | AR | AR | AIS, HUMARA, "Kennedy disease", NR3C4, SMAX1, "testicular feminization" | 10720 | 56 |
| Xq13 | XIST | XIST | DXS1089, LINC00001, "long intergenic non-protein coding RNA 1", NCRNA00001, swd66 | 40108 | 89 |



LABEL SMARTER using ULS™

- SPEED - nucleic acid labeling in just 30 mins
- PRECISION - perfect for FFPE when labeling degraded DNA
- ACCURACY - no enzymatic bias



**ONE universal labeling system
compatible with all your
biological samples.**
Enquire today!




Phone: +31 (0)20 6919181
Fax: +31 (0)20 630 42 47


Email: customerservice@kreatech.com
Website: www.kreatech.com

Kreatech Diagnostics
Vlierweg 20
1032 LG Amsterdam
The Netherlands


REPEAT-FREE™ POSEIDON™




KBI-10068
1q21 deletion
Material kindly provided by Dr. Balogh, Budapest



KBI-10102
13q14 deletion
Image kindly provided by Dr. Dastgheib, Toulouse



KBI-10209
7p11 deletion
Image kindly provided by Dr. Mohr, Dresden















KBI-10306
13q14 deletion
Image kindly provided by Dr. Stevens-Kroef, Nijmegen



KBI-10404
13q14 deletion
Image kindly provided by Prof. Christine J. Harrison, Newcastle



KBI-10611
19p13 deletion
Image kindly provided by Prof. Siebert, Kiel



For availability of the products highlighted in orange inquire at customerservice@kreatech.com

Kreatech offers the broadest commercial portfolio available of Fluorescent in Situ Hybridization (FISH) probes: more than 400 REPEAT-FREE™ POSEIDON™ probes. In addition, Kreatech has a product portfolio for microarray labeling applications using Kreatech's proprietary non-enzymatic Universal Linkage System (ULS™) labeling technology.

International
KREATECH Diagnostics
Vlierweg 20, 1032 LG Amsterdam
The Netherlands

Phone: +31 (0)20 691 9181
Fax: +31 (0)20 630 4247
E-mail: customerservice@kreatech.com
www.kreatech.com



PROBES — NEXT GENERATION FISH

CAN'T FIND THE FISH PROBE YOU'RE LOOKING FOR?
INQUIRE ABOUT A CUSTOM-MADE FISH PROBE AT FISH4U@KREATECH.COM

Cat. # KBI- Page

| | | |
|---------------|-------|-------|
| 4pter | 40207 | 92 |
| FGFR3 4p16 | 10602 | 38,47 |
| WHSC1 4p16 | 40107 | 88 |
| PPARGC1A 4p15 | 20004 | 94 |
| SE4 D4Z1 | 40061 | 77 |
| FIP1L1-CHIC2 | 10003 | 14 |
| PDGFRA 4q12 | 10007 | 15 |
| KIT-KDR 4q12 | | |
| AFF1 4q21-22 | 10404 | III |
| TET2 4q24 | | |
| 4qter | 40208 | 92 |

4

| | | |
|-------------|-------|----|
| 10pter | 40219 | 92 |
| DGSII 10p14 | 40105 | 83 |
| KIF5B 10p11 | | |
| SE10 D10Z1 | 20010 | 94 |
| RET 10q11 | 40062 | 77 |
| | 10753 | V |
| PTEN 10q23 | 10718 | 57 |
| GOT1 10q24 | | |
| FGFR2 10q26 | | |
| 10qter | 40220 | 92 |

10

| | | |
|------------|-------|----|
| 16pter | 40228 | 92 |
| FUS 16p11 | 10715 | 64 |
| SE16 D16Z2 | 20016 | 94 |
| | 40050 | 77 |
| | 40063 | 77 |
| CBFB 16q22 | 10304 | 32 |
| MAF 16q23 | 10610 | 37 |
| 16qter | 40229 | 92 |

16

| | | |
|---------------|-------|-------|
| NOR | 20033 | 95 |
| SE14/22 | 20028 | 94 |
| DGCR 'N25' | 40102 | 81,83 |
| 22q11 | | |
| DGCR 'Tuplet' | 40103 | 81,83 |
| HIRA 22q11 | | |
| TBX1 22q11 | 40104 | 82,83 |
| BCR 22q11 | 10005 | 11,35 |
| | 10006 | 11,35 |
| | 10008 | 12,35 |
| | 10009 | 12,35 |
| | 40050 | 77 |
| | 40061 | 77 |
| EWSR1 22q12 | 10750 | 62 |
| PDGFB 22q13 | 10751 | VII |
| SHANK3 22q13 | 40102 | 81,83 |
| | 40103 | 81,83 |
| | 40104 | 82,83 |
| 22qter | 40239 | 92 |

22

Cat. # KBI- Page

| | | |
|--------------|-------|----|
| 5pter | 40209 | 92 |
| TERT 5p15 | 10208 | 23 |
| | 10210 | 24 |
| | 10709 | 55 |
| | 40113 | 84 |
| CTNND2 5p15 | 40106 | 88 |
| SE1/5/19 | 20026 | 94 |
| Sq11 Control | | |
| CDC25C 5q31 | 10208 | 23 |
| EGR1 5q31 | 10209 | 24 |
| | 10210 | 24 |
| | 10709 | 55 |
| CSF1R 5q33 | 10209 | 24 |
| | 10210 | 24 |
| PDGFRB 5q33 | 10004 | 15 |
| | 40064 | 77 |
| RPS14 5q33 | 10209 | 24 |
| NPM1 5q35 | 10210 | 24 |
| FGFR4 5q35 | | |
| NSD1 5q35 | 40113 | 84 |
| 5qter | 40210 | 92 |

5

| | | |
|--------------|-------|-------|
| 11pter | 40221 | 92 |
| NUP98 11p15 | 10311 | III |
| SE11 D11Z1 | 20011 | 94 |
| | 40062 | 77 |
| MYEOV 11q13 | 10605 | 36,47 |
| CCND1 11q13 | 10604 | 43 |
| | 10609 | 46 |
| | 10734 | 69 |
| ATM 11q22 | 10103 | 18 |
| | 10108 | 21 |
| | 10114 | 22 |
| ZNF145 11q23 | 10502 | 38 |
| | 10307 | 30 |
| | 10308 | 30 |
| | 10309 | 31 |
| MLL 11q23 | 10303 | 29,35 |
| | 10307 | 30 |
| | 10308 | 30 |
| | 10309 | 31 |
| | 10404 | III |
| | 10711 | 60 |
| | 40222 | 92 |

11

| | | |
|---------------|-------|----|
| 17pter | 40230 | 92 |
| LIS 17p13 | 40101 | 87 |
| USP6 17p13 | | |
| TP53 17p13 | 10011 | 13 |
| | 10112 | 17 |
| | 10113 | 20 |
| | 10114 | 22 |
| | 10509 | 39 |
| | 10738 | 72 |
| AURKB 17p13 | 10722 | 68 |
| RAI1 17p11 | 40101 | 87 |
| 17p11.2 break | | |
| SE17 D17Z1 | 20017 | 94 |
| | 40062 | 77 |
| NF1 17q11 | 40114 | 84 |
| ERBB2 17q12 | 10701 | 49 |
| | 10735 | 50 |
| RARA 17q21 | 10302 | 29 |
| | 10305 | 31 |
| TOP2A 17q21 | 10724 | 49 |
| BRCA1 17q21 | 10735 | 50 |
| COL1A1 17q21 | 10742 | VI |
| MPO 17q22 | 10011 | 13 |
| TBX2 17q23 | 40114 | 84 |
| FASN 17q25 | | |
| 17qter | 40231 | 92 |

17

| | | |
|-------------|-------|----|
| XYpter | 40240 | 92 |
| SHOX Xp22 | 40112 | 89 |
| STS Xp22 | 40115 | 90 |
| KAL1 Xp22 | 40115 | 90 |
| TFE3 Xp11 | 10741 | 72 |
| | 40051 | 77 |
| SEX DXZ1 | 20023 | 94 |
| | 20030 | 94 |
| | 20032 | 94 |
| | 40005 | 79 |
| | 40008 | 79 |
| | 40063 | 77 |
| AR Xq12 | 10720 | 56 |
| XIST Xq13 | 40108 | 89 |
| MAGEC1 Xq27 | | |
| XYqter | 40241 | 92 |

X

Cat. # KBI- Page

| | | |
|------------------|-------|----|
| 6pter | 40211 | 92 |
| IRF4 DUSP22 6p25 | | |
| RREB1 6p24 | | |
| DEK 6p22 | 10306 | 32 |
| CCND3 6p21.1 | | |
| SE6 D6Z1 | 20006 | 94 |
| | 40061 | 77 |
| SEC63 6q21 | 10105 | 18 |
| | 10117 | 21 |
| | 10504 | 40 |
| ROS1 6q22 | 10752 | V |
| MYB 6q23 | | |
| TNFAIP3 6q23 | | |
| MLLT4 6q27 | 10309 | 31 |
| 6qter | 40212 | 92 |

6

| | | |
|-------------|-------|-------|
| 12pter | 40223 | 92 |
| NOL1 12p13 | | |
| ETV6 12p13 | 10401 | 33 |
| KRAS 12p12 | 10403 | 34 |
| SE12 D12Z3 | 20012 | 22,94 |
| ATF1 12q13 | 40062 | 77 |
| HER3 12q13 | | |
| GLI1 12q13 | 10104 | 20 |
| | 10108 | 21 |
| DDIT3 12q13 | 10714 | 63 |
| CDK4 12q14 | 10725 | 66 |
| MDM2 12q15 | 10717 | 65 |
| 12qter | 40224 | 92 |


12

| | | |
|-------------|-------|----|
| 18pter | 40232 | |
| TYMS 18p11 | | 92 |
| SE18 D18Z1 | 20018 | 94 |
| | 20032 | 94 |
| | 40005 | 94 |
| | 40050 | 77 |
| | 40051 | 77 |
| | 40063 | 77 |
| SS18 18q11 | 10713 | 63 |
| MALT1 18q21 | 10608 | 45 |
| DCC1 18q21 | 10731 | 71 |
| BCL2 18q21 | 10606 | 43 |
| | 10612 | 44 |
| | 10745 | 70 |
| | 10755 | 70 |
| 18qter | 40233 | 92 |


18

| | | |
|-------------|-------|----|
| XYpter | 40240 | 92 |
| SHOX YP11.3 | 40112 | 89 |
| SEY DYZ3 | 20024 | 94 |
| | 20030 | 94 |
| | 20032 | 94 |
| | 40005 | 79 |
| | 40008 | 79 |
| | 40063 | 77 |
| AZFc Yq11 | 40051 | 77 |
| SEY DYZ1 | 20025 | 94 |
| XYqter | 40241 | 92 |


Y




Amplification of Her2/Neu (ERBB2) & SE 17




KBI-10747
ALK 2p23 rearrangement
Image kindly provided by Prof. B. Tarris, Dr. P.A. Just, Paris




KBI-10752
ROS1 6q22 rearrangement



KBI-12754/14754
Amplification of FGFR1 3p11
Image kindly provided by Dr. Marcus Otto, ORDIS Biomarkers, Vienna



KBI-10717
Amplification of MDM2 12q15
Image kindly provided by W. Grifflinger, Amstern



KBI-10750
EWSR1 22q12 rearrangement

Patents and Trademarks

These products or the use of these products is subject to intellectual property rights. FISH and CISH probes are produced under an exclusive license from Veridex LLC for its REPEAT-FREE™ technology. FISH probes are labeled with the Kreatech's proprietary Universal Linkage system (ULS™). The fluorophores used in the PlatinumBright™ Nucleic Acid Labeling Kits (-415; -547; -590;

647) are subject to patents and know-how, owned or controlled, and manufactured by DYOMICS GmbH. The fluorophores used in the ULS™-Cys and ULS™-Cys3 labeling compounds are subject to patents, owned or controlled, and manufactured by GE Healthcare. The ULS™ technology and products are covered by one or more of the following US patents and corresponding patent claims outside the US owned by KREATECH: US 5,580,990; 5,714,327; 5,985,566;

6,133,038; US RE 40,557E; 6,797,818; 7,217,813. KREATECH is a trademark of Kreatech Biotechnology B.V.; ULS™, PlatinumBright™, KreaPure™, KreaBlock™ and POSEIDON™ are trademarks of KREATECH Biotechnology B.V. REPEAT-FREE™ is a trademark of Veridex LLC. ThermoBrite™ and CytoFuge™ 2 are trademarks of StatSpin, a division of IRIS Sample Processing.

Nothing disclosed herein is to be construed as a recommendation to use any of KREATECH products in violation of any patents. Furthermore, KREATECH does not warrant that the product, either alone, or in combination with other products, is immune from charges of patent infringement. KREATECH will not be held responsible for patent infringements or other violations that may occur from the use of any of its products.


INDEX BY PROBE NAME

| Description | Cat# | Page |
|---|-----------|----------|
| 6q21 / MYC (8q24) | KBI-10117 | 21 |
| 6q21 / SE 6 | KBI-10105 | 18 |
| ALK (2p23) Break | KBI-10747 | 53,73 |
| ALK/EML4 t(2;2); inv(2) Fusion | KBI-10746 | 54 |
| All Mouse Centromere (AMC) | KBI-30500 | 100 |
| AML/ETO t(8;21) Fusion | KBI-10301 | 28 |
| Angelman UBE3A (15q11) / PML(15q24) | KBI-40110 | 86 |
| AR (Xq12) / SE X | KBI-10720 | 56 |
| Arm Specific DNA Probes | | 89 |
| ATM (11q22) / GLI (12q13) | KBI-10108 | 21 |
| ATM (11q22) / SE 11 | KBI-10103 | 18 |
| AURKA (20q13) / 20q11 | KBI-10721 | 68 |
| AURKB (17p13) / SE 17 | KBI-10722 | 67 |
| Band Specific DNA Probes | | 89 |
| BCL1/IGH t(11;14) Fusion | KBI-10604 | 43 |
| BCL2 (18q21) Break | KBI-10612 | 44 |
| BCL2 (18q21) Break (tissue) | KBI-10745 | 70 |
| BCL2/IGH t(14;18) Fusion | KBI-10606 | 43 |
| BCL2/IGH@ t(14;18) Fusion (tissue) | KBI-10755 | IX |
| BCL6 (3q27) Break | KBI-10607 | 45 |
| BCL6 (3q27) Break (tissue) | KBI-10730 | 71 |
| BCR/ABL t(9;22), DC, S-Fusion, ES | KBI-10008 | 12,35 |
| BCR/ABL t(9;22) DC, S-Fusion | KBI-10009 | 12,35 |
| BCR/ABL t(9;22) Fusion | KBI-10005 | 35 |
| BCR/ABL t(9;22) Fusion | KBI-12005 | 11 |
| BCR/ABL t(9;22) TC, D-Fusion | KBI-10006 | 11,35 |
| C-MET (7q31) / SE 7 | KBI-10719 | 55 |
| C-MYC (8q24) / SE 8 | KBI-10106 | 19 |
| CBFB t(16;16), inv(16) Break | KBI-10304 | 32,35 |
| CC hTERC (3q26) C-MYC (8q24) / SE 7 TC | KBI-10704 | 52 |
| CCND1 (11q13) / SE 11 | KBI-10734 | 69 |
| CCND1 (BCL1;11q13) Break | KBI-10609 | 46 |
| CDK4 (12q13) / SE 12 | KBI-10725 | 66 |
| CHOP (12q13) Break | KBI-10714 | 63 |
| COL1A1/PDGFB t(17;22) DC, S-Fusion | KBI-10742 | VII |
| Cri-Du-Chat CTNND (5p15) / 5q31 | KBI-40106 | 88 |
| DEK / NUP214 t(6;9) Fusion | KBI-10306 | 32 |
| DiGeorge "N25" (22q11) / 22q13 (SHANK3) | KBI-40102 | 81,85 |
| DiGeorge II (10p14) / SE 10 | KBI-40105 | 83 |
| DiGeorge T-box1 (22q11) / 22q13 (SHANK3) | KBI-40104 | 82,85 |
| DiGeorge Tuple (22q11) / 22q13 (SHANK3) | KBI-40103 | 81,85 |
| DLEU (13q14) / 13qter | KBI-10102 | 17 |
| DLEU (13q14) / p53 (17p13) | KBI-10113 | 20 |
| EGFR, Her-1 (7p11) / SE 7 | KBI-10702 | 54 |
| ERBB2, Her-2/Neu (17q12) / SE 17 | KBI-10701 | 49 |
| ERBB2, Her-2/Neu (17q12) / SE 17 | KBI-14701 | 49 |
| ERCC1 (19q13) / ZNF443 (19p13) | KBI-10739 | 66 |
| ETV6 (TEL) (12p13) Break | KBI-10403 | 34 |
| EVI t(3;3); inv(3) (3q26) Break | KBI-10204 | 26 |
| EVI t(3;3); inv(3) (3q26) Break, TC | KBI-10205 | 27 |
| EWSR1 (22q12) Break | KBI-10750 | 62 |
| EWSR1 (22q12) Break | KBI-10750 | VIII |
| EWSR1/NFATC2 t(20;22) DC, S-Fusion | KBI-10751 | VIII |
| FGFR1 (8p11) / SE 8 (D8Z1) | KBI-12754 | VII |
| FGFR1 (8p12) Break | KBI-10737 | 16,51 |
| FGFR3/IGH t(4;14) Fusion | KBI-10602 | 38,47 |
| FIP1L1-CHIC2-PDGFR4 (4q12) Del, Break | KBI-10003 | 14 |
| FIP1L1-CHIC2-PDGFR4 (4q12) Del, Break, TC | KBI-10007 | 15 |
| FKHR (13q14) Break | KBI-10716 | 64 |
| FUS (16p11) Break | KBI-10715 | 64 |
| GLI (12q13) / SE 12 | KBI-10104 | 20 |
| hTERC (3q26) / 3q11 | KBI-10110 | 19 |
| hTERT (5p15) / 5q31 | KBI-10208 | 23 |
| hTERT (5p15) / 5q31 (tissue) | KBI-10709 | 55 |
| IGF1R (15q26) / 15q11 | KBI-40116 | 91 |
| IGH (14q32) Break | KBI-10601 | 22,41,44 |
| IGH (14q32) Break (tissue) | KBI-10729 | 70 |
| JAK2 (9p24) Break | KBI-10012 | 13 |
| MAF/IGH t(14;16) Fusion | KBI-10610 | 37 |
| MAFB/IGH@ t(14;20) Fusion | KBI-10510 | V |
| MALT (18q21) Break | KBI-10608 | 45 |
| MALT (18q21) Break (tissue) | KBI-10731 | 71 |
| MD GATA4 (8p23) / SE 8 | KBI-40118 | X |
| MDM2 (12q15) / SE 12 | KBI-10717 | 65 |
| MDM4 (1q32) / SE 1 | KBI-10736 | 61 |
| MDS 20q- (PTPRT 20q12) / 20q11 | KBI-10203 | 26 |
| MDS 5q- (5q31; 5q33) | KBI-10209 | 24 |
| MDS 5q- (5q31; 5q33) / hTERT (5p15) TC | KBI-10210 | 24 |
| MDS 7q- (7q22; 7q36) | KBI-10202 | 25 |

| Description | Cat# | Page |
|---|-------------|-------|
| MDS 7q- (7q22; 7q36) / SE 7 TC | KBI-10207 | 25 |
| MECOM / RUNX1 t(3;21) Fusion | KBI-10310 | II |
| Miller-Dieker LIS (17p13) / Smith-Magenis RAI (17p11) | KBI-40101 | 87 |
| MLL (11q23) / SE 11 | KBI-10711 | 60 |
| MLL (11q23) Break | KBI-10303 | 29 |
| MLL/AFF1 t(4;11) Fusion | KBI-10404 | IV |
| MLL/MLLT1 t(11;19) Fusion | KBI-10307 | 30 |
| MLL/MLLT3 t(9;11) Fusion | KBI-10308 | 30 |
| MLL/MLLT4 t(6;11) Fusion | KBI-10309 | 31 |
| MM 11q23 / DLEU (13q14) | KBI-10502 | 38 |
| MM 15q22 / 6q21 | KBI-10504 | 40 |
| MM 15q22 / 9q34 | KBI-10508 | 40 |
| MM 19q13 / p53 (17p13) | KBI-10509 | 39 |
| MM 1q21 / 8p21 | KBI-10503 | 39 |
| MM 1q21 / SRD (1p36) | KBI-10507 | 41 |
| Mm-BCR/ABL t(9;22), DC, S-Fusion, ES | KBI-10013 | II |
| MYC (8q24) Break, TC | KBI-10611 | 47 |
| MYC (8q24) Break, TC (tissue) | KBI-10749 | IX |
| MYC/IGH t(8;14) Fusion | KBI-10603 | 42 |
| MYCN (2p24) / LAF (2q11) | KBI-10706 | 59 |
| MYEOV/IGH t(11;14) Fusion | KBI-10605 | 36,47 |
| NF1 (17q11) / MPO (17q22) | KBI-40114 | 84 |
| NSD1 (5q35)/ hTERT (5p15) | KBI-40113 | 84 |
| NUP98 (11p15) Break | KBI-10311 | III |
| p16 (9p21) / 9q21 | KBI-10402 | 34 |
| p16 (9p21) / 9q21 (tissue) | KBI-10710 | 52 |
| p53 (17p13) / ATM (11q22) | KBI-10114 | 22 |
| p53 (17p13) / MPO (17q22) "ISO 17q" | KBI-10011 | 13 |
| p53 (17p13) / SE 17 | KBI-10112 | 17 |
| p53 (17p13) / SE 17 | KBI-12112 | 17 |
| p53 (17p13) / SE 17 (tissue) | KBI-10738 | 72 |
| PDGFRB (5q33) Break | KBI-10004 | 15 |
| PloidyScreen (21, X, Y) | KBI-40008 | 79 |
| PML/RARA t(15;17) Fusion | KBI-10302 | 29 |
| PML/RARA t(15;17) Fusion | KBI-12302 | 29 |
| PN 13 (13q14) | KBI-40001 | 78 |
| PN 13 (13q14) / 21 (21q22) | KBI-40003 | 79 |
| PN 21 (21q22) | KBI-40002 | 78 |
| PPARg (3p25) Break | KBI-10707 | 58 |
| Prader-Willi SNRNP (15q11) / PML(15q24) | KBI-40109 | 85 |
| PrenatScreen (13/21, X/Y18) | KBI-40005 | 79 |
| PrenatScreen (13/21, X/Y18) | KBI-40006 | 79 |
| PrenatScreen (13/21, X/Y18) | KBI-40007 | 79 |
| PTEN (10q23) / SE 10 | KBI-10718 | 57 |
| RAB9B (XqF1) / DSCR (16qC4) | KBI-30503 | 100 |
| RAB9B (XqF1) / WC Y | KBI-30505 | 100 |
| RARA (17q21) Break | KBI-10305 | 31 |
| RET (10q11) Break | KBI-10753 | VI |
| ROS1 (6q22) Break | KBI-10752 | VI |
| Satellite Enumeration DNA Probes | | 94,95 |
| SE (X,Y,18) | KBI-20032 | 79 |
| SE 18 (D18Z1) 5x conc | KBI-20018-B | 79 |
| SE 7 (D7Z1) / SE 8 (D8Z1) | KBI-20031 | 79 |
| SE X (DXZ1) / SE Y (DYZ3) | KBI-20030 | 79 |
| Short Stature (Xp22) / SE X | KBI-40112 | 89 |
| SRD (1p36) / SE 1(1qh) | KBI-10712 | 60 |
| STS (Xp22) / KAL (Xp22) / SE X TC | KBI-40115 | 90 |
| Sub-Telomere DNA Probes | | 92 |
| SYT (18q11) Break | KBI-10713 | 63 |
| TEL/AML t(12;21) Fusion | KBI-10401 | 33 |
| TFE3 (Xp11) Break | KBI-10741 | 72 |
| TK (11qE1) / AurKa (2qH3) | KBI-30501 | 100 |
| TK (11qE1) / WC Y | KBI-30502 | 100 |
| TPRSS2-ERG (21q22) Del, Break, TC | KBI-10726 | 57 |
| TOP2A (17q21) / Her-2/neu (17q12) / SE 17 TC | KBI-10735 | 50 |
| TOP2A (17q21) / SE 17 | KBI-10724 | 49 |
| TwinStar C-MET (7q31) | KBI-60719 | 74 |
| TwinStar CISH Detection Kit | KBI-60010 | 74 |
| TwinStar EGFR (7p11) | KBI-60702 | 74 |
| TwinStar Her2/neu (17q12) | KBI-60701 | 74 |
| UniStar C-MET (7q31) | KBI-50719 | 74 |
| UniStar CISH Detection Kit | KBI-50001 | 74 |
| UniStar EGFR (7p11) | KBI-50702 | 74 |
| UniStar Her2/neu (17q12) | KBI-50701 | 74 |
| Whole Chromosome DNA Probes | | 96 |
| Williams-Beuren ELN (7q11) / 7q22 | KBI-40111 | 86 |
| Wolf-Hirschhorn WHSC1 (4p16) / SE 4 | KBI-40107 | 88 |
| X-Inactivation XIST (Xq13) / SE X | KBI-40108 | 89 |
| ZNF217 (20q13) / 20q11 | KBI-10733 | 50 |

REPEAT-FREE™ POSEIDON™ DNA FISH PROBES: ADVANCING FISH


REPEAT-FREE™ POSEIDON™ FISH PROCEDURE



REPEAT-FREE™ (RF) POSEIDON™ FISH PROBES are the latest advancements in DNA FISH probes. The RF technology used to construct the RF POSEIDON™ FISH probes is based on subtractive hybridization that specifically removes all repetitive elements that are dispersed throughout the target area of interest. Eliminating these repeat sequences leads to more specific binding kinetics and eliminates the need to use Cot1 DNA during pre-annealing.

This results in FISH assays with reduced backgrounds and brighter signals. Traditional FISH is time consuming with standard hybridizations being 16 hours on average. With the use of RF POSEIDON™ probes, hybridization may be shortened to as little as 4 hours depending on the combination of probe and sample type used. RF technology can have a positive impact on the readability, sensitivity and specificity of your FISH assays.

CONVENTIONAL FISH PROCEDURE



ORDERING INFORMATION

How to place an order

Orders are accepted by phone, fax, e-mail or mail. When placing an order please include the following information:

- Your account number (required)
- Your name and phone number
- Complete billing address
- Complete shipping address
- Purchase order number
- Catalogue number and description of product
- Quantity needed
- VAT number (where applicable)

If it is necessary to send a written order confirmation by mail or fax, please indicate clearly that your order is marked "CONFIRMATION". This will assist in avoiding shipment duplications. KREATECH cannot be held responsible in cases where the latter procedure is not followed.

We accept Visa and MasterCard

If you wish to pay with a credit card please use the fax or e-mail option to send your order and please provide the name of your credit card, account number, expiration date and security code when placing credit card orders.

International

KREATECH Diagnostics
Vlierweg 20
1032 LG Amsterdam
The Netherlands

Phone: +31 (0)20 691 9181
Fax: +31 (0)20 630 4247

E-mail: customerservice@kreatech.com
Website: www.kreatech.com

Benelux

KREATECH Diagnostics
Vlierweg 20
1032 LG Amsterdam
The Netherlands

Phone: +31 (0)6 4850 0107
Fax: +31 (0)35 656 4806

E-mail: customerservice@kreatech.com
Website: www.kreatech.com

France

KREATECH Diagnostics
20 Avenue de la Paix
67080 Strasbourg Cedex
France

Phone: +33 (0)1 4372 0079
Fax: +33 (0)1 4348 8244

E-mail: customerservice@kreatech.com
Website: www.kreatech.com

Germany

KREATECH Diagnostics
Vlierweg 20
1032 LG Amsterdam
The Netherlands

Phone: +49 (0)223 3713 5979
Fax: +31 (0)20 630 4247

E-mail: customerservice@kreatech.com
Website: www.kreatech.com

United Kingdom

KREATECH Diagnostics
52 New Town, Uckfield
East Sussex, TN22 5DE
United Kingdom

Phone: +44 (0)208 350 543
Fax: +44 (0)208 711 3132

E-mail: customerservice@kreatech.com
Website: www.kreatech.com
